

USER MANUAL

HS898 MOBILITY SCOOTER



TABLE OF CONTENTS

INTRODUCTION1
IMPORTANT PRECAUTIONS2
ELECTROMAGNETIC INTERFERENCE AND WARNINGS3
SAFETY WARNING AND INSTRUCTION LABELS5
IDENTIFICATION OF PARTS6
CHARGING THE BATTERIES11
DISASSEMBLING YOUR SCOOTER12
ASSEMBLING YOUR SCOOTER13
OPERATING YOUR SCOOTER13
CARE AND MAINTENANCE15
OTHER INFORMATION16
TROUBLE SHOOTING17
TECHNICAL SPECIFICATIONS18

INTRODUCTION

Thank you and congratulation on purchasing your new Aspire mobility Scooter. It is designed to provide you with transportation ability indoors and outdoors.

We pride ourselves on providing safe and comfortable products. Our goal is to ensure your complete satisfaction. We sincerely hope you enjoy your Aspire mobility scooter.

Please read and observe all warning and instruction provided in owner's manual before you operate with various convenient function of this scooter. Also, please retain this booklet for future reference.


MANUFACTURED AND IMPORTED BY:

AIDACARE PTY LTD.

Building 3A, 1 Moorebank Avenue, Moorebank NSW Australia 2170 1300 133 120 | www.aspirecare.com.au

IMPORTANT PRECAUTIONS

- Only one person at a time could ride an Aspire Mobility Scooter.
- Maximum load is 181 kgs / 400 pounds.
- Turn key off before getting on or off.
- Always drive carefully and be aware of others using the same area.
- Always use pedestrian crossings wherever possible. Take extreme care when crossing roads.
- Do not drive on slope exceeding 10 degree, and take extreme care when turning on slope.
- Do not use full power when turning to sharp corner.
- Take great care and drive in low speed when backing up, riding downhill or on uneven surface, and climbing curb.
- Please use the lowest speed when driving through the descending road or uneven terrain. If speed is too fast, leave your hand off the handle bar, let the scooter stop. Make sure safety and start again.
- The weight capacity limit at different ramp degree (please refer to following picture).



• The grade climbing degree will be affected by weight capacity, driving speed, and ramp degree, and scooter parameter.

• To prevent any danger from motor defected; please avoid to drive on long ramp or any uneven terrain.

- A slow speed must always be used when ascending, descending or traversing aslope or incline and also on uneven terrain, ramps and soft or loose surfaces, such as gravel or grass.
- To prevent any danger, do not turn around at high speed on ascending, descending ramp.
- Scooter may not operate well in high humidity.
- Do not leave the powered scooter in a rain storm of any kind.
- Do not use the powered scooter in a shower.
- Direct exposure to rain or dampness will cause the scooter to malfunction electrically and mechanically; may cause the powered scooter to prematurely rust.
- Never put scooter in neutral when staying on slopes.
- Follow traffic laws when riding outside.
- When scooter on moving transport vehicles, do not sit or stay on scooter.

ELECTROMAGNETIC INTERFERENCE AND WARNINGS

CAUTION: It is very important that you read this information regarding the possible effects of Electromagnetic Interference on your motorized scooter.

Powered wheelchairs and motorized scooters may be susceptible to electromagnetic interference (EMI), which is interfering electromagnetic energy (EM) emitted from sources such a radio stations, TV stations, amateur radio (HAM) transmitters, two-way radios, and cellular phones. The interference (from radio wave sources) can cause the motorized scooter to release its brakes, move by itself, or move in unintended directions. It can also permanently damage the motorized scooter control system. The intensity of the interfering EM energy can be measured in volts per meter (V/m). Each motorized scooter can resist EMI up to a certain intensity. This is called its "immunity level." The higher the immunity level, the greater the protection. At this time, current technology is capable of achieving at least a 20 V/m immunity level, which would provide useful protection from the more common sources of radiated EMI. The immunity level of this motorized scooter model is not known.

There are a number of sources of relatively intense electromagnetic fields in the everyday environment. Some of these sources are obvious and easy to avoid. Others are not apparent and exposure is unavoidable. However, we believe that by following the warnings listed below, your risk to EMI will be minimized.

The sources of radiated EMI can be broadly classified into three types :

1.Hand-held portable transceivers (transmitters-receivers) with the antenna mounted directly on the transmitting unit. Examples include: citizens band (CB) radios, "walkie talkie," security, fire, and police transceivers, cellular telephones, and other personal communication devices

Some cellular telephones and similar devices transmit signals while they are ON, even when not being used

- 2.Medium-range mobile transceivers, such as those used in police cars, fire trucks, ambulances, and taxis. These usually have the antenna mounted on the outside of the vehicle; and
- 3.Long-range transmitters and transceivers such as commercial broadcast transmitters (radio and TV broadcast antenna towers) and amateur (HAM) radios

Other types of hand-held devices, such as cordless phones, laptop computers, AM/FM radios, TV sets, CD players, and cassette players, and small appliances, such as electric shavers and hair dryers, so far as. we know, are not likely to cause EMI problems to your motorized scooter.

Motorized Scooter Electromagnetic Interference :

Because EM energy rapidly becomes more intense as one moves closer to the transmitting antenna (source), the EM fields from hand-held radio wave sources (transceivers) are of special concern. It is possible to unintentionally bring high levels of EM energy very close to the motorized scooter control system while using these devices. This can affect motorized scooter movement and braking. Therefore, the warnings listed below are recommended to prevent possible interference with the control system of the motorized scooter.

Warnings :

Electromagnetic interference (EMI) from sources such as radio and TV stations, amateur radio (HAM) transmitters, two-way radios, and cellular phones can affect motorized scooters. Following the warnings listed below should reduce the chance of unintended brake release or motorized scooter movement which could result in serious injury.

- Do not operate hand-held transceivers (transmitters-receivers), such as citizens band (CB) radios, or turn ON personal communication devices, such as cellular phones, while the motorized scooter is turned ON;
- 2.Be aware of nearby transmitters, such as radio or TV stations, and try to avoid coming close to them;
- 3.If unintended movement or brake release occurs, turn the motorized scooter OFF as soon as it is safe;
- 4.Be aware that adding accessories or components, or modifying the motorized scooter, may make it more susceptible to EMI; and


There is no easy way to evaluate their effect on the overall immunity of the motorized scooter.

5.Report all incidents of unintended movement or brake release to the distributor listed on the inside front cover of this manual. Note whether there is a source of EMI nearby.

Important Information:

- 1.20 volts per meter (V/m) is a generally achievable and useful immunity level against EMI (as of May 1994). The higher the level, the greater the protection.
- 2. The immunity level of this product is at least 20/Vm.

SAFETY WARNING AND INSTRUCTION LABELS



IDENTIFICATION OF PARTS

Before attempting to drive this scooter on your own, it is important that you familiarize yourself with the controls, and how to operate






Figure 1 - HS-898 Front View







FUNCTION OF PARTS :

MAIN SWITCH (A)

1.By turning the key clockwise to NO to turn on the power

2.By turning the key anticlockwise to OFF to turn off the power.

Speed Dial Knob (B)

The Rabbit means fast and Turtle is slow.

By turning this you control total speed transferred to thumb controls.

Do not set the speed at maximum when operating the scooter. Adjust the speed at maximum when inclining a slope, and adjust to minimum speed when declining a slope.

Horn Button (C)

Pressing down the button to sound buzzer.

Battery Gauge (D)

There is a meter shows batteries capacity status.


When battery gauge's light trends to 🚝, it means sufficient power capacity.


When battery gauge's light trends to E, it means insufficient power capacity.

Check battery's capacity before driving, and charging batteries when at low power.

Hi/Lo (1)

High & Low Speed Switch allows speed to be changed from the normal maximum to half for pavement use. Light means activated Hi/Lo speed.





Panel Control

Pressing down (2) 🔁 to activate left turn signal.
Pressing down (3) 🚺 to activate headlight.
Pressing down (4) \square to activate right turn signal.
Pressing down (5) 🛕 to activate warning light.

%Right-hand/Left hand direction indicator (switches itself off automatically after 30 seconds).



Battery Gauge

Capacity Display	LED	D1	D2	D3	D4	D5	D6	D7	Battery Capacity%
	7	☆	☆	☆	☆	☆	☆	☆	>80
	6	☆	☆	☆	☆	☆	☆		<80
	5	☆	☆	☆	☆	☆			<65
	4	☆	☆	☆	☆				<50
	3	☆	☆	☆					<35
	2	☆	☆						<25
	1	☆							<20
	D1	D4 : A	mber	LED	LED				
Operation Characteristic	LED blocks will only decrease gradually during driving.								
Usage Condition	(1)When turn off power (2)When charging, function will shut off.								

Low Voltage Warning Sound

Operation Characteristic	Warning sound will signal, when battery capacity's lower than 25%, after turn on the power, 2 short "Be" sounds will signal 3 times (Be Be-Be Be-Be Be)
Remark	Can be set as on or off.

Warning Sound (On, Off) Setting

Operation Characteristic	Sound of turn signal, stop, low voltage, reversing can be set as on or off
	1.Ensure power is turned off
	2.Pressing left and right turn signal buttons downwards simultaneously.
	3.Turn on the power.
Description of Volume Setting Steps	4.After 2 seconds, LED (D7) flashes to indicate the setting start, if hold on buttons for 10 seconds without releasing (left and right turn signal buttons) will automatically leave setting mode.
	5.Pressing volume setting button (As Stop (Warning) light button) to access volume setting mode, by pressing button (Ex: Stop (Warning) light button) to adjust volume.
	6.Waiting for 10 seconds without action, or pressing left and right turn signal buttons downwards simultaneously for 2 seconds to return to normal operation mode.

Table of Warning Sound Settings

Function	Volume Button	Volume Adjustment	
Reverse Warning Sound	Head Light Button	Flash	
Turn Signal Warning	Left or Right Turn Signal	■ [] [] [] []] Off	
Sound	Button	🛯 🗶 🛛 🗶 🗍 🖉 🖉 Minimum	
Stop (Warning) Warning Sound	Stop (Warning) Light Button	Low	
Low Voltage Warning		Mid	
Low Voltage Warning Sound	Hi & Lo Button	Loud	
Horn Volume	Horn Button	Except horn volume can not be set as off	

P.S: Pressing volume button (ex:Stop (Warning) to access volume adjustment mode, the LED displays to volume blocks and the volume (ex:Stop (Warning) Sound signal, press (Stop (Warning) Light button again, the volume and LED will decrease gradually, by pressing button once will increase or decrease one block:

Loud >> Mid >> Low >> Minimum >> Off

	Reversing Warning Sound Volume Adjustment
	1.Turn off the power (Key off)
	2.Pressing left and right turn signal buttons downwards simultaneously.
	3.Turn on the power(Key on)
	After 2 seconds, LED (D7) flashes (meaning accessed in volume setting mode),
Examples	Release the buttons within 10 seconds (left and right turn signal buttons).
	4.Pressing Head Light button to access Reverse Warning Sound.
	5.Pressing Head Light button to adjust volume
	6.After setting, waiting for 10 seconds without action, or pressing left and right turn signal buttons downwards simultaneously for 2 seconds to return to normal operation mode.
	Turn Signal Warning Sound=Mid
	Stop (Warning) Light Warning Sound=Mid
Default	Low Voltage Warning Sound=Mid
	Reversing Warning Sound=Mid
	Horn Volume=Loud

Throttle

Pulling right Throttle (E) moves scooter forward. Pulling left Throttle (F) moves the scooter backward. (This can be reversed if required by local dealer.) Releasing both, engages automatic brake. These are also your accelerator. The further you depress them, the faster you go. (Subject to the position of the Rabbit/Turtle control).



Always ensure that the scooter is switched off before getting on or off the scooter and before removing any items of the scooter.

Tiller Angle Adjustment

Pulling the tiller adjustment (G) downwards to adjust tiller's angle and release to lock at a desired comfortable position.


Seat Swivel Adjustment


Setting at a comfortable angle by pulling lever (H) upwards to rotate seat (a) left and right.

Seat fore-aft Adjustment

Setting at a comfortable position by lifting lever (I) upwards to adjust the seat (a).

When driving the scooter, set the seat at foremost position to prevent tip over.






Armrest Adjustment Thumbscrews

- 1.Loosen the two thumbscrews (K) to adjust arm width tighten again to lock in at desired position.
- 2.Adjust the attached screw (J)'s height, to control armrest's angle.

Pull the armrest up when get on or off the scooter. Do not hang heavy parts on the armrests to cause tip over.



Anti-tip wheels (M)


Helps keep scooter from tipping over.

Free-Wheeling Lever (L)

When lever is in N (Neutral) position, scooter can be moved manually without power. When lever is in D (Drive) position, scooter can be driven. Normal position is D.



Freewheel operation is only recommended on flat surfaces, never on gradients. Never leave your scooter on a gradient with its motors disengaged. Always re-engage the motors immediately after pushing the scooter.



CHARGING THE BATTERIES

Batteries must be charged before using the scooter for the first time and should be recharged after each day use. You will need the scooter and the battery charger.

Each country may supply different charger. The charging procedure may be different from below. If you require more details, please contact your authorized dealer.

Be sure the scooter key is in the OFF position. Do not use extension cord to connect the charger. Do not use other brand charger to charge the scooter.

- 1.Insert battery charger cord into the charger connector on the charger output. Refer to above figure for correct Figure 12
- position
- 2.Plug the other end of the battery charger cable into a standard electrical outlet.
- 3. Turn the power on. Normally, The LED (Power) Light will turn on when electric current passes.
- 4. Charging starts. During charging, LED (Charge) will indicate orange light, when it turns to green light, that means well-charged.
- **5.LED** indication
 - LED(POWER) GREEN LIGHT ON : Power On.
 - LED(CHARGE) ORANGE LIGHT ON : Charging / GREEN LIGHT ON : Fully Charged
- 6.Charger Trouble Shooting
 - (A)If LED (POWER) light is off
 - Check the input voltage (115V/230V) is the same as you adjust.
 - If light is still off, please check and repair the battery charger.
 - (B)If LED (CHARGE) light is off
 - · Check to see clips connection is correct.
 - If the battery is fully charged, the LED (CHARGE) light will be off.
 - If light is still off, the battery may be defective.
 - (C)If ORANGE light can turn to GREEN
 - The battery can not be charged. Please check and recover it.
 - (D)If ORANGE light turns to GREEN immediately
 - Check to see the battery is fully charged, if not, The battery may be defective Check and recover it.

The time needed to recharge will vary depending on the depletion of the batteries. Charging for longer than necessary will not harm the batteries. They can not be overcharged.

Keep in mind these rules:

- Fully charge batteries at least once a month, or more if you use scooter regularly. Charge after each trip exceeding 3 kilometers.
- If storing your scooter for some time (1 month or more), make sure that batteries are fully charged, and on returning, charge them again before using scooter.
- Batteries will only give maximum performance after scooter has been used, and batteries have been recharged up to 10 times. A bit like breaking in a new car.

Please be aware that the travelling range of your mobility scooter is impacted by how fast the batteries are discharged. This will depend on many circumstances, such as ambient temperature, condition of the surface of the road, tyre pressure, weight of the driver, driving environment (inclines etc.) and utilisation of your lighting system if fitted. We recommend that you test your local ride with a family member to ensure a safe journey.

HS898

DISASSEMBLING YOUR SCOOTER

Disassembling Basket


By lifting the basket (O) upwards.

Disassembling Seat


By lifting lever (H) and seat (a) upward.

Rear Compartment Cover

Open the Compartment Cover (P) according to arrow's direction.








Disassembling Batteries

Tear off battery straps (Q) and disconnect battery connectors (R) to remove the batteries (b).

Be careful when removing batteries, as batteries are heavy. Do not contact battery terminal **O** with any metals to avoid danger.



ASSEMBLING YOUR SCOOTER

To assemble scooter, you can repeat disassembly directions in reverse. Abbreviated directions are given below. Refer to Figures on pages 12 to locate parts.

- 1.Use the tiller adjustment to move tiller up and out of the way.
- 2.Place front basket.
- 3.Place battery pack in battery compartment.
- 4.Place seat on seat post and lift lever (H) to insert the seat into the seatpost, then swivel the seat to lock at a desired position.

OPERATING YOUR SCOOTER

You could make the following adjustments to increase your comfort when driving.

- adjust seat position.
- adjust armrest width to comfortable position.
- adjust tiller's angle.

1.Before operating your scooter, check the following:

- free-wheeling lever in D
- speed dial knob is at turtle picture.
- 2.Sit on scooter and turn on key, Battery Gauge meter should be indicate at Full. The Self-Diagnostic Warning Light should not be blinking.
- 3.When your hands rest comfortably on handlebars, the Throttle should be within easy

reach. The right lever moves scooter forward, the left one moves it backward. When you release both levers, scooter will stop.

This scooter has automatic braking system. Release the Throttle and brake will stop scooter.

- 4. Steer scooter by turning tiller toward the way you want to go.
- 5.Practice driving where there are is No obstacle. Start at the slowest speed and drive forward and backward; make some turns. As you get more comfortable, you can increase speed by turning speed dial toward picture of rabbit.
- 6.If Battery Gauge indicates 🛅, you should plan to recharge batteries very soon.
- 7. When you are finished riding, turn off the key before getting off.
- 8.If you are finished riding for the day, immediately recharge batteries.

See CHARGING THE BATTERIES, page 11.

Keep in mind these rules \colon

- Release Throttle and allow scooter to stop completely before changing from forward to reverse, or reverse to forward.
- When turning to corner, swing front wheels widely, so back wheels will turn more tightly.
- Use scooter only where it is safe to walk.
- Drive in low speeds when reversing, riding dowhill, over rainp or curb, or on uneven sukpace, downhills, ramps, curbs, or uneven surfaces.

Other Operating Information:

Hill climbing: You may need to use a higher speed. Turn to lower speed before going downhill.

- **Down slopes:** proceed with downward slopes slowly, and set speed control in proximity of turtle. The closer speed control is set toward turtle, the slower scooter travels. However, this scooter will not self accelerate down hills due to the automatic braking, taking effect should you attempt to drive too fast.
- **Kerb climbing:** Approach slowly from right angles to curb. A slight angle is permissible with a 4-wheel scooter. Do not attempt greater than a 7cm curb. If Self-Diagnostic Warning Light starts to blink, identify problem from

chart on page 16 and take action.

If the scooter breaks down and must be moved, get off scooter, engage Free-Wheeling Lever to N, push scooter slowly to a safe location, and push lever back to D.

In unlikely event of a panel display error, you need to re-set the display system by cycling the on/off main switch. The display circuitry is independent of the motor control system. A display console error does not affect scooter speed control.

CARE AND MAINTENANCE

DAILY CHECK

Please always check your scooter before you start your every journey.

Check point	Inspection	What to do if the inspection is failed
N-D lever	Check for correct function	Contact your dealer.
Horn	Check for correct function	Contact your dealer.
Throttles	Pull the wigwag to test the scooter movement	Contact your dealer.
Electro- mechanical brake and Emergency hand brake	Pull the wigwag a little bit and release it to test if brake works. If your scooter comes with emergency handbrake, please check it as well.	Contact your dealer.
Battery Gauge	Check if the battery gauge is displayed and whether it is at low power.	 Contact your dealer if battery gauge is not working. Recharge the battery immediately if low
Rear mirror (s)	Check if the parts are clean and firmly tighten to the scooter and do not wobble.	 Clean up the dirt by damp cloth. Tighten the screw or clamping stem that holds the mirror(s)
Lighting	Check if all lights, such as head lights, rear lights, and turn signal are functioning correctly.	Contact your dealer.

WEEKLY CHECK

Check point	Inspection	What to do if the inspection is failed
Speed Dial Knob	Check for correct function	Contact your dealer.
Armrests	Check if the parts are clean and firmly tighten to the scooter and do not wobble. Tighten the screw knob that holds the armrest.	Contact your dealer.
Wheels/Tires	 Inflate the tire to the correct pressure and check that 1.Drive wheels rotate without wobbling. 2.Tire tread depth is not less than 0.5mm. 3.No foreign objects in tires. 	Contact your dealer.
Motor	Check for correct function	Contact your dealer.
Battery Charger	Check if the charger is functioning correctly and the batteries are charged.	Contact your dealer.

MONTHLY CHECK

Check point	Inspection	What to do if the inspection is failed
Seat / Upholstery	Check for movement and if it's worn	Contact your dealer.
Electronics	Check if all the battery cables and connectors are firmly tighten to the scooter	Contact your dealer.

CLEANING YOUR SCOOTER

- Do not use any abrasive or scouring liquids for cleaning. Only use a damp cloth and gentle detergent.
- Do not use hose pipe or splash water directly onto the scooter as this may cause damage to electronics.

MAINTENANCE

- User should inspect the scooter regularly to keep scooter in good running order.
- Check if the electrical cable connectors are fully connected.
- All maintenance and repair of scooter should be done by an authorized dealer.

Seat Upholstery :

Only use damp cloth and a little soap to wipe the seat. Do not use abrasive cleaners as this will damage the seat.

Storage :

- Please store the scooter in a dry location. If store the scooter in long time, please disconnect the battery terminals.
- Do not store your scooter where it will be exposed to source of direct heat, damp, oil. acid, alkaline, or where Ozone could be possibly generated. All of the above will minimize scooter / tire cycle and shorten its lifetime.

Remark :

- Obstacle height approx. 7cm (Figure.17)
- Feasible groove width approx. 20cm, scooter must go straight line to across the groove. (Figure.18) Suggest, not to across the obstacle over 5cm.
- When driving scooter on ramp, adjust body center of gravity to keep scooter more safety.







General driving posture



On ramp forward your body will let scooter more safety.

OTHER INFORMATION

RECYCLING & DISPOSAL

- The equipment wrapping is potentially recyclable.
- The metal parts are used for scrap metal recycling. The plastic parts are used for plastic recycling.
- Electric components and printed circuit boards are disposed of as electronic scrap.
- Exhausted or damaged batteries can be returned to your medical equipment supplier.
- Disposal must be carried out in accordance with the respective national legal provisions.
- Ask your city or district council for details of the local waste management companies.

SERVICE LIFE

We estimate a service life of five years for this product, provided it is used in strict accordance with the intended use as set out in this document and all maintenance and service requirements are met. The estimated service can be exceeded if the product is carefully used and properly maintained, and provided technical and scientific advances do not result in technical limitations. The service life can also be considerably reduced by extreme or incorrect usage. The fact that we estimate a service life for this product does not constitute an additional warranty.

TROUBLE SHOOTING

Here are some suggestions about solving problems you may have with your scooter. There is a self-diagnostic warning light on the control panel. To check the self-diagnostic warning light, turn on the key and count the number of blinks on the warning light.

SCOOTER WON'T MOVE WHEN KEY IS TURNED ON

Check point	Solution	
Check if the power is off	Turn the power on.	
Check if the N-D lever is in Neutral position	Switch to D (drive) position. Turn off the power and turn on again.	
Check if the battery power is enough. (Battery gauge is under 25%)	Recharge the battery and then retest.	
Check if the charger power cord is still plugged in scooter	Unplug the charger power cord.	
1.Check Power Reserve Indicator on control panel, it should be lighting in green, yellow, and red zones.		

2. Check Self-Diagnostic Warning Light, it should be steady. If it is flashing, see chart below for problem identification.

3. Check all electrical connections to be sure they are tight.

4.If none of above correct problem, contact your authorized dealer.

ERROR CODE

Flash	Description	Initial check points		
1	Battery Low	The batteries are running low. • Recharge the batteries.		
2	Low Battery Fault	The batteries have run out of charge.Recharge the batteries.Check the battery and associated connections and wiring.		
	The Low Battery Fault Flash code, described above, is a requirement of various safety standards. The scooter will output a visible and audible low battery warning if the battery voltage drops below 90% of its cut-off voltage. The warning will be two short flashes, and will take priority over all other flash codes in the system.			
3	High Battery Fault	 Battery voltage is too high. This may occur if overcharged and/or traveling down a long slope. If traveling down a slope, reduce your speed to minimize the amount of regenerative charging. 		
4	Current limit time-out or controller overheat	 The motor has been exceeding its maximum current rating for too long. The scooter may have stalled. Turn the controller off, leave for a few minutes and turn back on again. The motor may be faulty. Check the motor and associated connections and wiring. 		
5	Park Brake	Either a park brake release switch is active or the park brake is faulty. • Check the park brake and associated connections and wiring. • Ensure any associated switches are in their correct positions.		
6	Drive Inhibit	 Either a stop function is active or charger inhibits or OONAPU condition has occurred. Release the stop condition (seat raised etc.) Disconnect the battery charger. Ensure the throttle is in Neutral when turning the controller on. The throttle may require re-calibration. 		
7	Speed Pot	The throttle, speed limit pot. SRW or their associated wiring may be faulty. • Check the throttle and speed pot and associated connections and wiring.		
8	Motor Voltage	The motor or its associated wiring is faulty. • Check the motor and associated connections and wiring.		
9	Other Error	The controller may have an internal fault. • Check all connections and wiring.		

OTHER PROBLEMS

• TIRE : Low tire pressure: pump up tires to 30~35 psi.

• Charger : During charging, light on charger doesn't change to green. Contact your authorized dealer.

TECHNICAL SPECIFICATIONS

Overall Length	1440 mm	
Overall Width	720 mm	
Overall Height	1360 mm	
Front Wheels	330 mm	
Rear Wheels	330 mm	
Weight W/ Batteries	145 kg	
Max. Speed	10 kmph	
Weight Capacity	180 kg	
Ground Clearance	100 mm	
Grade Climbable	10 degree	
Curb Climbable	80 mm	
Turning Radius	1600 mm	
Suspension	Front & Rear	
Brake	Hand Brake & Electro-Mechanical	
Seat Type	Swivel Captain Seat With Seat Sliding Mechanism & Back Angle Adjustment	
Seat Width	520 mm	
Motor Size	850W, 3450 r.p.m.	
Battery Size	(2) . 75Ah	
Weight of Battery	55 kg	
Travel Range	40 km	
Battery Charger	8A Off Board	
Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Knob	

*Subject to change without notice. (Rev. 5, 2020/06/29)

SERVICE LOG BOOK

YEAR 1	Date:
Service Type:	
Condition Report:	
Action taken:	
Inspected By:	
Signature:	
YEAR 2	Date:
Service Type:	
Condition Report:	
Action taken:	
Signature:	
YEAR 3	Deter
Commiss Transm	Date:
Service Type:	
Action taken:	
Charles 1	
Signature:	
YEAR 4	
	Date:
Service Type:	
Condition Report:	
Action taken:	
Inspected By:	
Signature:	
YEAR 5	
	Date:
Service Type:	
Condition Report:	
Inspected By:	
Signature:	

YEAR 6	Date:
Service Type:	
Condition Report:	
Action taken:	
Inspected By:	
Signature:	
YEAR 7	Data
Service Type:	Date:
Condition Report:	
Action taken:	
Inspected By	
Signature:	
YEAR 8	Date:
Service Type:	
Condition Report:	
Action taken:	
Inspected By:	
Signature:	
YEAR 9	Date:
Service Type:	
Condition Report:	
Action taken:	
Signature:	
VEAD 10	
YEAR 10	Date:
Service Type: Condition Report:	
Action taken:	
Inspected By:	
Signature:	



MANUFACTURED AND IMPORTED BY:

AIDACARE PTY LTD. Building 3A, 1 Moorebank Avenue, Moorebank NSW Australia 2170 1300 133 120 | www.aspirecare.com.au