

# **COMPACT**

# **USER MANUAL**

**Exclusively Distributed by:** 



© 2023 Life Assistant LTD

Manual identification: LAPLAB.000.040-R01, March. 2023

© 2023 *Life Assistant* Ltd, all rights reserved. No part of this manual may be reproduced in any form without the written permission of *Life Assistant* Ltd.

**Life Assistant** Ltd. reserves the right to change or improve its products and accompanying technical literature without specific notice to changes or improvements.



# Manufacture

Life Assistant LTD,

8 Harod St.Afula, 1871208, ISRAEL email info@sitnstand.com web www.lifeassistant.co.il







#### **Authorized Representative in the EU:**

Advena Ltd.

Malta Office Tower Business Centre, 2nd Flr, Tower Street, Swatar, BKR 4013, Malta UK Office Pure Offices, Plato Close, Tachbrook Park, Warwick, CV34 6WE, UK





# TABLE OF CONTENTS

INTRODUCTION	E
About This Manual	
Life Assistant Sitnstand Product Overview	
Contents	8
Special messages: Warnings, Cautions and Notes	8
Conditions for Use	
SAFETY	10
General Warnings and Cautions	1
LABELS AND SYMBOLS	12
USING THE SITNSTAND	14
Setup Instructions	15
Following the Appropriate Seat Selection	17
Controller indicators modes	19
Sitting Down	20
Standing Up	20
Replacing the Battery	2
Charging the Battery	2
MAINTENANCE	22
TROUBLESHOOTING	24
Basic Troubleshooting	25
To Obtain Information	27
OPERATING SPECIFICATIONS	
Environmental	29
1.2. Electrical	
1.3. Physical and Operational	32
LIMITED WARRANTEE	32
EMC STATEMENT	34
Warnings	
Electromagnetic Emissions	36





#### **ABOUT THIS MANUAL**


This manual provides the information necessary to operate **Life Assistant SitnStand Compact** in a safe and efficient manner. Please review and understand this manual before operating the system. If any part of this manual is unclear, contact your supplier or **Life Assistant** Customer Support for clarification.

#### **Life Assistant Sitnstand Product Overview**

**The SitnStand Compact** is a lightweight, portable inflatable system that is safe, simple and easy to operate. It comprises an electronic control system, multiple air valves, multiple air sealed layers, and a portable air pump that supplies the required pressure and air capacity to assist standing up and sitting down with minimum/no help required.







#### **Contents**

- 1 **SitnStand Compact** unit
- 1 SitnStand Compact compatible case
- 1 Li-Ion 14.4VDC rechargeable battery
- 1 Wall battery charger with 1m cable and connector
- · 1 User manual

## **Warnings, Cautions and Notes**

Three types of special messages appear in this User Manual:



**WARNING** indicates the possibility of injury to the patient or operator.



**CAUTION** indicates a condition that may lead to equipment malfunction.



**INFORMATION** provides other important information.



#### **Conditions for Use**

#### **INTENDED USE**

The SitnStand Compact is a powered inflatable cushion intended for assisting patients suffering from mobility difficulties in rising from a seated position to standing, as well as in sitting down from a standing position. The SitnStand Compact is intended to be used with any decent seating area.

The SitnStand Compact is intended to be used in the homecare environment while performing daily activities, as well as in medical care facilities.

#### **INTENDED USERS**

The SitnStand Compact is intended to be used by elderly people and/or people with compromised mobility or difficulties.

The SitnStand Compact is intended to be used by people with chronic or acute mobility difficulty.

#### **CONTRAINDICATIONS**

- Do not use **The SitnStand Compact** while under the influence of narcotics or alcohol.
- Do not use The SitnStand Compact while holding a child or infant.
- Do not use **The SitnStand Compact** inside a vehicle.
- Do not use The SitnStand Compact without clothing.
- In the case of Cognitive impairment, use **The SitnStand Compact** only by a caregiver.





## **General Warnings and Cautions**



WARNING: DO NOT USE BEFORE READING AND UNDERSTANDING THIS MANUAL.



**WARNING:** To avoid the potential for electrical shock, use only the designated charger provided with the device. Using any other charger may cause an electrical and/or electromagnetic hazard.



**WARNING:** To avoid the potential for electric shock, the power supply must be disconnected from the power mains prior to cleaning.



**WARNING:** To avoid the potential for electric shock, contact your seller or *Life Assistant Ltd*. when a malfunction occurs.



**WARNING:** Do not inflate under the bottom of a table or other item that interferes with the cushion inflation.



**INFORMATION:** The only user-serviceable components in the system are the battery and cushion jacket.



**CAUTION:** Changes or modifications not expressly approved by *Life Assistant Ltd* can affect the safety and effectiveness of the system, and will invalidate the system's warranty.



**CAUTION:** Do not operate with damaged cords or plugs. If damaged, please contact the company.



**CAUTION**: At the end of its useful life, the system and battery must be disposed of in accordance with local law and/or code concerning electrical and electronic equipment.





# LABELS AND SYMBOLS

SYMBOL	DESCRIPTION
MD	Medical Dvice
SN	Serial Number
REF	Catalog Number
$\bigcap$ i	User Manual
	Keep Dry
	Manufactured by
	Not for general waste Special Requirements for Waste of Electrical and Electronic Equipment (WEEE Directive)
$\triangle$	Caution
i	Information
CE	CE Mark
EC REP	Authorized Representative





# THE SITNSTAND COMPACT

#### **Setup Instructions**

Setting up the SitnStand Compact is performed by the user.

The setting up process consists of three steps:

- A. Charging the battery for the first time.
- B. Placing *the SitnStand Compact* on an appropriate seat.
- C. Connecting the battery.


#### A – CHARGING THE BATTERY FOR THE FIRST TIME-FIGURE 2

Charge the battery for at least 12 hours before use:

- 1. Remove the battery from the box.
- 2. Plug in the battery charger into an electric wall socket.
- 3. Connect the battery charger connector to the battery charging socket.
- 4. When the green light is on, battery is fully charged.
- 5. When red light is on, battery is not fully charged.



**NOTE:** To prevent falling over any tangled cable, place the battery close to the charger.





14



# B – PLACING THE SITNSTAND COMPACT ON A SUITABLE SEAT – FIGURE 3

Selection of an Appropriate Seat

When selecting the seat on which **the SitnStand Compact** will be fitted, it is advisable to select a seat that is steady, stable and cannot be easily moved. An armchair is preferable.


An appropriate seat size for **the SitnStand Compact** should be at least the minimum size of **the SitnStand** Compact seating area in its deflated state, with at least one arm rest.



## Following the Appropriate Seat Selection

Place the seating area of **the SitnStand** Compact directly onto the chair seat.

Place the strip on to the backrest of the chair sit as shown in the picture bellow.










#### C - PLACING THE BATTERY IN ITS PLACE - FIGURE 6

Slide the charged battery into its place within the air supply unit, as illustrated in **Figure 6**.

**The SitnStand Compact** is now in standby mode, and is ready for use.




# Sit Stand

18

#### **Controller Indicator Modes - Figure 7**

- 1. **Standby Mode** the Green light flashes slowly (about one flash per second) *The SitnStand Compact* is now ready for use and is waiting for your action.
- 2. **Working Mode –** the Green light flashes quickly (about 5 times a second) *The SitnStand Compact* is operating (while in this mode, the buttons will not respond).
- 3. **Charge Battery** the Orange light is on while the Standby/ Working mode green light is also flashing as described in 1 or 2 above indicating that the battery needs recharging.
- 4. **Charge Battery** the Orange light is on and Standby/Working mode Green light is off battery is depleted. In this case, the SitnStand Compact cushion can only be deflated.
- 5. **Error –** the Red and Green lights flash quickly **The SitnStand Compact** has identified a system fault See troubleshooting chapter.





#### **Sitting Down**

To sit down – inflate the cushion using the controller to 1,2,3 or 4 layers (whichever is most comfortable for you), then sit down on the pillow facing forward.

When you feel steady, start deflating the cushion by pressing on the Down Arrow button on the controller.

## **Standing Up**

To stand up – First, disconnect your upper-back from the back of the chair by bending slightly forward, and only then, start inflating the cushion one layer at a time by pressing the Up-Arrow button on the controller.

**The SitnStand Compact** has 4 layers: however, Inflation can be ended at any time when the level of elevation is sufficient.

Do not attempt to stand up before the layer is totally inflated and there is no sound from the air pump.

**The SitnStand Compact** may remain inflated until the next time it is used for sitting down.

### Replacing the Battery - Figure 8


When the Charge Battery indicator turns Orange, battery should be charged.

Open the zipper covering the air supply unit and remove battery by sliding it towards you as illustrated in **Figure 8**.

The battery is connected by a magnet to the air supply unit.




**NOTE:** Don't remove the battery while The **SitnStand compact** is at "Working Mode" (see controller indicator modes)



## **Charging the Battery-Figure 9**

- 1. Plug the battery charger into an electric wall socket.
- 2. Connect the battery charger connector to the battery charging socket.
- 3. When the Red light on the charger changes to Green, the battery is fully charged.
- 4. When the Red light is on, the battery is not fully charged.







# **MAINTENANCE**

22

In order to maintain *The SitnStand Compact* system, please note the following:

- 1. Keep **The SitnStand Compact** away from dust or dirt
- 2. Use a soft dry brush to clean the fabric surfaces (e.g. seat area).
- 3. For cleaning **The SitnStand Compact** use A damp cloth
- 4. Do not iron the cloth jacket
- 5. Use a damp cloth to clean the plastic surfaces (e.g. air supply unit)
- 6.Do not use solvents or detergents to clean any part
- 7. Do not allow liquids near the system openings
- 8. Keep dry

User Manual | Maintenance



**WARNING:** To avoid potential for electrical shock, the power supply must be disconnected from the power mains prior to cleaning.



**WARNING:** To avoid the potential for electric shock, contact your seller or *Life Assistant Ltd*. when a malfunction occurs.

- (i)
  - **INFORMATION:** The only user-serviceable components in the system are the battery and cushion jacket.
- (i)
  - **INFORMATION:** No portion of the system should be immersed in water or other fluids
- (i)

**INFORMATION:** Never spray or otherwise apply any cleaner fluid directly into the openings or seams of the system. Always apply the cleaner fluid to the cloth instead.





# **TROUBLESHOOTING**

Before addressing the troubleshooting table, please check and confirm the following:

- 1. The battery is attached to *the SitnStand compact*.
- 2. The battery is charged. (low battery indicator is off, standby/ready indicator is on)

## **Basic Troubleshooting**

When addressing a malfunction, as an initial action please follow STEP 1 described at the table below. In case of unsolved malfunction follow STEP 2.

MALFUNCTION	STEP1	STEP 2
Battery not charging	Charge the battery	Contact your
No light on the controller		seller or <b>Life</b> <b>Assistant Ltd</b> .
Product is not working when pressing the controller		
Does not deflate		
Does not inflate		
Air leaking	Contact your seller or <b>Life Assistant Ltd.</b>	
Red triangle light and the Green light on the controller are flashing	Press the two remote control buttons (up and down arrows) simultaneously for 10 seconds to reset system	Contact your seller or <b>Life Assistant Ltd.</b>



24



## **To Obtain Information:**

Life Assistant LTD.

26

8 Harod St. Afula, 1871208, ISRAEL

email: info@sitnstand.com

web: www.sitnstand.com





# **OPERATING SPECIFICATIONS**



**NOTE:** Unless otherwise indicated, all specifications are subject to change without notice. Specifications and test methods will be made available upon request.

### 1.1. Environmental

User Manual | Operating Specifications

#### 1.1.1. TRANSPORT AND STORAGE ENVIRONMENT

-25°C to +70°C (-13°F to 158°F) **Temperature range:** 

Without relative humidity – 93% **Relative humidity range:** 

#### 1.1.2. OPERATION ENVIRONMENTAL CONDITIONS

+5°C to +40°C (41°F to 104°F) **Temperature range:** Relative humidity range: 15% to 93%, with no condensation 700 – 1060 hPa (10.15 - 15.37 psi) **Atmospheric pressure: Transport altitude:** 0-2000 meters (0 - 6562 Feet)



28



## 1.2. Electrical

Lithium-Ion Rechargeable, 14.4 V, 2.6 Ah
DC 14.4V @ 10A
100 to 240VAC
16VDC, 0.5A (8W)
50 – 60Hz
~500 charging cycles / one year under warranty that the battery should be over 80% of its capacity



**CAUTION:** At the end of its useful life, the system and battery must be disposed of in accordance with local law and/or code concerning electrical and electronic equipment. Do not discard to standard trash bin.

## 1.3. Physical and Operational

Seat footprint – Deflated:	40 x 53 x 1 (cm)   (16 x 21 x 0.4 (ln))
Seat footprint – Inflated:	40 x 48 x 20 (cm)   16 x 19 x 8 (In)
Carrying package dimensions:	43 x 28 x 15 (cm)   17 x 11 x 6 (ln)
Product Weight:	3 kg 6.6 lb
Max supported	<b>USA</b> 180 Kg /400 Lbs
weight:	<b>Non USA</b> 121 Kg / 19 st
Nameiral Oneration	~25 sec
Nominal Operation time (inflate/deflate):	~25 Sec
•	Up / down pushbutton ↑↓
time (inflate/deflate):	
time (inflate/deflate):	Up/down pushbutton ↑↓
time (inflate/deflate):	Up / down pushbutton ↑↓  Green LED: system ready
time (inflate/deflate):	Up / down pushbutton ↑↓  Green LED: system ready  Yellow LED: low bat.
time (inflate/deflate):  Control mechanism:	Up / down pushbutton ↑↓  Green LED: system ready  Yellow LED: low bat.  Red LED: Error





# LIMITED **WARRANTY**

32

Life Assistant Ltd. guarantees this product for a period of two years from the date of purchase.

User Manual | Limited Warranty

The guarantee covers defects in materials and workmanship, comprising the system, assemblies, components and labor, but excluding battery and textile/fabric.

Life Assistant Ltd. guarantees the battery and the textile/fabric for one year from the date of purchase. The guarantee covers defects in materials and workmanship comprising the system, assemblies, components and labor.

The guarantee is not valid in the case of abnormal wear and tear, willful damage, neglect, misuse or unauthorized repair or alteration of the product, or any of its parts or accessories.

The following statements do not affect the consumer's statutory rights. This product is sold on the understanding that, in the event of any defect of manufacture or material appearing within two years of the date of receipt of this product, the defect will be amended free of charge on condition that evidence of date of ownership is provided. Such evidence may be an invoice, delivery note, etc. and, evidence that the defect is not attributable to accidental damage either in transit or otherwise, misuse or unauthorized repair or alteration.

Product Serial Number:
Battery Serial Number:
Location of Purchase:
Date of Purchase:
Delivery Note/Invoice No.: ————————————————————————————————————





# EMC STATEMENT

#### Warnings



**WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided as it may result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.



**WARNING:** Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment, may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment, resulting in improper operation.



**WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of *the SitnStand*, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment may result.

## **Electromagnetic Emissions**

- **SitnStand** is intended for use in the electromagnetic environment specified in the following tables. This is not a life-sustaining device.
- The user and/or installer of the unit must ensure that it is used in such an environment.





DECLARATION - ELECTROMAGNETIC EMISSIONS	CTROMAGNET	IC EMISSIONS
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group1 Class B	<b>The SitnStand compact</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference nearby electronic equipment.
Harmonic emissions IEC 61000-3-2	N/A	<b>The SitnStand compact</b> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings
Voltage Fluctuations And Flicker IEC 61000-3-3:2013	Complies	used for domestic purposes





User Manual







#### **Exclusively Distributed by:**



© 2023 **Life Assistant LTD**, all rights reserved. No part of this manual may be reproduced in any form without the written permission of **Life Assistant Ltd.**