



# Instruction manual

## **W-NCB Wireless Nurse Call Button**

(inc wall mount, wrist strap and safety lanyard)



UK Version

The **Alerta Wireless Nurse Call Button** is for use with Alerta Wall Point Receiver. It gives the patient / user the assurance that help is at their fingertips. It uses radio frequency technology to connect to Alerta Wall Point Receiver with a wireless range of up to 95 meters. Once paired with Alerta Wall Point Receiver it will only transmit to that Wall Point Receiver and will not interfere with other Alerta Wall Point Receivers in use nearby.

The **Alerta Wireless Nurse Call Button** is supplied with a wall mountable holder that can be mounted to a wall or bedside unit for holding the **Wireless Nurse Call Button** meaning it is always in one place if the user is bed bound, also making it easy for those affected by dementia to remember its location.

The **Alerta Wireless Nurse Call Button** is also supplied with a safety lanyard for use over the neck, and a wrist strap enabling it to be wrist worn. The button can be mounted in the wall holder with the wrist strap still attached.

**IMPORTANT:** It is your responsibility to assess the suitability of the patient for whether or not the lanyard should be used. Potential danger of asphyxiation or restriction to blood flow may occur.

## Set-up Instructions

Prior to use the **Alerta Wireless Nurse Call Button** must be paired to the Alerta Wall Point Receiver. It is supplied with a battery.

1. Slightly unscrew the battery cover on the back with a screwdriver.
2. Pull the battery pull tag to remove.
3. Secure the battery cover firmly into place and tighten the screw before operation.



Figure 1

## Replacing the Battery

1. Unscrew the battery cover on the back with a screwdriver and remove the cover. Figure 2
2. Replace with a new CR2450 battery, taking care to seat the battery correctly.
3. Replace the battery cover, ensure firmly in place, and tighten the screw before operation.



Figure 2

## Low Battery

The lower LED will flash slowly in red to indicate battery level is low, and trigger the alarm on the Alerta Wall Point Receiver and also nurse call system if connected. Please replace a new CR2450 as soon as possible.

## Pairing

The **Alerta Wireless Nurse Call Button** must be paired to the Alerta Wall Point Receiver. Please refer to the manual for the W-WPR Wall Point Receiver and see 'Pairing' section of the manual.

Once paired, your **Alerta Wireless Nurse Call Button** is now ready for use.

To activate, press and release the '+' button in the middle of the device once. The outer LED will flash yellow and trigger the alarm on the Alerta Wall Point Receiver and also nurse call system if connected. Please test prior to use with patient / user.

|                                                                                                                                                           |                                                                                     |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <b>CAUTION</b>                                                                                                                                            |  |
| Improper installation of the battery or the installation of the incorrect battery type can result in damage to the contacts and cause improper operation. |                                                                                     |



EU 2017/745



**CS Lifesciences Europe Ltd**  
The Black Church  
St. Mary's Place  
Dublin 7  
D07 P4AX  
Ireland

[eurep@cslifesciences.com](mailto:eurep@cslifesciences.com)

Rev: 2024.04

# Making **cutting-edge** medical equipment **affordable**

Distributed by  
**aidacare**<sup>™</sup>  
HEALTHCARE EQUIPMENT  
1300 133 120 | [aidacare.com.au](http://aidacare.com.au)

