

Bracing & Supports - Neck & Shoulder

I ndis For NDIS Item Code, refer to page: 500



Cosmac Soft Foam Collar

- Washable polyurethane foam covered with stockinette
- Velcro fastening
- Sizing: (1) measure depth (cm) from chin to clavicle (2) measure neck circumference for length



Philadelphia Collar

- Anterior piece cradles mandible to enhance rotation control
- Made of light, strong foam, easy to trim
- Bilateral velcro fastening
- Sizing: (1) measure depth (cm) from chin to clavicle (2) measure neck circumference (S, M or L)

Code	Size	Depth (cm)	Length (cm)	DVA:
FBS301201	Small	7.5	28 - 35	AR18
FBS301202	Small/Long	7.5	35 - 43	AR18
FBS301203	Medium	8.0	28 - 35	AR18
FBS301204	Medium/Long	8.0	35 - 43	AR18
FBS301205	Large	10.0	35 - 43	AR18
FBS301206	XLarge	11.0	38 - 46	AR18

Depth (cm)	Small (25 - 30cm)	Medium (33 - 40cm)	Large (40 - 46cm)	DVA:
6	FBS301207	FBS301207AM	FBS301207AL	AR18
8	FBS301207BS	FBS301207BM	FBS301207BL	AR18
11	FBS301207CS	FBS301207CM	FBS301207CL	AR18
13	FBS301207DS	FBS301207DM	FBS301207DL	AR18



Exo 101 Firm Collar

- Designed for symptomatic relief of neck pain or injury
- Made of polyethylene plastic, with vinyl-covered foam padding
- Height adjustable from 12.5 to 14.5cm
- Air holes for ventilation and hook-and-loop closure



Oppo 4190 Rigid Adj Collar w/ chin

- Completely height-adjustable to provide desired flexion or extension positioning
- Vinyl-covered foam padding for patient's comfort and easy cleaning
- Ventilated to allow air circulation
- Hook and loop closure for quick, easy application and removal

Code	Size	Neck (cm)	DVA:
FBS301208	Small	31 - 37	AR18
FBS301208M	Medium	35 - 41	AR18
FBS301208L	Large	40 - 45	AR18

Code	Size	Neck (cm)	DVA:
FBS301210	Small	34.3 - 40	AR18
FBS301210M	Medium	40 - 46.4	AR18
FBS301210L	Large	46.4 - 52.7	AR18

Note: Measure around the circumference of the neck.

Note: Measure around the circumference of the neck.

DVA Contact Details

P 1300 888 052 | F 1300 787 052 | E dva@aidacare.com.au

I ndis For NDIS Item Code, refer to page: 500



Aspen Collar - Adult

- Designed to optimise support and comfort, the key components for better patient outcomes
- All contact surfaces of the collar are cushioned with cotton-lined, breathable foam padding
- Provides substantial motion restriction without producing painful pressure points that can lead to skin breakdown or poor patient compliance

Code			Size	Top of Shoulders to Chin (cm)	Neck (cm)	DVA:
Collar (1 set of pads)	Collar Set (2 sets of pads)	Replacement Pads				
FBS301211	FBS301215	FBS301219	Adult Short	2	33 - 53	AR18
FBS301212	FBS301216	FBS301219AR	Adult Regular	2 - 3	33 - 53	AR18
FBS301213	FBS301217	FBS301219AT	Adult Tall	3 - 5	33 - 53	AR18
FBS301214	FBS301218	FBS301219AXT	Adult X/Tall	5+	33 - 53	AR18

Note: Measure (1) top of shoulders to chin and (2) around the circumference of the neck.



Miami-J Collar - Adult

- A phenotype-driven sizing system - with unique fronts and backs that are bioengineered to work together for superior immobilisation - makes size selection simple, while drastically reducing the margin of error in fitting a patient
- Adjustable mandibular and occipital support enables customisation to disperse pressure and avoid skin breakdown
- Sorbatex™ pads that are specifically engineered to be anti-decubitus - protecting skin integrity during extended wear - and are removable for easy cleaning
- Sternal pad disperses pressure and is removable for comfort and improved compliance when supine and/or swallowing
- Patient assist strap allows for flexibility and comfort without sacrificing immobilisation

	Code			Size	Sizing Guidelines	DVA:
	Collar (1 set of pads)	Collar Set (2 sets of pads)	Replacement Pads			
Standard Sizes: > 70% of Population	FBS301221	FBS301226	FBS301232	Small	• Standard for mature men • Short-necked women	AR18
	FBS301222	FBS301227	FBS301232M	Medium	• Standard for mature women • Thin, mature men	AR18
	FBS301223	FBS301228	FBS301232L	Large	• Long, tall neck • 'Swan' neck • Young women • Adolescents	AR18
Specialty Sizes: < 30% of Population	FBS301224	FBS301229	FBS301233	Stout	• Very large neck circumference • Obese • 'No neck' • Massive shoulders	AR18
	FBS301225	FBS301230	FBS301232SS	Super, Short	• Kyphotic (chin-on-chest) neck • Geriatric • Osteoporitic and/or ankylosing spondylitic	AR18
	FBS301220	FBS301231	FBS301234	X-Small	• Very short, thin neck principally female • More prevalent in Asian population • Rare	AR18



BodyAssist N71 Thermal Shoulder Brace

- Provides soft tissue support & compression
- Increases elasticity of muscles, ligaments & tendons
- Assists in prevention of injury or re-injury
- Elastic support strap fits comfortably under the arm pit & assists to generate better than average compression & support to the shoulder
- Can be worn while playing sports
- Multi-directional compression assists to reduce soft tissue swelling
- Promotes greater blood circulation, provides a more speedy recovery
- Contoured panels fit most body shapes
- Completely washable & hygienic
- Stabiliser strap is optional

Code - Shoulder Brace		Size	(cm)	DVA:	Code - Stabiliser Strap	Size	To Fit	(cm)	DVA:
Left	Right								
FBS281398	FBS281403	XS	80 - 90	AR03	FBS285281	Regular	To Fit Brace XS to M	80 - 90	AR02
FBS281399	FBS281405	S	90 - 100	AR03	FBS285282	Large	To Fit Brace M to XL	90 - 100	AR02
FBS281400	FBS281406	M	100 - 110	AR03					
FBS281401	FBS281407	L	110 - 120	AR03					
FBS281402	FBS281408	XL	120 - 130	AR03					

Note: Measure circumference around chest, around outside of arm on affected side (refer illustration).



Oppo 1072 Shoulder Support

- Provides compression and support to shoulder and aids healing
- Therapeutic warmth for the shoulder and upper arm
- Offers support and protection
- Minimises potential for future injuries
- Fits left or right

Code	Size	Chest (cm)	DVA:
FBS301236	Small	83.2 - 93.3	AR03
FBS301236M	Medium	93.3 - 103.5	AR03
FBS301236L	Large	103.5 - 113.7	AR03
FBS301236XL	X-Large	113.7 - 126.4	AR03

Note: Measure around the chest (cm).



Oppo 2172 Reinforced Shoulder Support

- Two reinforced elastic straps on upper arm provide extra support for the shoulder
- Allows the patient gradually return to full range of motion
- Fits left or right

Code	Size	Chest cm	DVA:
FBS301237	Small	81.6 - 91.8	AR03
FBS301237M	Medium	91.8 - 101.9	AR03
FBS301237L	Large	101.9 - 114.6	AR03
FBS301237XL	X-Large	114.6 - 127.3	AR03

Note: Measure around the chest (cm).



Oppo 2075 Clavicle Brace

- Excellent aid in correcting poor posture and shoulder slump
- Adjustable hook and loop closure facilitates proper fit

Code	Size	(cm)	DVA:
FBS301239	S	53.7 - 63.8	AR03
FBS301239M	M	63.8 - 76.5	AR03
FBS301239L	L	76.5 - 89.2	AR03
FBS301239XL	XL	89.2 - 104.5	AR03
FBS301239XXL	XXL	104.5 - 117.2	AR03
FBS301239XXXL	XXXL	117.2 - 129.9	AR03

Note: Measure around scapula under arm.



Oppo 4072 Humerus Brace

- Padded webbing strap fits under the opposite axilla
- Hook and loop closure for quick, easy, adjustable application and removal

Code	Size	Mid Upper Arm (cm)	DVA:
FBS301240	S	24 - 29	AR03
FBS301240M	M	29 - 34	AR03
FBS301240L	L	34 - 39	AR03

Note: Measure mid upper arm circumference (cm).



Rehband 5055 Acro Comfort Shoulder Support

- Stabilises the shoulder joint
- Stretch flex system is anatomically designed for optimal function
- Excellent comfort through the temperature regulating features of comfortemp®
- Fits left or right

Code	Size	(cm)	DVA:
FBS301257	X-Small	24 - 28	AR03
FBS301257S	Small	28 - 32	AR03
FBS301257M	Medium	32 - 35	AR03
FBS301257L	Large	35 - 38	AR03
FBS301257XL	X-Large	38 - 42	AR03

Note: Measure circumference around the shoulder joint going under the arm (cm).



Rehband 7126 Shoulder Support

- Soft and flexible shoulder support with excellent fit
- Provides warmth and pain relief to the shoulder joint

Code		Size	Chest (cm)	DVA:
Left	Right			
FBS301235	FBS301235XSR	X-Small	24 - 28	AR03
FBS301235SL	FBS301235SR	Small	28 - 32	AR03
FBS301235ML	FBS301235MR	Medium	32 - 35	AR03
FBS301235LL	FBS301235LR	Large	35 - 38	AR03
FBS301235XLL	FBS301235XLR	X-Large	38 - 42	AR03

Note: Measure around the chest (cm).

DVA Contact Details

P 1300 888 052 | F 1300 787 052 | E dva@aidacare.com.au



For NDIS Item Code, refer to page: 500





Rolyan MFC II Shoulder Orthosis

- Lycra/spandex fabric offers excellent conformability to body contours for a snug fit
- Easy to put on and take off
- Can be worn under clothing
- Hand wash, air dry

Code		Size	Mid Upper Arm (cm)	DVA:
Left	Right			
FBS301251	FBS301252	X-Small	23 - 28	AR03
FBS301253	FBS301254	Small	28 - 33	AR03
FBS301255	FBS301256	Medium	33 - 38	AR03

Note: Measure circumference of mid upper-arm.



Rolyan Hemi Arm Sling II

- Positions the head of the humerus in the glenoid fossa and the scapula in external rotation
- Optional universal cuff adapter A653-11 available - positions the hand above heart to treat or prevent oedema. Attaches to a d-ring on the unaffected side
- Length 18cm, colour white
- Elastic straps are dynamic, allowing a patient to actively move returning muscles and perform self-ranging exercises
- Approximately 3.8cm longer than the original model

Code		Size	(cm)	Chest (cm)	DVA:
Left	Right				
FBS301241	FBS301242	Small	36 - 46	-	AR03
FBS301243	FBS301244	Medium	46 - 61	78 - 89	AR03
FBS301245	FBS301246	Medium - Long	46 - 61	89 - 102	AR03
FBS301247	FBS301248	Long	61 - 71	-	AR03
FBS301249	FBS301250	X-Large	71 - 86	-	AR03

Note: To fit, measure the distance from acromion to acromion posteriorly across the midpoint of the scapula. To determine if patient is medium or medium long, measure the chest circumference at the axilla.





Bracing & Supports - Trunk & Back

I ♥ ndis For NDIS Item Code, refer to page: 500



Oppo 2060 Abdominal Binder - Firm

- Made of a combination of elastic and soft nylon pile bonded to 6mm foam for comfort
- Elastic panel construction provides even support; prevents "roll over"
- Contoured to fit between waist and hips



Oppo 2265 Lumbar Sacral Support - Tapered

- Four removable posterior flexible aluminum stays support the lower back
- Two flexible stays in front help reduce rolling of the support when sitting
- Allows limited flexion and extension
- 25cm high at back, tapered to 17.5cm at front

Code	Size	Waist (cm)	DVA:
FBS301269	Small	61.6 - 81.3	AR14
FBS301269M	Medium	81.3 - 94	AR14
FBS301269L	Large	94 - 104.1	AR14
FBS301269XL	X-Large	104.1 - 114.3	AR14
FBS301269XXL	XX-Large	114.3+	AR14

Note: Measure around the waist.

Code	Size	Waist (cm)	DVA:
FBS301278	Small	58.4 - 71.1	AR08
FBS301278M	Medium	71.1 - 83.8	AR08
FBS301278L	Large	83.8 - 96.5	AR08
FBS301278XL	X-Large	96.5 - 109.2	AR08
FBS301278XXL	XX-Large	109.2 - 121.9	AR08

Note: Measure around the waist.



Rehband 1068 EQ Posture Support

- A flexible aid for training a better posture
- Bendable stays for additional support
- Elastic and well-fitting
- Wear with singlet or t-shirt underneath



Oppo 2175 Posture Aid

- Triangle configuration correctly distributes pressure to tip of shoulder
- Suitable for adults or children
- Adjustable hook and loop closure facilitates proper fit

Code	Size	Waist (cm)	DVA:
FBS301258	X-Small	56 - 64	AR03
FBS301258S	Small	64 - 72	AR03
FBS301258M	Medium	72 - 80	AR03
FBS301258L	Large	80 - 90	AR03
FBS301258XL	X-Large	90 - 100	AR03

Note: Measure waist circumference (cm).

Code	Size	Waist (cm)	DVA:
FBS301259	Small	64.1 - 76.5	AR03
FBS301259M	Medium	76.5 - 89.2	AR03
FBS301259L	Large	89.2 - 104.5	AR03
FBS301259XL	X-Large	104.5 - 117.2	AR03
FBS301259XXL	XX-Large	117.2 - 129.9	AR03

Note: Measure around scapula under arm.

DVA Contact Details

P 1300 888 052 | F 1300 787 052 | E dva@aidacare.com.au

I ♥ ndis For NDIS Item Code, refer to page: 500



Dale Abdominal Binder

- Single ply stretch loop elastic (latex-free)
- Velcro closure anywhere on belt for perfect fit
- Provides full, even pressure and aids abdominal muscle
- Colour: white
- Sizing: available in either 23cm (pictured) or 30cm widths, both in S, M and L sizes



Dale Abdominal Lumbar Support

- Single ply stretch loop elastic (latex-free)
- Two removable aluminium supports
- Supports can be hand-shaped to hip curvature
- Velcro closure anywhere on belt for perfect fit
- Colour: white
- Sizing: available in either 23cm (pictured) or 30cm widths, both in S, M and L sizes

Code		Size	Waist/Hips (cm)	DVA:
23cm Width	30cm Width			
FBS301263	FBS301266	S	Up to 114	AR14
FBS301264	FBS301267	M	114 - 157	AR14
FBS301265	FBS301268	L	152 - 191	AR14

Note: Measure circumference of male waist/female hips.

Code		Size	Waist/Hips (cm)	DVA:
23cm Width	30cm Width			
FBS301270	FBS301273	S	Up to 114	AR08
FBS301271	FBS301274	M	114 - 157	AR08
FBS301272	FBS301275	L	152 - 191	AR08

Note: Measure circumference of male waist/female hips.



Oppo 2263 Lumbar Sacral Support

- Built-in flexible stays provide additional support and prevent rolling
- Double layer design for adjustable compression
- Ventilated material for breath-ability and comfort
- Flexible stays provide stabilising support to the abdominal muscles
- Two flexible stays keep the belt smooth in the front for comfort when sitting
- Removable pad provides extra compression to the lower back

Code	Size	Waist (cm)	DVA:
FBS301276	Small	58.4 - 71.1	AR08
FBS301276M	Medium	71.1 - 83.8	AR08
FBS301276L	Large	83.8 - 96.5	AR08
FBS301276XL	X-Large	96.5 - 109.2	AR08
FBS301276XXL	XX-Large	109.2 - 121.9	AR08

Note: Measure around the waist.



Oppo 2267 Lumbar Sacral Support - Light

- Extra stays provide effective back support
- Elastic reinforced straps for a perfect fit
- Breathable elastic mesh fabric ensures comfort
- May be worn under clothing and will not hinder daily activities



Oppo 2264 Lumbar Sacral Support - Foam

- Designed with straps crossing at the back and fastened in front
- Comfortable to wear and provides excellent support
- Four pliable support splints at the back for extra stability
- Breathable elastic/cotton and soft pile bonded to foam
- 22.8cm back depth

Code	Size	Waist (cm)	DVA:
FBS301277	Small	58.4 - 71.1	AR08
FBS301277M	Medium	71.1 - 83.8	AR08
FBS301277L	Large	83.8 - 96.5	AR08

Note: Measure around the waist.

Code	Size	Waist (cm)	DVA:
FBS301279	Small	58.4 - 71.1	AR08
FBS301279M	Medium	71.1 - 83.8	AR08
FBS301279L	Large	83.8 - 96.5	AR08
FBS301279XL	X-Large	96.5 - 109.2	AR08
FBS301279XXL	XX-Large	109.2 - 121.9	AR08
FBS301279XXXL	XXX-Large	122 - 135	AR08

Note: Measure around the waist.



Oppo 2061 Lumbar Sacral Belt

- Provides support and variable compression to the sacrum and sacroiliac region of the spine
- 50% cotton blend
- Used for therapeutic or prophylactic support
- Ideal for sports medicine



Oppo 2065 Lumbar Sacral Support

- Designed with straps crossing at the back and fastened in front
- Comfortable to wear with 59% cotton content, and provides excellent support
- Four pliable support splints at the back for extra stability
- Choice of 2 back heights, 23cm and 30cm (pictured)

Code	Description	Size	Waist (cm)	DVA:
FBS301280	Oppo 2061 Lumbar Sacral Belt	One size	66 - 91.4	AR08

Note: Measure around the waist.

Code		Size	Waist (cm)	DVA:
23cm Width	30cm Width			
FBS301282	FBS301283	S	58.4 - 71.1	AR08
FBS301282M	FBS301283M	M	71.1 - 83.8	AR08
FBS301282L	FBS301283L	L	83.8 - 96.5	AR08
FBS301282XL	FBS301283XL	XL	96.5 - 109.2	AR08
FBS301282XXL	FBS301283XXL	XXL	109.2 - 121.9	AR08

Note: Measure around the waist.



Oppo 2164 Lumbar Sacral Support - Multiband

- Designed with straps crossing at the back and fastened in front
- Comfortable to wear, and provides excellent support
- Four pliable support splints at the back for extra stability
- Breathable elastic/cotton blend and soft pile bonded to foam



Oppo 2261 Light Industrial Back Support

- Two extra stays on the back provide effective back support
- May be worn under clothing and will not hinder daily activities
- Breathable elastic/cotton mesh fabric ensures comfort
- Elastic reinforced straps for a custom fit
- Back height 25.4cm

Code	Size	Waist (cm)	DVA:
FBS301284	Small	58.4 - 71.1	AR08
FBS301284M	Medium	71.1 - 83.8	AR08
FBS301284L	Large	83.8 - 96.5	AR08
FBS301284XL	X-Large	96.5 - 109.2	AR08
FBS301284XXL	XX-Large	109.2 - 121.9	AR08

Note: Measure around the waist.

Code	Size	Waist (cm)	DVA:
FBS301285	Small	58.4 - 71.1	AR08
FBS301285M	Medium	71.1 - 83.8	AR08
FBS301285L	Large	83.8 - 96.5	AR08
FBS301285XL	X-Large	96.5 - 109.2	AR08
FBS301285XXL	XX-Large	109.2 - 121.9	AR08

Note: Measure around the waist.



Oppo 2169 Back Support/Shoulder Straps

- Helps prevent and relieve lower back pain
- Easy to apply and remove
- Shoulder straps enable user to affix when needed for lifting or bending
- May be used for both the lumbar and lumbosacral regions
- Four pliable support splints in the back for extra stability



Rehband 1069 Dosi EQX Straight Fit

- Classic back support for straighter figure
- 22cm width
- Four stabilising splints in back
- Adjustable cross-strap for individual support

Code	Size	Waist (cm)	DVA:
FBS301286	Small	73.7 - 83.8	AR08
FBS301286M	Medium	83.8 - 96.5	AR08
FBS301286L	Large	96.5 - 109.2	AR08
FBS301286XL	X-Large	109.2 - 124.5	AR08
FBS301286XXL	XX-Large	124.5 - 137.2	AR08

Note: Measure around the waist.

Code	Size	Waist (cm)	DVA:
FBS301287	X-Small	56 - 64	AR08
FBS301287S	Small	64 - 72	AR08
FBS301287M	Medium	72 - 80	AR08
FBS301287L	Large	80 - 90	AR08
FBS301287XL	X-Large	90 - 100	AR08
FBS301287XXL	XX-Large	100 - 110	AR08

Note: Measure around the waist.

DVA Contact Details

P 1300 888 052 | F 1300 787 052 | E dva@aidacare.com.au



For NDIS Item Code, refer to page: 500





Rehband 1079 Dosi EQ SI Shaped Fit

- Classic back support for curved (female) figure
- 22cm width
- Four stabilising splints in back
- Adjustable cross-strap for individual support

Code	Size	Waist (cm)	DVA:
FBS301288	X-Small	56 - 64	AR08
FBS301288S	Small	64 - 72	AR08
FBS301288M	Medium	72 - 80	AR08
FBS301288L	Large	80 - 90	AR08
FBS301288XL	X-Large	90 - 100	AR08
FBS301288XXL	XX-Large	100 - 110	AR08

Note: Measure around the waist.



Rehband 1070 Dosi EQ SI High Fit

- Classic back support, in high fit
- 30cm width
- Covers the lower back and half of the thoracic
- Four stabilising splints in back
- Adjustable cross-strap for individual support

Code	Size	Waist (cm)	DVA:
FBS301289	Small	64 - 72	AR08
FBS301289M	Medium	72 - 80	AR08
FBS301289L	Large	80 - 90	AR08
FBS301289XL	X-Large	90 - 100	AR08
FBS301289XXL	XX-Large	100 - 110	AR08

Note: Measure around the waist.



BodyAssist 109 Sacro Cinch

- Cool, breathable lightweight lycra and elastic construction
- Incorporates two layers of closure, giving firm but comfortable support to the lower back
- Incorporates pliable support stays for extra stability

Code		Size	Waist (cm)	DVA:
Beige	Black			
FBS283620	FBS283662	X-Small	65 - 75	AR08
FBS283630	FBS283663	Small	75 - 85	AR08
FBS283640	FBS283664	Medium	85 - 95	AR08
FBS283650	FBS283665	Large	95 - 112	AR08
FBS283660	FBS283666	X-Large	112 - 130	AR08

Note: Measure around the waist.



BodyAssist 105 Deluxe Sacro Cinch Support

- Elastic / lycra support with two layers of enclosure
- Front split design for perfect fit of waist/hips
- The upper side pull closure, wraps neatly around the waist, while the lower straps close perfectly over the hips

Code	Size	Waist (cm)	DVA:
FBS280529	Small	75 - 85	AR08
FBS280530	Medium	85 - 95	AR08
FBS280531	Large	95 - 105	AR08
FBS280532	X-Large	105 - 115	AR08

Note: Measure around the waist.



Elastic Sacro Cinch Support

- Made of firm elastic with webbing cinch straps
- Reduces intradiscal pressure by increasing intra-abdominal pressure
- Velcro closure with cinch pull strap
- Two posterior metal stays, which need to be shaped at time of fitting
- Back length 25cm, front length 19cm
- Straighter male fit, and shaped female fit

Code		Size	Waist (cm)	DVA:
Male	Female			
FBS281610M	FBS281610	Small	76 - 86	AR08
FBS281620M	FBS281620	Medium	87 - 96	AR08
FBS281630M	FBS281630	Large	97 - 106	AR08
FBS281640M	FBS281640	X-Large	107 - 116	AR08
FBS281650M	FBS281650	XX-Large	117 - 126	AR08

Note: Measure around the waist.



Elastic Sacro Cinch Support - Long Fit

- Made of firm elastic with webbing cinch straps
- Reduces intradiscal pressure by increasing intra-abdominal pressure
- Velcro closure with cinch pull strap
- Two posterior metal stays, which need to be shaped at time of fitting
- Back length 36cm, front length 22cm

Code	Size	Waist (cm)	DVA:
FBS301290	Small	76 - 86	AR08
FBS301290M	Medium	87 - 96	AR08
FBS301290L	Large	97 - 106	AR08
FBS301290XL	X-Large	107 - 116	AR08
FBS301290XXL	XX-Large	117 - 126	AR08

Note: Measure around the waist.



Oppo 1063 Neoprene Back Belt - Univeoprene

- Provides maximum support to mid and lower back without restricting normal movement
- Three flexible stays on the back provide extra support
- Retains body heat for constant moist heat to lower back muscle group
- Hook and loop closure allows for a wide range of adjustable compression

Code	Description	Size	Waist (cm)	DVA:
FBS301293	Oppo 1063 Neoprene Back Belt - Univ (Beige)	One size	71.1 - 116.8	AR08

Note: Measure around the waist.



Oppo 1064 Lumbar Support - Neoprene

- Provides support for lower back injuries
- Maximum support to mid and lower back without restricting normal movement
- Three flexible stays on the back provide extra support
- Hook and loop closure allows for a wide range of adjustable compression
- Two extra straps provide more support and stability
- 25cm width back shape provides extra support

Code	Size	Waist (cm)	DVA:
FBS301294	Small	61.0 - 78.7	AR08
FBS301294M	Medium	78.7 - 89.5	AR08
FBS301294L	Large	89.5 - 99.7	AR08
FBS301294XL	X-Large	99.7 - 110.5	AR08
FBS301294XXL	XX-Large	110.5 - 121.3	AR08

Note: Measure around the waist.



McDavid A491 Waist Belt - Thermal

- Adjustable thermal waist belt that delivers various levels of compression with therapeutic warmth
- The shaped design and adjustable closures assures an optimum fit
- One size fits up to a 127cm or 50" waist

Code	Description	Size	DVA:
FBS285050	McDavid A491 Waist Belt - Thermal Black	Universal size	AR08



McDavid A493 Sacro Cinch Support - Thermal

- Made of lightweight lycra, thermal rubber & orthopaedic elastic that delivers firm support
- Two layer closure. The first layer consists of heat retaining medical neoprene across the lower back, attached to lycra panels
- The outer straps deliver firm, powerful support when attached
- When closed the first layer delivers mild compression so the brace to be worn without effect
- 5 pliable support stays

Code	Size	Waist (cm)	DVA:
FBS283565	Small	60 - 80	AR08
FBS283570	Medium	67 - 105	AR08
FBS283580	Large	90 - 132	AR08

Note: Measure around the waist.



BodyAssist 103 Lumbar Sacral Support – Neoprene

- Thermal neoprene delivers therapeutic warmth & protection to soft tissue
- Anatomically designed to fit perfectly
- Includes 4 rigid medical splints & 4 flexible spiral side stays
- 2 layers of closure
- Outer straps deliver firm, orthopaedic support
- Inner layer consists of thermal neoprene delivering minimal-to-mild compression

Code		Size	Waist (cm)	DVA:
Beige	Black			
FBS283600	FBS283595	Small - Medium	Up to 85	AR08
FBS283610	FBS283590	Large - X-Large	85 - 105	AR08

Note: Measure around the waist.



Aspen TLSO - Thoraco Lumbar Sacral Bracing System

- The Aspen TLSO is a prefabricated, custom fitted orthosis
- Constructed of three rigid plastic shells including the rigid posterior panel and two rigid side panels
- Together with the padded chest panel, shoulder straps and torso control straps, the Aspen TLSO provides triplanar motion restriction
- Requiring several specific anatomical measurements and characterised by a structurally rigid layer of high density polyethylene that encircles the trunk, the Aspen TLSO is designed to be trimmed, bent, moulded or modified for use by a specific patient

Code		Size	Back Panel Ht. (cm)	Cut-Out Ht. (cm)	Circum. (cm)	DVA:
TLSO	Replacement Pads					
FBS301304	FBS301308	S/S - Short/Small	33	20	66 - 99	AR08
FBS301305	FBS301309	S/L - Short/Large	33	20	91 - 124	AR08
FBS301306	FBS301310	T/S - Tall/Small	41	25	81 - 119	AR08
FBS301307	FBS301311	T/L - Tall/Large	41	25	112 - 150	AR08



Aspen LSO - Lumbar Sacral Bracing System

- The Aspen LSO is a prefabricated, custom fitted body jacket style orthosis made from high density polyethylene and rigid aluminum rods with a cotton-lined foam interface material

Code		Size	Back Panel Ht. (cm)	Cut-Out Ht. (cm)	Circum. (cm)	DVA:
LSO	Replacement Pads					
FBS301312	FBS301316	S/S - Short/Small	33	20	66 - 99	AR08
FBS301313	FBS301317	S/L - Short/Large	33	20	91 - 124	AR08
FBS301314	FBS301318	T/S - Tall/Small	41	25	81 - 119	AR08
FBS301315	FBS301319	T/L - Tall/Large	41	25	112 - 150	AR08



Bracing & Supports - Arm & Elbow

I ♥ ndis For NDIS Item Code, refer to page: 500



BodyAssist 95 Elbow Sleeve

- Provides a basic level of compression and orthopaedic support

Code	Description	Size	Elbow (cm)	DVA:
FBS284640	95 Elbow Sleeve	Small	20 - 25	AR03
FBS284650	95 Elbow Sleeve	Medium	25 - 28	AR03
FBS284660	95 Elbow Sleeve	Large	28 - 30	AR03
FBS284670	95 Elbow Sleeve	X-Large	30 +	AR03

Measure elbow circumference.



McDavid A488 Thermal Elbow Wrap

- Multi-purpose elbow wrap provides warmth, firm compression and some protection
- Contoured and darted for a better fit
- Velcro closure ensures the best fit
- Fully adjustable

Code	Description	Colour	DVA:
FBS285040	A488 Thermal Elbow Wrap	Beige	AR03
FBS285038	A488 Thermal Elbow Wrap	Black	AR03

Measure elbow circumference.



Oppo 1085 Elbow Sleeve

- Provides full support of elbow region in conjunction with reducing of the insertion of the extensor carpi radialis tendon
- Provides support and protection during activities which require extensive elbow use

Code	Description	Size	Elbow (cm)	DVA:
FBS301321	1085 Elbow Sleeve	Small	23.2 - 25	AR03
FBS3013212	1085 Elbow Sleeve	Medium	25.7 - 28	AR03
FBS3013213	1085 Elbow Sleeve	Large	28.3 - 30	AR03
FBS3013214	1085 Elbow Sleeve	X-Large	30.8 - 33	AR03

Measure elbow circumference.



Oppo 2080 Elbow Sleeve

- Provides support and compression to the elbow joint without restricting flexibility
- Retains heat to relieve pain

Code	Description	Size	Elbow (cm)	DVA:
FBS301322	1085 Elbow Sleeve	Small	23.2 - 26	AR03
FBS3013222	1085 Elbow Sleeve	Medium	25.7 - 31.1	AR03
FBS3013223	1085 Elbow Sleeve	Large	31.1 - 37.5	AR03
FBS3013224	1085 Elbow Sleeve	X-Large	37.5 - 43.5	AR03

Measure elbow circumference.

DVA Contact Details

P 1300 888 052 | F 1300 787 052 | E dva@aidacare.com.au

I ♥ ndis For NDIS Item Code, refer to page: 500

aidacare
HEALTHCARE EQUIPMENT



Rehband 7124 Hyper-Extension Elbow Brace

- Prevents hyper-extension of the elbow
- Constructed of warming neoprene
- Metal splints
- Cross-wise straps prevent hyper extension

Code	Description	Size	Forearm (cm)	DVA:
FBS301323	7124 Hyper-Ex Elbow Brace	X-Small	21 - 23	AR03
FBS301323S	7124 Hyper-Ex Elbow Brace	Small	23 - 25	AR03
FBS301323M	7124 Hyper-Ex Elbow Brace	Medium	25 - 28	AR03
FBS301323L	7124 Hyper-Ex Elbow Brace	Large	28 - 31	AR03
FBS301323XL	7124 Hyper-Ex Elbow Brace	X-Large	31 - 33	AR03

Measure forearm circumference.



Rolyan DynaPro Flex Elbow

- Abnormal reflex arc therapy provided by the device assists in neuro-rehabilitation
- Axial hinge with legible goniometer facilitates brace fitting
- Elbow cuff extenders available
- New fabric reduces shear/friction
- Gel elbow pad included

Code	Description	Size	Upper Arm (cm)	Lower Arm (cm)	DVA:
FBS301324	DynaPro Flex Elbow	Small	20.3 - 25.4	12.7 - 20.3	AR03
FBS301324M	DynaPro Flex Elbow	Medium	22.8 - 29.2	19.0 - 25.4	AR03
FBS301324L	DynaPro Flex Elbow	Large	27.9 - 35.5	22.8 - 30.5	AR03

Upper Arm: circ (cm) around mid upper-arm
Lower Arm: circumference (cm) around mid - lower arm.



Rolyan Progressive Elbow/Knee Splint

- Provides soft, comfortable positioning with secure support while gently increasing the range of motion of elbow or knee
- 3-part Ezeform stay needs to be formed to patient and can be reformed as needed
- Latex free, hand wash, air dry

Code	Description	Size	(cm)	DVA:
FBS301325	Progressive Elbow/Knee Splint	X-Small	Up to 27	AR03
FBS301326	Progressive Elbow/Knee Splint	Small	27 - 38	AR03
FBS301327	Progressive Elbow/Knee Splint	Medium	39 - 44	AR03
FBS301328	Progressive Elbow/Knee Splint	Large	45 - 50	AR03
FBS301329	Progressive Elbow/Knee Splint	X-Large	51 - 56	AR03

Measure arm circ 5cm above elbow.
Min 5cm clearance from axillary necessary.



Oppo 1080 Tennis Elbow Brace

- Provides full support of elbow region in conjunction with reducing of the insertion of the extensor carpi radialis tendon
- Circular pad for extra pressure on the muscle in the arm or wrist
- Adjustable compressive strap reduces tension in wrist and finger extensors during gripping
- Tennis elbow, golfer's elbow

Code	Description	Size	Elbow (cm)	DVA:
FBS301335	1080 Tennis Elbow Brace	S	23.2 - 25.7	AR03
FBS3013352	1080 Tennis Elbow Brace	M	25.7 - 28.3	AR03
FBS3013353	1080 Tennis Elbow Brace	L	28.3 - 30.8	AR03
FBS3013354	1080 Tennis Elbow Brace	XL	30.8 - 33.3	AR03

Measure elbow circumference.



Oppo 1086 Tennis Elbow Support

- Relieves pain from strained muscles in the forearm
- Circular pad for extra pressure on the muscle in the arm or wrist
- Adjustable hook and loop strap for greater tension if needed
- One size fits most
- Tennis elbow

Code	Description	Size	Elbow (cm)	DVA:
FBS301336	1086 Tennis Elbow Support	Universal	17.8 - 33	AR03

Measure elbow circumference.



Oppo 1486 Tennis Elbow Support With Silicone Pad

- Specially shaped silicone inserts surrounding the epicondyles and the adjoining radius
- Oedema and haematoma are quickly reabsorbed
- Rapidly improves the joint function
- Relieves the strain on the periarticular tissues and compresses the soft tissue area
- One size fits most

Code	Description	Size	Forearm (cm)	DVA:
FBS301337	1486 Tennis Elbow Support with Silicone Pad	Universal	17.8 - 33	AR03

Measure forearm circumference.



Oppo 3086 Tennis Elbow Brace

- Detachable inner pad may be removed for cleaning
- Raised soft foam may be positioned to reduce muscle loading
- Clasp plastic provides pressure evenly on musculature around the elbow
- Adjustable hook and loop closure for a comfortable fit and stabilisation during daily activities

Code	Description	Size	Elbow (cm)	DVA:
FBS301338	3086 Tennis Elbow Brace	S	24.1 - 25.4	AR03
FBS3013382	3086 Tennis Elbow Brace	M	25.4 - 27.3	AR03
FBS3013383	3086 Tennis Elbow Brace	L	27.3 - 29.8	AR03
FBS3013384	3086 Tennis Elbow Brace	XL	27.3 - 29.8	AR03

Measure elbow circumference.



Rehband 4801 Epiflex Tennis Elbow Brace

- Neoprene/elastic elbow support provides relief of the tendon insertion
- Easy fastening buckles allow adjustment of compression

Code	Description	Size	Elbow (cm)	DVA:
FBS301339	4801 Epiflex Tennis Elbow Brace	S	24 - 25	AR03
FBS301339M	4801 Epiflex Tennis Elbow Brace	M	26 - 27	AR03
FBS301339L	4801 Epiflex Tennis Elbow Brace	L	28 - 30	AR03
FBS301339XL	4801 Epiflex Tennis Elbow Brace	XL	31 - 33	AR03

Measure elbow circumference.



Rehband 7122 Epi Ulna Tennis Elbow Support

- Warming elbow support which provides relief of the tendon insertion
- Warming neoprene
- Adjustable width for the arm

Code	Description	Size	Elbow (cm)	DVA:
FBS301340	7122 Epi Ulna Tennis Elbow Support	XS	21 - 23	AR03
FBS301340S	7122 Epi Ulna Tennis Elbow Support	S	23 - 25	AR03
FBS301340M	7122 Epi Ulna Tennis Elbow Support	M	25 - 28	AR03
FBS301340L	7122 Epi Ulna Tennis Elbow Support	L	28 - 31	AR03
FBS301340XL	7122 Epi Ulna Tennis Elbow Support	XL	31 - 33	AR03
FBS301340XXL	7122 Epi Ulna Tennis Elbow Support	XXL	33 - 36	AR03

Measure elbow circumference.



BodyAssist 970 Tennis Elbow Band

- Rigid, non - stretch strap
- Heavy duty D - loop closure for snug fit
- Incorporates pressure pad
- One size fits all

Code	Description	Size	DVA:
FBS280880	970 Tennis Elbow Band	Universal	AR03

Measure elbow circumference.



McDavid A486 Thermal Tennis Elbow Strap

- Applies pressure over the tendons
- Rigid velcro cinch-strap closes neatly over a neoprene inner wrap
- Shaped design with adjustable touch tape closure

Code	Description	Size	DVA:
FBS285030	A486 Thermal Tennis Elbow Strap	Universal	AR03

Measure elbow circumference.



McDavid A485 Thermal Tennis Elbow Support

- Heavy-duty neoprene sleeve with rigid cinch strap
- Applies pressure over the tendon
- Delivers therapeutic heat and protection to the elbow joint

Code	Description	Size	Elbow (cm)	DVA:
FBS285018	A485 Thermal Tennis Elbow Support	S	20 - 24	AR03
FBS285019	A485 Thermal Tennis Elbow Support	M	24 - 28	AR03
FBS285020	A485 Thermal Tennis Elbow Support	L	28 - 30.5	AR03
FBS285021	A485 Thermal Tennis Elbow Support	XL	30.5 - 35.5	AR03

Measure elbow circumference.



Rolyan Takeoff Forearm Band

- Removable gel pad helps dampen vibration at the elbow
- Rinse after use to extend band life
- Latex free
- TensoWrap™ forearm band changes force's direction to bypass the injured portion of the muscle/tendon unit

Code	Description	Size	DVA:
FBS301341	Takeoff Forearm Band	Universal	AR03



Oppo 3089 Arm Sling

- Supports the arm and carries the weight across the back and shoulder
- Universal design fits right or left arm
- Adjustable slide buckle for proper height adjustment
- Comfortable lightweight material

Code	Description	Size	Forearm (cm)	DVA:
FBS301342	3089 Arm Sling	S	25.4 - 29.2	AR03
FBS3013422	3089 Arm Sling	M	29.2 - 36.8	AR03
FBS3013423	3089 Arm Sling	L	36.8 - 43.2	AR03
FBS3013424	3089 Arm Sling	XL	43.2 - 48.3	AR03

Measure the length of forearm.



Oppo 3289 Mesh Arm Sling

- Offers immobilisation and support for mild sprains/strains of the shoulder area
- Netting material for greater comfort, particularly for long-term use
- Adjustable padded straps for height adjustment

Code	Description	Size	Forearm (cm)	DVA:
FBS301347	3289 Mesh Arm Sling	S	23.5 - 25.4	AR03
FBS3013472	3289 Mesh Arm Sling	M	25.4 - 30.5	AR03
FBS3013473	3289 Mesh Arm Sling	L	30.5 - 35.6	AR03
FBS3013474	3289 Mesh Arm Sling	XL	35.5 - 40.6	AR03

Measure the length of forearm.



Oppo 3189 Arm Sling/Immobiliser

- Helps immobilise and support the arm
- Comfortably supports a sprained or broken arm or for post-surgical support
- Extra strap offers more stabilisation

Code	Description	Size	Forearm (cm)	DVA:
FBS301348	3189 Arm Sling/Immobiliser	S	25.4 - 29.2	AR03
FBS3013482	3189 Arm Sling/Immobiliser	M	29.2 - 36.8	AR03
FBS3013483	3189 Arm Sling/Immobiliser	L	36.8 - 43.2	AR03
FBS3013484	3189 Arm Sling/Immobiliser	XL	43.2 - 48.3	AR03



Rolyan Deluxe Shoulder Immobiliser

- Shoulder pad for added comfort
- Extra-deep pocket improves arm positioning
- Waist strap limits shoulder movement

Code	Size	Forearm (cm)	DVA:
FBS301349	X-Small	Up to 28	AR03
FBS301350	Small	28	AR03
FBS301351	Medium	36	AR03
FBS301352	Large	43	AR03
FBS301353	X-Large	51	AR03

Measure the length of the forearm from the elbow to the knuckles.



Collarcare Padded Arm Sling/Immobiliser

- Designed to provide support to the clavicle area or for full shoulder immobilisation
- The padded fabric provides superior body conformity and weight-bearing comfort
- Suitable for CVA patients
- 'D' buckle design allows adjustment for comfort
- Soft padded arm pouch can be folded to fit

Code	Description	Size	DVA:
FBS301360	Padded Arm Sling/Immobiliser	Child	AR03
FBS301359	Padded Arm Sling/Immobiliser	Adult	AR03



Collarcare Light Arm Sling/Immobiliser

- Designed to provide support to the clavicle area or for full shoulder immobilisation
- Suitable for CVA patients
- 'O' buckle design allows adjustment for comfort
- Circumferential waist strap allowing arm to be fixed to the body

Code	Description	Size	(cm)	DVA:
FBS301355	Light Arm Sling/Immobiliser	S	16 x 32	AR03
FBS301356	Light Arm Sling/Immobiliser	M	19 x 39	AR03
FBS301357	Light Arm Sling/Immobiliser	L	21 x 44	AR03
FBS301358	Light Arm Sling/Immobiliser	XL	24 x 50	AR03

Sizing: depth x forearm length.



Rolyan Shoulder Support Abduction Sling

- Shoulder abduction positioning with support
- Pillow portion rests on waist to support arm weight
- Y-tabs provide easy strap length adjustability
- Snap-open front panel facilitates putting support on and taking it off
- Sling positioning allows air exchange at axilla

Code	Description	Size	Forearm (cm)	DVA:
FBS301361	Shoulder Support Abduction Sling	S	28	AR03
FBS301362	Shoulder Support Abduction Sling	M	36	AR03
FBS301363	Shoulder Support Abduction Sling	L	43	AR03

Measure the length of the forearm from the elbow to the knuckles.



Rolyan Figure Of 8 Sling

- Soft strapping system distributes pressure
- Supportive for orthopaedic and neurological conditions
- Hook and loop construction provides wide adjustability

Code	Description	Size	Forearm (cm)	DVA:
FBS301364	Figure-of-8 Sling	S/M	Up to 34.5	AR03
FBS301365	Figure-of-8 Sling	L/XL	Up to 37.5	AR03

Measure elbow crease to base of the index finger.



Rolyan CVA Sling

- Machine washable
- Latex free

Code	Description	Size	DVA:
FBS301366	CVA Sling	Universal	AR03

Bracing & Supports - Wrist & Hand

I ♥ ndis For NDIS Item Code, refer to page: 500



Oppo 2083 Wrist Wrap w/Thumb Loop

- Universal design fits left and right hand
- Provides comfortable compression and support

Code	Description	Size	Wrist (cm)	DVA:
FBS301370	Oppo 2083 Wrist Wrap w/Thumb Loop	Universal Size	14 - 23	AR03

Measure wrist circumference.



Oppo 2481 Wrist Wrap w/Silicone

- Provides localised support and compression, and avoids unnecessary compression to areas where it is not required
- Tension is fully adjustable to provide from light support to total compression
- Cohesive fabric, causing no pain or irritation to the skin whilst in use, comfortable to wear and painless to remove
- 7.5 cm width
- Secure the hook closure

Code	Description	DVA:
FBS301371	Oppo 2481 Wrist Wrap w/Silicone	AR03



Oppo 1081 Wrist Wrap

- Helps relieve pain through pressure on injured wrists
- Adjustable hook and loop strap ensures comfortable fit and stabilisation during activities
- Provides firm support around weak or strained area
- Fits wrist circumference 12 - 23 cm

Code	Description	Size	DVA:
FBS301367	Oppo 1081 Wrist Wrap	Universal Size	AR03



Oppo 1083 Wrist Wrap w/Thumb Loop

- Provides comfortable compression and support for weak or stressed wrists
- Helps treatment and prevention of injury
- Retains warmth, which may aid and speed up recovery period
- Fits wrist circumference 12 - 23cm

Code	Description	Size	DVA:
FBS301368	Oppo 1083 Wrist Wrap W/Thumb Loop	Universal Size	AR03



McDavid A451 Thermal Wrist Wrap

- Multi-purpose wrist wrap providing therapeutic warmth, firm compression and some protection
- Thumb hole, fully adjustable
- Velcro closure
- Colour black
- One size fits all
- Fits wrist circumference 12 - 23cm

Code	Description	Size	DVA:
FBS284970	McDavid A451 Thermal Wrist Wrap	Universal Size	AR03



BodyAssist 250 Deluxe Wrist Support

- Compressive wrist brace
- Adjustable thumb strap

Code	Description	Size	Wrist (cm)	DVA:
FBS283910	BodyAssist 250 Deluxe Wrist Support	S	Up to 16.5	AR03
FBS283920	BodyAssist 250 Deluxe Wrist Support	M	16.5 - 19	AR03
FBS283930	BodyAssist 250 Deluxe Wrist Support	L	19 +	AR03

Measure wrist circumference.



Oppo 2081 Wrist Wrap

- Helps relieve pain through pressure of injured wrists
- Elastic hook and loop strap insures comfortable fit and stabilisation during activities
- Provides firm support around weak or strained area
- Fits wrist circumference 12 - 23cm

Code	Description	Size	DVA:
FBS301369	Oppo 2081 Wrist Wrap	Universal Size	AR03



Oppo 2288 Wrist Splint With Strap

- Semi-rigid thumb spica portion is contoured to provide support, immobilisation and thumb abduction
- Provides wrist stability with moderate motion restriction
- An elastic wrist strap gives greater wrist stability
- Aluminum stay can easily be reversed to fit right or left hand

Code	Description	Size	Wrist (cm)	DVA:
FBS301373	Oppo 2288 Wrist Splint with Strap	S	14 - 15.9	AR03
FBS3013732	Oppo 2288 Wrist Splint with Strap	M	15.9 - 18.4	AR03
FBS3013733	Oppo 2288 Wrist Splint with Strap	L	18.4 - 21	AR03
FBS3013734	Oppo 2288 Wrist Splint with Strap	XL	21 - 23.5	AR03

Measure wrist circumference.



Oppo 3082 Wrist Support With D-Rings

- Hook and loop closures ensure proper fit and easy application
- Removable aluminum splint helps relieve pain
- Specify left or right when ordering

Code		Description	Size	Wrist (cm)	DVA:
Left	Right				
FBS301374L1	FBS301374	Oppo 3082 Wrist Support with D-Rings	S	14 - 15.9	AR03
FBS301374L2	FBS301374R2	Oppo 3082 Wrist Support with D-Rings	M	15.9 - 18.4	AR03
FBS301374L3	FBS301374R3	Oppo 3082 Wrist Splint with D-Rings	L	18.4 - 21	AR03
FBS301374L4	FBS301374R4	Oppo 3082 Wrist Splint with D-Rings	XL	21 - 23.5	AR03

Measure wrist circumference.



Oppo 2182 Wrist Splint With Strap

- An elastic wrist strap gives greater wrist stability
- Reversible aluminum splint may be easily be relocated to fit right or left hand

Code	Description	Size	Wrist (cm)	DVA:
FBS301375	Oppo 2182 Wrist Splint with Strap	S	14 - 15.9	AR03
FBS3013752	Oppo 2182 Wrist Splint with Strap	M	15.9 - 18.4	AR03
FBS3013753	Oppo 2182 Wrist Splint with Strap	L	18.4 - 21	AR03
FBS3013754	Oppo 2182 Wrist Splint with Strap	XL	21 - 23.5	AR03

Measure wrist circumference.



BodyAssist 260 Deluxe Wrist/Hand Brace

- Contoured support for weak, injured wrists
- Pliable rigid splint to minimise wrist movement

Code	Description	Size	Wrist (cm)	DVA:
FBS283940	BodyAssist 260 Deluxe Wrist/Hand Brace	S	Up to 16.5	AR03
FBS283950	BodyAssist 260 Deluxe Wrist/Hand Brace	M	16.5 - 19	AR03
FBS283960	BodyAssist 260 Deluxe Wrist/Hand Brace	L	19 - 21.5	AR03
FBS283970	BodyAssist 260 Deluxe Wrist/Hand Brace	XL	22 +	AR03

Measure wrist circumference.



McDavid A454 Thermal Carpal Tunnel Wrist Splint

- Contoured thermal neoprene for weak injured wrists
- Includes a plantar splint to minimise wrist movement
- Specify left or right
- One size fits all

Code		Description	DVA:
Left	Right		
FBS285000	FBS285010	McDavid A454 Thermal Carpal Tunnel Wrist Splint	AR03



Rehband 4024 Manu Rheuma Flex Wrist

- Wrist support made of soft, skin caring and breathable material
- Supports the wrist through a flexible palmar splint
- Velcro closure simplifies applying the support for rheumatic patients
- A soft and pre-shaped plastic splint provides light support
- Straight model to fit the arthritic hand

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301377LXS	FBS301377	XS	14 - 15	AR03
FBS301377LS	FBS301377RS	S	15 - 16.5	AR03
FBS301377LM	FBS301377RM	M	16.5 - 18	AR03
FBS301377LL	FBS301377RL	L	18 - 19.5	AR03
FBS301377LXL	FBS301377RXL	XL	19.5 - 21	AR03
FBS301377LXXL	FBS301377RXXL	XXL	21 - 22.5	AR03

Measure wrist circumference.



Rehband 4055 Manu Comfort T Wrist

- Stabilising wrist brace
- Easily applied
- Excellent comfort through the temperature regulating features of ComforTemp
- Includes plastic and steel splints
- 3 cm shorter than Rehband 4058
- Specify left or right when ordering

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301378L1	FBS301378	XS - S	13 - 16	AR03
FBS301378L2	FBS301378R2	M - L	16 - 19	AR03
FBS301378L3	FBS301378R3	XL - XXL	19 - 23	AR03

Measure wrist circumference.



Rehband 4074 Diagonal Lycra Short Wrist

- A thumb support with splint to stabilise the CMC and MCP joints of the thumb
- The material offers optimal fit and comfort
- Pain relieving and unloading effects
- Specify left or right when ordering

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301382LXS	FBS301382	XS	13 - 14	AR03
FBS301382LS	FBS301382RS	S	14 - 16	AR03
FBS301382LM	FBS301382RM	M	16 - 18	AR03
FBS301382LL	FBS301382RL	L	18 - 20	AR03
FBS301382LXL	FBS301382RXL	XL	20 - 22	AR03

Measure wrist circumference.



Rehband 4088 Wrist

- Provides exceptional stability while immobilising the wrist
- Made of thermoprene material for optimal comfort
- A firm wrist stay supports the palmar side, while a dorsal pocket holds a thermoplastic insert
- Both the stay and insert can be heat-moulded
- Specify left or right when ordering

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301383LS	FBS301383	S	14 - 16	AR03
FBS301383LM	FBS301383RM	M	16 - 18	AR03
FBS301383LL	FBS301383RL	L	18 - 20	AR03
FBS301383LXL	FBS301383RXL	XL	20 - 22	AR03

Measure wrist circumference.



Rehband 4142 Manu 3D

- High quality 3 tab wrist brace with flexion stay
- Provides a light support to the basal thumb joint
- Specify left or right when ordering

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301384LXS	FBS301384	XS	14.5 - 15.5	AR03
FBS301384LS	FBS301384RS	S	15.5 - 16	AR03
FBS301384LM	FBS301384RM	M	16 - 17	AR03
FBS301384LL	FBS301384RL	L	17 - 18	AR03
FBS301384LXL	FBS301384RXL	XL	18 - 19	AR03
FBS301384LXXL	FBS301384RXXL	XXL	19 - 20	AR03

Measure wrist circumference.



Rehband 4143 Manu 3D Basis

- High quality 3 tab wrist brace with flexion stay
- Allows free thumb motion
- Six sizes left and right

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301385LXS	FBS301385	XS	14.5 - 15.5	AR03
FBS301385LS	FBS301385RS	S	15.5 - 16	AR03
FBS301385LM	FBS301385RM	M	16 - 17	AR03
FBS301385LL	FBS301385RL	L	17 - 18	AR03
FBS301385LXL	FBS301385RXL	XL	18 - 19	AR03
FBS301385LXXL	FBS301385RXXL	XXL	19 - 20	AR03

Measure wrist circumference.



Rehband 4058 Manu Comfort T Stable Wrist

- Extra stabilising wrist brace
- Easily applied
- Excellent comfort through the temperature regulating features of ComforTemp
- Includes plastic and steel splints
- 3 cm longer than Rehband 4055
- Specify left or right when ordering

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301379L1	FBS301379	XS - S	13 - 16	AR03
FBS301379L2	FBS301379R2	M - L	16 - 19	AR03
FBS301379L3	FBS301379R3	XL - XXL	19 - 23	AR03

Measure wrist circumference.



Rolyan D-Ring Wrist Long Fit - Firm

- Volar and dorsal stays for extra wrist support
- Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers and thumb
- Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation
- D-ring closure system for simple application, removal and customised tension control
- Double layer of stockinette under closures cushions straps

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301390	FBS301386	XS	Up to 14.6	AR03
FBS301391	FBS301387	S	14.6 - 16.5	AR03
FBS301392	FBS301388	M	16.6 - 19.1	AR03
FBS301393	FBS301389	L	16.6 - 19.1	AR03

Measure wrist circumference.



Rolyan D-Ring Wrist Brace

- Non-elastic construction does not restrict circulation
- Circumferential design holds the brace in position during application
- D-ring closures provide tension control to adjust amount of support
- Double layer of cotton stockinette under closures cushions straps
- Forearm length is proportional to brace size
- Removable, adjustable aluminium support bar
- Splints available in regular and long length
- Regular length splint is also available with a larger thumb hole, which accommodates larger thumbs and fluctuating oedema.
- Fabricated from a polyester/cotton fabric laminated to 1cm foam

Code		Size	Wrist (cm)	Length (cm)	DVA:
Left	Right				
Style - Regular					
FBS301394LXS	FBS301394	XS	Up to 14.6	15.9	AR03
FBS301394LS	FBS301394RS	S	14.6 - 16.5	16.5	AR03
FBS301394LM	FBS301394RM	M	16.6 - 19.1	17.8	AR03
FBS301394LL	FBS301394RL	L	19.2 - 21.6	19	AR03
FBS301394LXL	FBS301394RXL	XL	Over 21.6	19.7	AR03
Style - Regular - Large Thumb					
FBS301395LXS	FBS301395	XS	Up to 14.6	15.9	AR03
FBS301395LS	FBS301395RS	S	14.6 - 16.5	16.5	AR03
FBS301395LM	FBS301395RM	M	16.6 - 19.1	17.8	AR03
FBS301395LL	FBS301395RL	L	19.2 - 21.6	19	AR03
FBS301395LXL	FBS301395RXL	XL	Over 21.6	19.7	AR03
Style - Long					
FBS301396LXS	FBS301396	XS	Up to 14.6	21	AR03
FBS301396LS	FBS301396RS	S	14.6 - 16.5	21.6	AR03
FBS301396LM	FBS301396RM	M	16.6 - 19.1	22.2	AR03
FBS301396LL	FBS301396RL	L	19.2 - 21.6	22.6	AR03
FBS301396LXL	FBS301396RXL	XL	Over 21.6	23.2	AR03

Measure wrist circumference.



Rolyan Take Off Universal Wrist Splint

- Semi-rigid foam support with plastic stay
- Soft, breathable Fabrifoam material provides friction to keep splint from migrating and helps remove moisture
- One size fits most adults

Code		Size	DVA:
Left	Right		
FBS301405	FBS301406	Universal	AR03



Rolyan Alignrite Wrist Support w/Wrist Strap

- Provides excellent support for the wrist while allowing fingers and thumb free movement
- Woven cotton-coated rubber fabric has open weave for ventilation
- Wrap - around wrist strap provides extra support
- Hand or machine wash, air dry
- Latex free

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301407	FBS301408	XS	11.4 - 14.0	AR03
FBS301409	FBS301410	S	14.0 - 16.5	AR03
FBS301411	FBS301412	M	16.6 - 18.4	AR03
FBS301413	FBS301414	L	18.5 - 20.3	AR03
FBS301415	FBS301416	XL	20.4 - 22.2	AR03

Measure wrist circumference.



Rolyan Wrist Splint

- Unlike neoprene or perforated neoprene, breathoprene allows air through, and wicks perspiration away from the skin for enhanced comfort and compliance
- The two-way stretch 3.2 mm thick, soft material provides compression and a perfect fit
- Removable metal dorsal and palmar stays and wrist D-ring strap for added support
- Wraparound design
- Latex free

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301417	FBS301418	S	14 - 16.5	AR03
FBS301419	FBS301420	M	16.6 - 18.4	AR03
FBS301421	FBS301422	L	18.5 - 20.3	AR03
FBS301423	FBS301424	XL	20.4 - 22.2	AR03

Measure wrist circumference.



Oppo 1289 Wrist/Thumb Support w/Strap

- Three removable metal splints in thumb area for improved stability
- Additional strap for stabilisation during activities
- Supports the wrist and thumb, leaving the palm free
- Ideal for sprains of thumb, ligament support, and minor fractures of the thumb

Code	Description	Size	Wrist (cm)	DVA:
FBS301425	Oppo 1289 Wrist/Thumb Support W/Strap	S	13 - 15.2	AR03
FBS3014252	Oppo 1289 Wrist/Thumb Support W/Strap	M	15.2 - 17.2	AR03
FBS3014253	Oppo 1289 Wrist/Thumb Support W/Strap	L	17.2 - 19.1	AR03
FBS3014254	Oppo 1289 Wrist/Thumb Support W/Strap	XL	19.1 - 21	AR03

Measure wrist circumference.



Oppo 2184 Wrist/Thumb Sleeve

- Pull-on sleeve in four-way stretch tubular elastic
- Not splinted, light compression
- Universal design fits left and right hand

Code	Description	Size	Wrist (cm)	DVA:
FBS301426	Oppo 2184 Wrist/Thumb Sleeve	S	14 - 15.9	AR03
FBS3014262	Oppo 2184 Wrist/Thumb Sleeve	M	15.9 - 18.4	AR03
FBS3014263	Oppo 2184 Wrist/Thumb Sleeve	L	18.4 - 21	AR03
FBS3014264	Oppo 2184 Wrist/Thumb Sleeve	XL	21 - 23.5	AR03

Measure wrist circumference.



Oppo 3182 Wrist Support With Thumb

- Wraparound design and hook and loop straps for easy fit
- Supports and limits the motion of the thumb and wrist joint
- Elastic fabric provides compression
- Elastic hook and loop strap may be used to limit wrist motion
- Specify left or right

Code		Description	Size	Wrist (cm)	DVA:
Left	Right				
FBS301427L1	FBS301427	Oppo 3182 Wrist Support with Thumb	S	14 - 15.9	AR03
FBS301427L2	FBS301427R2	Oppo 3182 Wrist Support with Thumb	M	15.9 - 18.4	AR03
FBS301427L3	FBS301427R3	Oppo 3182 Wrist Support with Thumb	L	18.4 - 21	AR03
FBS301427L4	FBS301427R4	Oppo 3182 Wrist Support with Thumb	XL	21 - 23.5	AR03

Measure wrist circumference.



Oppo 1089 Wrist Thumb Support

- Thumb orthosis especially suitable for rheumatics
- Relieves pain and stabilises CMC joint
- Reinforced metal stick stabilises thumb
- Adjustable hook and loop fasteners insure comfortable fit and stabilisation during activities
- Thermal wrist support for tenosynovitis, carpal tunnel syndrome and general wrist support

Code	Description	Size	Wrist (cm)	DVA:
FBS301428	Oppo 1089 Wrist Thumb Support	S	14 - 15.9	AR03
FBS3014282	Oppo 1089 Wrist Thumb Support	M	15.9 - 18.4	AR03
FBS3014283	Oppo 1089 Wrist Thumb Support	L	18.4 - 21	AR03
FBS3014284	Oppo 1089 Wrist Thumb Support	XL	21 - 23.5	AR03

Measure wrist circumference.



Oppo 1188 Mouldable Wrist/Thumb Support

- Relieves pain
- Thermoplastic support stabilises thumb
- Thermoplastic insert may be removed as healing progresses
- Available in 3 lengths: 14 cm, 16.5 cm and 21 cm

Code			Size	Wrist (cm)	DVA:
Short (14cm)	Medium (16.5cm)	Long (21cm)			
FBS301429	FBS3014296-1	FBS3014298-1	S	14 - 15.9	AR03
FBS3014295-2	FBS3014296-2	FBS3014298-2	M	15.9 - 18.4	AR03
FBS3014295-3	FBS3014296-3	FBS3014298-3	L	18.4 - 21	AR03
FBS3014295-4	FBS3014296-4	FBS3014298-4	XL	21 - 23.5	AR03
FBS3014295-5	FBS3014296-5	FBS3014298-5	XXL	23.5 - 26	AR03

Measure wrist circumference.



Rehband 4025 Rheuma Wrist & Thumb Support

- Wrist/thumb support made of soft and breathable material
- Supports the wrist through a flexible palmar splint
- Velcro closure simplifies applying the support for rheumatic patients
- Straight model to fit the arthritic hand
- The placement of the velcro is chosen to avoid any pressure on caput ulnae

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301430LXS	FBS301430	XS	14 - 15	AR03
FBS301430LS	FBS301430RS	S	15 - 16.5	AR03
FBS301430LM	FBS301430RM	M	16.5 - 18	AR03
FBS301430LL	FBS301430RL	L	18 - 19.5	AR03
FBS301430LXL	FBS301430RXL	XL	19.5 - 21	AR03
FBS301430LXXL	FBS301430RXXL	XXL	21 - 22.5	AR03

Measure wrist circumference.



Rehband 4146 Manu 3D Pollex

- High quality 3 tab wrist brace with flexion stay
- Aluminium thumb stay
- Six sizes left and right

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301431LXS	FBS301431	XS	14 - 15	AR03
FBS301431LS	FBS301431RS	S	15.5 - 16	AR03
FBS301431LM	FBS301431RM	M	16 - 17	AR03
FBS301431LL	FBS301431RL	L	17 - 18	AR03
FBS301431LXL	FBS301431RXL	XL	18 - 19	AR03
FBS301431LXXL	FBS301431RXXL	XXL	19 - 20	AR03

Measure wrist circumference.



Rolyan D-Ring Wrist/Thumb Spica

- Made of cotton/polyester fabric lined with cotton stockinette
- Circumferential design with D-ring closures makes the brace easy to put on and comfortable to wear
- Two metal stays support wrist and thumb
- Specify left or right

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS282670	FBS282680	XS	Up to 14.6	AR03
FBS282690	FBS282700	S	14.6 - 16.5	AR03
FBS282710	FBS282715	M	16.6 - 19.1	AR03
FBS282720	FBS282730	L	19.2 - 21.6	AR03
FBS282725	FBS282731	XL	Over 21.6	AR03

Measure wrist circumference.



Rolyan D-Ring Wrist/Thumb Spica - Thermal

- Made from a neoprene blend with a Trioxon lining
- Maintains body heat, yet allows skin to breathe
- Circumferential design with D-ring closures makes the brace easy to put on and comfortable to wear
- Adjust the amount of wrist extension with aluminium stay

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS280191	FBS280196	XS	U to 14.6	AR03
FBS280192	FBS280197	S	14.6 - 16.5	AR03
FBS280193	FBS280199	L	19.2 - 21.6	AR03
FBS280194	FBS280201	XL	Over 21.6	AR03

Measure wrist circumference.



BodyAssist N29 Thumb Abduction Splint

- Minimises thumb movement while leaving palm and fingers free

Code	Description	Size	Wrist (cm)	DVA:
FBS301432	Body Assist N29 Thumb Abduction Splint	S/M	Up to 19	AR03
FBS3014323	Body Assist N29 Thumb Abduction Splint	L/XL	19 +	AR03

Measure wrist circumference.



BodyAssist 280 Thumb Splint

- Designed to minimise wrist and thumb movement
- Supplied with two interchangeable splints - a rigid splint, and pliable metal stay

Code	Description	Size	Wrist (cm)	DVA:
FBS280680	BodyAssist 280 Thumb Splint	S	Up to 16.5	AR03
FBS2806802	BodyAssist 280 Thumb Splint	M	16.5 - 19	AR03
FBS2806803	BodyAssist 280 Thumb Splint	L	19 - 21.5	AR03
FBS2806804	BodyAssist 280 Thumb Splint	XL	22+	AR03

Measure wrist circumference.



Rehband 4064 EQ Diagonal Lycra Thumb

- Thumb support made of durable lycra
- Features a splint to stabilise the CMC and MCP joints of the thumb
- Elasticity of the material is used diagonally for optimal fit and support

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301433LXS	FBS301433	XS	13 - 14	AR03
FBS301433LS	FBS301433RS	S	14 - 16	AR03
FBS301433LM	FBS301433RM	M	16 - 18	AR03
FBS301433LL	FBS301433RL	L	18 - 20	AR03
FBS301433LXL	FBS301433RXL	XL	20 - 22	AR03

Measure wrist circumference.



Rehband Thumboform Splint

- The thumboform thumb support comes with a heat-mouldable insert
- This splint stabilises the CMC and MCP joint of the thumb
- Available in short, long and Xlong styles

Code - Short		Code - Long		Code - XLong		Size	Wrist (cm)	DVA:
Left	Right	Left	Right	Left	Right			
FBS301434LS	FBS301434	FBS301435LS	FBS301435	FBS301436LS	FBS301436	S	14 - 16	AR03
FBS301434LM	FBS301434RM	FBS301435LM	FBS301435RM	FBS301436LM	FBS301436RM	M	16 - 18	AR03
FBS301434LL	FBS301434RL	FBS301435LL	FBS301435RL	FBS301436LL	FBS301436RL	L	18 - 20	AR03

Measure wrist circumference.



Rolyan Collum CMC Thumb Brace

- Made of quality leather with a soft suede lining for comfort and durability
- Restricts extremes of motion in the thumb joint
- Allows the user to perform routine hand functions
- The support fits over the thumb and wraps around the wrist
- It is held in place with two hook and loop fasteners

Code		Size	(cm)	DVA:
Left	Right			
FBS283201LXS	FBS283201	XS	Up to 9.5	AR03
FBS283201LS	FBS283201RS	S	9.8 - 10.8	AR03
FBS283201LM	FBS283201RM	M	11.1 - 12.1	AR03
FBS283201LL	FBS283201RL	L	over 12.4	AR03

Measure length from the top thumb-joint, along outside of thumb to the wrist-bone.



Rolyan Gel Thumb Spica Splint

- Cushions and helps protect the thumb
- Gel shell pad provides gentle pressure on the incision site to help reduce scar formation and hypersensitivity
- Optional polyform stay along radial side of thumb provides additional support to the thumb

Code		Size	(cm)	DVA:
Left	Right			
FBS301623	FBS301624	XS	Up to 9.5	AR03
FBS301625	FBS301626	S	9.5 - 10.8	AR03
FBS301627	FBS301628	M	10.9 - 12	AR03
FBS301629	FBS301630	L	Over 12.1	AR03

Measure length from the top thumb - joint, along outside of thumb to the wrist-bone.



Rolyan Take Off Thumb Splint

- Made using fabric form NuStimWrap
- The splint is available with either a standard plastic thumb stay for increased support, or a mouldable stay for custom fit
- Provides warmth, compression and support to reduce pain
- Easy to apply over joints and helps prevent overuse
- One size fits most

Code		Style	DVA:
Left	Right		
FBS301437	FBS301438	Regular	AR03
FBS301439	FBS301440	With Mouldable Insert	AR03



Rolyan Thermo-Form Thumb Splint

- These splints contain a pocket on the inside where the thermoplastic piece is inserted.
- The thermoplastic material can be moulded directly in the pocket by pouring in hot water
- The thermoplastic material becomes mouldable in approx. 1 minute, allowing the therapist to mould the splint to the thumb for proper positioning
- The splint comes with one strap

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301441	FBS301442	XS	Up to 14	AR03
FBS301443	FBS301444	S	14 - 16.5	AR03
FBS301445	FBS301446	M	16.5 - 19.1	AR03
FBS301447	FBS301448	L	19.1 - 21.6	AR03
FBS301449	FBS301450	XL	21.6 - 23	AR03

Measure wrist circumference.



Rolyan Polycentric Hinged Ulnar Deviation Splint

- This splint leaves the palm free for gripping objects
- The polycentric hinge allows full flexion and extension of the fingers
- Padded finger stays are malleable for ease of application and to provide optimal alignment
- Ideal for use following MP arthroplasty

Code		Size	(cm)	DVA:
Left	Right			
FBS301452L1	FBS301452	XS	5.7 - 7	AR03
FBS301452L2	FBS301452R2	S	7 - 7.6	AR03
FBS301452L3	FBS301452R3	M	7.6 - 8.6	AR03
FBS301452L4	FBS301452R4	L	8.6 - 9.2	AR03

Measure width from the 2nd to 5th MP to size.



Rolyan Radial Hinged Deviation Splint

- Soft finger straps permit each finger to be individually aligned for maximum control
- Lined in soft felt, dorsal based splint leaves palm free for function

Code		Size	(cm)	DVA:
Left	Right			
FBS301453L2	FBS301453	S	7 - 7.6	AR03
FBS301453L3	FBS301453R3	M	7.6 - 8.6	AR03
FBS301453L4	FBS301453R4	L	8.6 - 9.2	AR03

Measure width from the 2nd to 5th MP to size.



Rolyan Comforter Splint

- The Comforter splint holds the wrist and MP joints of the fingers in neutral to reduce stress on muscles and tendons
- Malleable padded stays allow each finger to be individually aligned
- IP joints are free to allow the wearer to manage bedding and night time activities
- The softly padded cover is removable for washing

Code		Size	(cm)	DVA:
Left	Right			
FBS301454L2	FBS301454	S	7 - 7.6	AR03
FBS301454L3	FBS301454R3	M	7.6 - 8.6	AR03
FBS301454L4	FBS301454R4	L	8.6 - 9.2	AR03

Measure width from the 2nd to 5th MP to size.



Rolyan Hand-Based in-Line Splint

- MCP support straps have reinforced neoprene to prevent over-stretching caused by continual use
- Finger straps support each finger individually
- 2.4 mm neoprene gives needed support without bulkiness
- Neoprene fabric retains warmth, increasing comfort

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301455	FBS301456	XS	Up to 14.6	AR03
FBS301457	FBS301458	S	14.6 - 16.5	AR03
FBS301459	FBS301460	M	16.6 - 19.1	AR03
FBS301461	FBS301462	L	Over 19.1	AR03

Measure wrist circumference.



Finger Flexor Splint

- Aids flexion at the proximal and distal interphalangeal joints on an individual finger

Code	Description	Size	Length	DVA:
FBS301463	Finger Flexor Splint	S	Child	AR03
FBS301463M	Finger Flexor Splint	M	Short	AR03



Finger Extension Splint

- A dynamic extension splint which resists flexion of the interphalangeal joints and assists in their extension

Code	Description	Size	Length	DVA:
FBS301464	Finger Extension Splint	S	Child	AR03
FBS301464M	Finger Extension Splint	M	Short	AR03
FBS301464L	Finger Extension Splint	L	Long	AR03



Oppo 4285 Finger Splint

- Lightweight aluminum protects injured finger
- Ventilated foam padding
- Folds over the end of the finger to provide protection and stabilisation

Code	Description	Size	Finger (cm)	DVA:
FBS301465	Oppo 4285 Finger Splint	S	5.5	AR03
FBS3014652	Oppo 4285 Finger Splint	M	7.5	AR03
FBS3014653	Oppo 4285 Finger Splint	L	9.5	AR03

Length of finger (cm).



Oppo 4280 Baseball Finger Splint

- Help stabilise the phalangeal joints in position
- Malleable aluminum splint can be reshaped to hold finger at an exact degree of flexion or extension
- Provides firm support while injury heals

Code	Description	Size	Finger (cm)	DVA:
FBS301466	Oppo 4280 Baseball Finger Splint	S	8	AR03
FBS3014662	Oppo 4280 Baseball Finger Splint	M	10	AR03
FBS3014663	Oppo 4280 Baseball Finger Splint	L	12	AR03

Length of finger (cm).



Oppo 4282 Curved Finger Splint

- Lined and padded with soft, ventilated material
- Curved finger splints support and protect dislocated fracture
- Upward curve at the tip to protect the end of the finger

Code	Description	Size	Finger (cm)	DVA:
FBS301467	Oppo 4282 Curved Finger Splint	S	7	AR03
FBS3014672	Oppo 4282 Curved Finger Splint	M	10	AR03
FBS3014673	Oppo 4282 Curved Finger Splint	L	13	AR03

Length of finger (cm).



Oppo 4283 Open Thimble Splint

- Malleable aluminum splint for immobilisation of finger
- Protects and stabilises fractures of the distal finger joint
- Encloses finger injuries and protects on four sides

Code	Description	Size	Finger (cm)	DVA:
FBS301469	Oppo 4283 Open Thimble Splint	S	5.5	AR03
FBS3014692	Oppo 4283 Open Thimble Splint	M	6.5	AR03
FBS3014693	Oppo 4283 Open Thimble Splint	L	7.5	AR03

Length of finger (cm).



Oppo 4182 Wrist/Hand Positioning Splint

- Two unstretchable straps give greater wrist stability
- Incorporating a moulded metal palmar stay which is used to position the wrist and hand
- Insert each two fingers, except thumb, into the two unstretchable straps
- Be sure the metal stay is on palm side of hand, then wrap the strap over hand and secure snugly

- Thread through the buckles and secure the two hook and loop straps to desired compression

Code		Description	Size	Wrist (cm)	DVA:
Left	Right				
FBS301470LS	FBS301470	OPPO 4182 Wrist/Hand Positioning Splint	S	13.1 - 17.1	AR03
FBS301470LM	FBS301470RM	OPPO 4182 Wrist/Hand Positioning Splint	M	17.1 - 21.6	AR03
FBS301470LL	FBS301470RL	OPPO 4182 Wrist/Hand Positioning Splint	L	21.6 - 30.5	AR03

Measure wrist circumference.



Rolyan Palm Protector With Finger Separator

- Made from closed-cell foam, lined with pile
- Soft finger separator strap weaves between fingers and can be cut to adjust fit
- Hand washable
- One size fits all

Code		Description	DVA:
Left	Right		
FBS282830	FBS282840	Rolyan Palm Protector with Finger Separator	AR03



Rolyan Palm Shield - Foam

- This foam support wraps around the user's wrist and hand and holds the thumb and fingers in place
- It is held in place with a Velcro closure at the wrist
- The support has four rounded foam sections that separate and position the fingers
- It also has a padded section that rests in the user's palm, and supports and positions the user's thumb

- This hand support is designed for users with moderate to severe hand tightness
- One size fits most adults - left and right fittings

Code		Description	DVA:
Left	Right		
FBS283010	FBS283020	Rolyan Palm Shield - Foam	AR03



Rolyan Palm Shield-Sof-Gel

- A palm shield with gel which helps prevent skin breakdown in the palm and between the fingers
- Gel in finger separators and palm won't bottom out giving the product a longer life
- Ideal for moderate to severe hand tightness
- Helps prevent finger contractures
- One size fits most adults - left and right fittings


Code		Description	DVA:
Left	Right		
FBS283030	FBS283040	Rolyan Palm Shield - Sof - Gel	AR03



Rolyan Palm Protector

- Holds fingers out to prevent digging into the palm of the hand
- Cuff fabricated from closed-cell foam lined with a pile fabric
- Pre-rolled and attached to simplify application
- Hand wash, air dry
- One size for each of left and right hand

Code		Description	DVA:
Left	Right		
FBS282860	FBS282880	Rolyan Palm Protector	AR03



Rolyan Carve-It Hand Positioning Splint

- For use with strokes, Alzheimer's disease, dementia and other neurological impairments
- Helps anatomically position the hand with mild to moderate finger spasticity
- Design allows use with left and right hand
- Easy to modify by cutting foam with an electric knife for a custom fit
- Flat bottom provides easy positioning in bed or on lap tray, arm tray etc
- Removable cover can be hand washed and air dried

Code		Size	(cm)	DVA:
Splint	Replacement Cover			
FBS283075	-	Sm/Med	Up to 7.6	AR03
FBS283070	FBS301472	Lge/XL	7.6 +	AR03

For size, measure width of 2nd to 5th MCP.



Rolyan Progressive Palm Protector Splint (With Optional Wrist Support)

Choice Guide			Code				Replacement Pads	DVA:
			Palm Splint		Palm Splint with Wrist Support			
Series	Foam Size (cm)	MCP Flexion	Left	Right	Left	Right		
1	2.5	90°	FBS301473	FBS301474	FBS282890	FBS282900	FBS301479	AR03
2	5.1	60°	FBS301475	FBS301476	FBS282910	FBS282920	FBS301480	AR03
2	7.6	30°	FBS301477	FBS301478	FBS282930	FBS282940	FBS301481	AR03

Indications:

Series 1	Tight hands with moderate to severe spasticity or contractures
Series 2	Hands with mild to moderate spasticity or contractures
Series 3	Hands with mild spasticity or contractures, or non - spastic hands, to provide a functional position in a soft splint

- Progressive soft splint for neurological conditions and arthritis
- Adjust splint as condition improves by removing foam palm piece and replacing it with a larger series
- Variable strapping system allows strapping thumb only, fingers only, or both thumb and fingers
- Each splint comes with one foam palm piece
- Splint base is large enough to accommodate all hand sizes and can be modified for a smaller hand by trimming the edges
- Wrist support style has a removable Ezeform stay that can be moulded for wrist extension and D-ring attachments for secure fit
- Available in universal adult size
- Hand wash splint, air dry





Rolyan Kwik-Form Functional Resting Splint

- Can be custom fitted without tools or heat - simply bend the wire core to the desired position
- Supports the wrist, thumb and fingers in a functional or resting position
- Features hook and loop closures and removable machine washable cover

Code		Size	Wrist (cm)	DVA:
Left	Right			
FBS301482	FBS301483	Sm/Med	Up to 18	AR03
FBS301484	FBS301485	Med/Lge	18cm +	AR03

Measure wrist circumference.



Rolyan Kwik-Form Progressive Hand Splint

- Can be custom fitted without tools or heat - simply bend the wire core to the desired position
- Supports the hand and wrist and provides progressive extension of the fingers
- Includes 3 hand rolls - as condition progresses, simply replace the hand roll
- For especially tight hands, can be used without hand roll

- Removable cover sewn from comfortable, cool soft material for maximum compliance
- Features hook and loop closures and removable machine washable cover

Splint Code	Replacement Cover Code	Size	Wrist (cm)	DVA:
FBS301490	FBS301492	Sm/Med	Up to 18	AR03
FBS301491	FBS301493	Med/Lge	Up to 18+	AR03

Measure wrist circumference.



Rolyan AirPro Air Graduate Splint

- Air Therapy-low load prolonged stretch
- Conforms to the shape of the hand accommodating different degrees of finger contracture
- Ideal for severely contracted fingers
- Removable cone shaped palmar air bladder can be used independently of wrist support
- Extra cover included
- Universal size

Code		Description	DVA:
Left	Right		
FBS301494L	FBS301494	Rolyan AirPro Air Graduate Splint	AR03

Universal size, left and right.



Rolyan SoffPro Functional Resting Hand

- Treats moderate contractures of the wrist/hand/fingers
- Easily adjustable for ulnar radial deviation
- Finger separator and ulnar drift strap available separately to allow custom fitting

Code		Size	Width (cm)	Length (cm)	DVA:
Left	Right				
FBS301495LS	FBS301495	Small	Up to 7.6	Up to 15.2	AR03
FBS301495LM	FBS301495RM	Medium	7.6 - 8.9	15.2 - 20.3	AR03
FBS301495LL	FBS301495RL	Large	8.9 - 10.2	20.3 - 24.1	AR03
FBS301495LC	FBS301495RC	Child	5 - 7.6	17.8	AR03
FBS301495LY	FBS301495RY	Youth	5 - 7.6	20.3	AR03

Width: MP joints (cm). Length: wrist to finger tip (cm).



Rolyan SoffPro Grip

- Treats moderate to severe contractures of the wrist/hand/fingers
- Easily moulded without tools yet sturdy enough to maintain shape during wear
- Finger separator and ulnar rift strap available separately to allow custom fitting

Code		Size	Width (cm)	Length (cm)	DVA:
Left	Right				
FBS301496LS	FBS301496	Small	Up to 7.6	Up to 15.2	AR03
FBS301496LM	FBS301496RM	Medium	7.6 - 8.9	15.2 - 20.3	AR03
FBS301496LL	FBS301496RL	Large	8.9 - 10.2	20.3 - 24.1	AR03

Width: MP joints (cm). Length: wrist to finger tip (cm).



- For use with SoftPro Grip and SoftPro functional resting hand splints
- Provides finger separation and mild abduction
- Protects finger skin integrity
- Universal size

Code	Description	Size	DVA:
FBS301497	Rolyan Softpro Finger Separator (Accessory to SoftPro Splint)	Universal	AR03



Rolyan SoftPro Ulnar Drift Strap (Accessory To SoftPro Splint)

- For use with SoftPro Grip and SoftPro functional resting hand splints
- Holds fingers in comfortable, neutral position

Code	Description	Size	DVA:
FBS301498	Rolyan Softpro Ulnar Drift Strap (Accessory to SoftPro Splint)	Universal	AR03



Cosmac Cock-Up Resting Hand Splint

- Pre-shaped plastic with velcro fastening straps
- Use for carpal tunnel syndrome, RSI, resting splint

Code		Size	Width (cm)
Left	Right		
FBS301499LM	FBS301499	S/M	7.5
FBS301499LL	FBS301499RL	L	9

Measurement: At the widest part of the knuckles



Cosmac Cock-Up Resting Hand Splint With Thumb

- Pre-shaped plastic (with thumb) with velcro fastening straps
- Use for carpal tunnel syndrome, RSI, resting splint

Code		Size	Width (cm)
Left	Right		
FBS301500LM	FBS301500	S/M	7.5
FBS301500LL	FBS301500RL	L	9

Measurement: At the widest part of the knuckles



IMAK Arthritis Gloves

- Design to help relieve aches, pains and stiffness associated with arthritis of the hands
- Provides mild compression for warmth and helps increase circulation
- Made of soft, breathable cotton and comfortable enough to wear all day and night
- Sold as a pair

Code	Size	Width (cm)	DVA:
FBS303522	XS	up to 6.9	AR03
FBS303519	S	up to 7.9	AR03
FBS303520	M	up to 8.9	AR03
FBS303521	L	up to 10	AR03
FBS303523	XL	up to 11.4	AR03

Measurement: At the widest part of the knuckles



IMAK Active Gloves

- Provide mild compression to enhance blood circulation
- Made of soft, breathable cotton and comfortable enough to wear all day and night
- Anti-slip design on fingers and palm provides excellent grip
- Sold as a pair

Code	Size	Width (cm)	DVA:
FBS303524	S	up to 7.9	AR03
FBS303525	M	up to 8.9	AR03
FBS303526	L	up to 10	AR03

Measurement: At the widest part of the knuckles



IMAK SmartGlove

- Helps prevent and relieve carpal tunnel syndrome
- Increases comfort and ensures an ergonomically correct wrist position
- Cushions and protects the underside of the wrist, including the pisiform
- Reversible to fit either the left or right hand
- Sold individually

Code	Size	Width (cm)	DVA:
FBS303511	XS	up to 6.9	AR03
FBS303512	S	up to 8.2	AR03
FBS303513	M	up to 9.5	AR03
FBS303514	L	up to 10.8	AR03

Measurement: At the widest part of the knuckles



IMAK SmartGlove with Thumb

- Helps prevent and relieve carpal tunnel syndrome, De Quervain's and other wrist and thumb pain
- Increases comfort and ensures an ergonomically correct wrist position
- Features a flexible support splint to support your thumb
- Sold individually

Code	Size	Width (cm)	DVA:
FBS303515	XS	up to 6.9	AR03
FBS303516	S	up to 8.2	AR03
FBS303517	M	up to 9.5	AR03
FBS303518	L	up to 10.8	AR03

Measurement: At the widest part of the knuckles



For Crutches
Refer to Page 261

Bracing & Supports - Knee



For NDIS Item Code, refer to page: 500



Oppo 1024 Universal Knee Support

- Designed for all knee sizes with adjustable support
- Open patella relieves pressure
- Provides compression for weak or overstressed knees
- Exclusive reinforced patella stabiliser to prevent displacement
- Retains body heat to warm knee joint and increases blood circulation



Oppo 1124 Knee Support Short

- Wrap - around design allows for adjustable fit and unrestricted range of motion of the knee
- Open patella design to relieve pressure
- Exclusive reinforced patella stabiliser to prevent displacement
- Heat retention and compression reduces potential for strain or injury

Code	Description	Size	DVA:
FBS301501	Oppo 1024 Universal Knee Support	Universal Size	AR02

Code	Description	Size	DVA:
FBS301502	Oppo 1124 Knee Support Short	Universal Size	AR02



Oppo 1021 Knee Support Open Patella

- Provides compression for weak or overstressed knees
- Open patella design relieves pressure
- Exclusive reinforced patella stabiliser to prevent displacement
- Helpful for post - operative knee problems, strains and arthritis



Oppo 1132 Contour Knee Support

- Unique contour shape to snugly fit the knee joint to prevent slipping due to excessive movement
- The adjustable design provides easy application and maximum support
- Breathable CoolPrene® material lined with COOLMAX® behind the knee for optimal comfort

Code	Description	Size	Knee (cm)	DVA:
FBS301508	Oppo 1021 Knee Support Open Patella	S	33.7 - 35.9	AR02
FBS301508M	Oppo 1021 Knee Support Open Patella	M	35.9 - 38.4	AR02
FBS301508L	Oppo 1021 Knee Support Open Patella	L	38.4 - 43.5	AR02
FBS301508XL	Oppo 1021 Knee Support Open Patella	XL	43.5 - 51.1	AR02

Note: Measure around the mid - section of the knee, while in a standing position.

Code	Description	Size	Knee (cm)	DVA:
FBS301505	Oppo 1132 Contour Knee Support	S	33 - 39.4	AR02
FBS301505M	Oppo 1132 Contour Knee Support	M	39.4 - 43.2	AR02
FBS301505L	Oppo 1132 Contour Knee Support	L	43.2 - 47	AR02
FBS301505XL	Oppo 1132 Contour Knee Support	XL	47 - 50.8	AR02

Note: Measure around the mid - section of the knee, while in a standing position.

General Enquiries

P 1300 133 120 | E sales@aidacare.com.au | W www.aidacare.com.au



McDavid A408 Thermal Knee Wrap

- Multi purpose knee wrap
- Provides warmth, compression, & protection
- Contoured & darted for a better fit
- Retains body warmth & can increase blood circulation
- Fully adjustable



McDavid A409 Thermal Knee Wrap Open Patella

- Multi purpose knee wrap
- Provides warmth, compression, & protection
- Contoured & darted for a better fit
- Retains body warmth & can increase blood circulation
- Fully adjustable

Code		Size	DVA:
Black	Beige		
FBS284780	FBS284780BE	Universal Size	AR02

Code		Size	DVA:
Black	Beige		
FBS284790	FBS284790BE	Universal Size	AR02



Oppo 2029 Open Knee Brace

- Open design for better fit
- For arthritic knee conditions, minor sprains and strains
- Provides firm support and gentle compression over the kneecap
- Improves healing and keeps joints and muscles flexible

Code	Description	Size	Knee (cm)	DVA:
FBS301506	Oppo 2029 Open Knee Brace	S	30.5 - 34.3	AR02
FBS301506M	Oppo 2029 Open Knee Brace	M	34.3 - 38.1	AR02
FBS301506L	Oppo 2029 Open Knee Brace	L	38.1 - 41.9	AR02
FBS301506XL	Oppo 2029 Open Knee Brace	XL	41.9 - 45.7	AR02
FBS301506XXL	Oppo 2029 Open Knee Brace	XXL	45.7 - 49.5	AR02

Note: Measure around the mid - section of the knee, while in a standing position.



Oppo 2120 Knee Wrap

- Provides comfortable and uniform compression
- Made of unique extra - strength elastic and transverse nylon fibers for long - lasting, heavy - duty support
- Wrap may be pre - positioned for controlled compression



Oppo 1125 Knee Support Short Donut

- Adjustable wrap - around design for custom fit
- Open patella design relieves pressure on the kneecap
- Donut shape around the patella for added comfort and support

Code	Description	Size	DVA:
FBS301507	Oppo 2120 Knee Wrap	Universal Size	AR02

Code	Description	Size	DVA:
FBS301503	Oppo 1125 Knee Support Short Donut	Universal Size	AR02



McDavid A402 Thermal Knee Support Open

- Warmth, firm compression & some protection
- Retains body warmth & can increase blood circulation
- Reversible & lined both sides (Heavy duty outer lining for durability)



Rehband 8361 Genu Therma

- Provides warmth and compression, which helps to reduce pain and improves blood circulation
- Elastic and supple made from a warming material
- Anatomic design with soft trimmings to provide excellent comfort and fit

Code		Size	(cm)	DVA:
Black	Beige			
FBS280970	FBS280975	S	30 - 33	AR02
FBS280971	FBS280976	M	33 - 38	AR02
FBS280972	FBS280977	L	38 - 43	AR02
FBS280973	FBS280978	XL	43 - 51	AR02
FBS280974	FBS280979	XXL	51 - 56	AR02

Note: Measure around the knee joint.

Code	Size	15 cm Below Mid-Patella (cm)	15 cm Above Mid-Patella (cm)	DVA:
FBS301509	XXS	29 - 32	36 - 40	AR02
FBS301509XS	XS	32 - 35	40 - 44	AR02
FBS301509S	S	35 - 38	44 - 48	AR02
FBS301509M	M	38 - 41	48 - 52	AR02
FBS301509L	L	41 - 44	52 - 56	AR02
FBS301509XL	XL	44 - 48	56 - 61	AR02
FBS301509XXL	XXL	48 - 51	61 - 67	AR02

Note: Measure around leg, 15cm above, and 15cm below middle of knee.

DVA Contact Details

P 1300 888 052 | F 1300 787 052 | E dva@aidacare.com.au



For NDIS Item Code, refer to page: 500





Rehband 8364 Genu Thermo Open

- Provides warmth and compression, which helps to reduce pain and improves blood circulation
- Supports sensory motor control
- Reduces pressure on the patella
- Anatomic design with soft trimmings to provide excellent comfort and fit
- Can be used for left or right

Code	Size	15 cm Below Mid-Patella (cm)	15 cm Above Mid-Patella (cm)	DVA:
FBS301510	XXS	29 - 32	36 - 40	AR02
FBS301510XS	XS	32 - 35	40 - 44	AR02
FBS301510S	S	35 - 38	44 - 48	AR02
FBS301510M	M	38 - 41	48 - 52	AR02
FBS301510L	L	41 - 44	52 - 56	AR02
FBS301510XL	XL	44 - 48	56 - 61	AR02
FBS301510XXL	XXL	48 - 51	61 - 67	AR02

Note: Measure around leg, 15cm above, and 15cm below middle of knee.



Rehband 8354 Genu Thermo Fit

- Provides warmth and compression, which helps reduce pain and improve blood circulation
- Supports sensory motor control
- Reduces pressure on the patella
- Straps to ensure optimal fit and comfort
- Anatomic design with soft trimmings to provide excellent comfort and fit
- Two elastic straps, can be used for left or right

Code	Size	15 cm Below Mid-Patella (cm)	15 cm Above Mid-Patella (cm)	DVA:
FBS301511	XXS	29 - 32	36 - 40	AR02
FBS301511XS	XS	32 - 35	40 - 44	AR02
FBS301511S	S	35 - 38	44 - 48	AR02
FBS301511M	M	38 - 41	48 - 52	AR02
FBS301511L	L	41 - 44	52 - 56	AR02
FBS301511XL	XL	44 - 48	56 - 61	AR02
FBS301511XXL	XXL	48 - 51	61 - 67	AR02

Note: Measure around leg, 15cm above, and 15cm below middle of knee.



Oppo 2021 Open Patella Knee Support

- For arthritic knee conditions, minor sprains and strains provides firms support for the maintenance of joint structures with gentle compression over the kneecap

Code	Description	Size	Knee (cm)	DVA:
FBS301512	Oppo 2021 Open Patella Knee Support	S	30.5 - 34.3	AR02
FBS301512M	Oppo 2021 Open Patella Knee Support	M	34.3 - 38.1	AR02
FBS301512L	Oppo 2021 Open Patella Knee Support	L	38.1 - 41.9	AR02
FBS301512XL	Oppo 2021 Open Patella Knee Support	XL	41.9 - 45.7	AR02

Note: Measure around the mid-section of the knee, while in a standing position.



BodyAssist 450 Knee Support

- Basic multi purpose slip-on knee support
- Provides a low level of compression & support
- Helps ease painful, injured knees
- Lightweight ventilated elastic for all day comfort
- Retains some body warmth & can increase blood circulation
- Length 18cm

Code	Description	Size	Knee (cm)	DVA:
FBS284340	Body Assist 450 Knee Support	S	Up to 32	AR02
FBS284350	Body Assist 450 Knee Support	M	32 - 37	AR02
FBS284360	Body Assist 450 Knee Support	L	37 - 42	AR02
FBS284370	Body Assist 450 Knee Support	XL	42 - 47	AR02

Note: Measure around the knee joint.



Oppo 1029 Patella Tendon Strap

- Extra pad provides compression on target area
- Non - slip backing keeps wrap from sliding out of position
- Great for use after an injury or to help prevent injuries
- Reduces vibrations in the femoral groove
- Fits knee circumference 32 - 40.5cm

Code	Description	Size	DVA:
FBS301513	Oppo 1029 Patella Tendon Strap	One size	AR02

Note: Measure around the knee joint.



Oppo 1429 Jumpers Strap W/Pad

- Specially shaped silicone insert offers extra pressure on the patella
- Adjustable hook and loop strap allows you to apply greater tension if needed
- Minimises potential for future injuries
- Fits knee circumference 32 - 40.5cm

Code	Description	Size	DVA:
FBS301515	Oppo 1429 Jumpers Strap W/Pad	One size	AR02

Note: Measure around the knee joint.



McDavid A414 Patella Knee Strap

- A semi-tubular buttress in a thermal strap with velcro closure
- For compression and pressure on the patellar tendon



BodyAssist PSB Patella Stabilising Brace

- Designed according to knee cap taping technique, locking the knee cap
- Anatomically shaped to provide orthopaedic support
- Retains some body warmth & can increase blood circulation
- A great brace for anyone who needs to tape their knee cap

Code		Size	DVA:
Black	Beige		
FBS301516	FBS301517	Universal Size	AR02

Code	Size	Knee (cm)	DVA:
FBS281440	Child	Up to 32	AR02
FBS285270	Regular	32 - 42	AR02
FBS285280	X Large	42+	AR02

Note: Measure around the knee joint.



Rehband 8360N Genu Carezza Patella Stabiliser

- Keeps patella in correct position
- Prevents patella from luxation
- Supports sensory motor control
- Flexible patella stabiliser
- C-shaped pad, which is adjusted by two separate straps
- Innovative material that combines compression and breathability
- Excellent fit underneath clothes

Code	Size	15 cm Below Mid-Patella (cm)	15 cm Above Mid-Patella (cm)
FBS301520	XXS	29 - 32	36 - 40
FBS301520XS	XS	32 - 35	40 - 44
FBS301520S	S	35 - 38	44 - 48
FBS301520M	M	38 - 41	48 - 52
FBS301520L	L	41 - 44	52 - 56
FBS301520XL	XL	44 - 48	56 - 61
FBS301520XXL	XXL	48 - 51	61 - 67

Note: Measure around leg, 15cm above, and 15cm below middle of knee.



BodyAssist 43C Elastic Cartilage Knee Brace

- Adjustable front closure, lateral/medial condyle pads & patella cutout
- Firm elastic orthopaedic soft tissue support with additional compression over the medial and lateral meniscus (cartilage)
- Three medial and three lateral spring steel stays provide greater support

Code	Size	Knee (cm)
FBS284255	S	Up to 32
FBS284256	M	32 - 37
FBS284257	L	37 - 42
FBS280765	XL	42 - 47
FBS280766	XXL	47+

Note: Measure around the knee joint.



Oppo 1030 Knee Stabiliser

- Neoprene circle pad around the patella provides padded protection and stabilisation to prevent patella displacement
- Two medial and lateral stays provide additional support and stability for the knee
- Helps relieve pain from tendonitis and arthritis
- Retains body heat to warm knee joint and increases blood circulation

- Two hook and loop straps insure comfortable fit and stabilisation during activities

Code	Size	Knee (cm)	DVA:
FBS301631	S	33.7 - 37.5	AR02
FBS301631M	M	37.5 - 40.6	AR02
FBS301631L	L	40.6 - 45.1	AR02
FBS301631XL	XL	45.1 - 50.2	AR02

Note: Measure around the knee joint.



Oppo 1028 Patella Tracker Support

- Prevents patella pain and weakness due to instability
- Helps maintain proper patella position
- Helps prevent Chondromalacia Patella
- The pad design ensures comfortable fit and provides even pressure to the side of the knee
- Does not interfere with normal motion, even when straps are crossed

Code	Size	Knee (cm)	DVA:
FBS301519	XS	30.5 - 33.0	AR02
FBS301519S	S	33.0 - 35.6	AR02
FBS301519M	M	35.6 - 38.1	AR02
FBS301519L	L	38.1 - 40.6	AR02
FBS301519XL	XL	40.6 - 43.2	AR02
FBS301519XXL	XXL	43.2 - 45.5	AR02



Oppo 1230 Adjustable Knee Stabiliser

- Two medial and lateral spiral stays provide additional support and stability for the knee
- The additional circle pad around the patella provides padded protection and stabilisation to prevent patella displacement
- Open patella design to relieve pressure
- Heat retention and compression reduces potential for strain or injury

Code	Description	Size	DVA:
FBS301524	Oppo 1230 Adjustable Knee Stabiliser	One size	AR02



Oppo 1130 Knee Stabiliser Adjustable

- Two medial and lateral spiral stays provide additional support and stability for the knee
- Heat retention and compression reduces potential for strain or injury
- Adjustable for custom fit
- Fits knee circumference 32 - 40.5cm

Code	Description	Size	DVA:
FBS301525	Oppo 1130 Knee Stabiliser	One Size	AR02



Oppo 1033 Knee Support Patella Stabiliser

- Neoprene circle pad around the patella provides protection and stabilisation to prevent patella displacement
- Four neoprene straps provide support directly to patella area
- Two medial and lateral spiral stays provide additional support and stability for the knee
- Helps relieve pain from tendinitis and arthritis
- Retains body heat to warm knee joint and increases blood circulation

Code	Size	Knee (cm)	DVA:
FBS301521	XXS	33.7 - 37.5	AR02
FBS301521M	XS	37.5 - 40.6	AR02
FBS301521L	S	40.6 - 45.1	AR02
FBS301521XL	M	45.1 - 50.2	AR02

Note: Measure around the mid-section of the knee, while in a standing position.



Rehband 8362N Genu Carezza Knee Support

- Support and stabilisation by coil springs
- Supports sensory motor control
- Wraparound design for easy application
- Lateral and medial coil springs for stability and to prevent the brace from rolling
- Anatomic design with soft trimmings to provide excellent comfort and fit
- Innovative material that combines compression and breathability

Code	Size	15 cm Below Mid-Patella (cm)	15 cm Above Mid-Patella (cm)	DVA:
FBS301526	XXS	29 - 32	36 - 40	AR02
FBS301526XS	XS	32 - 35	40 - 44	AR02
FBS301526S	S	35 - 38	44 - 48	AR02
FBS301526M	M	38 - 41	48 - 52	AR02
FBS301526L	L	41 - 44	52 - 56	AR02
FBS301526XL	XL	44 - 48	56 - 61	AR02

Note: Measure around leg, 15cm above, and 15cm below middle of knee.



BodyAssist 41S Cross - Cut Elastic Knee Brace

- Anatomically shaped by a cross-cut design that minimises irritation & pressure behind the knee
- Provides orthopaedic support for minor knee joint strains & sprains
- Removable medical support spring stays on sides of the knee joint to provide medio-lateral support
- Brace retains some body warmth & can increase blood circulation

- Lightweight, ventilated elastic for all day comfort

Code	Size	Knee (cm)	DVA:
FBS284180	S	Up to 32	AR02
FBS284190	M	32 - 37	AR02
FBS284200	L	37 - 42	AR02
FBS284210	XL	42 - 47	AR02
FBS284215	XXL	47+	AR02

Note: Measure around the knee joint.



BodyAssist 42S Knee Support

- An open patella with medial & lateral medical side stays for extra overall support
- Open Patella minimises pressure over the knee cap
- Quality lightweight elastic orthopaedic soft tissue support with additional side supports

Code	Size	Knee (cm)	DVA:
FBS284260	S	30.5 - 35.5	AR02
FBS284270	M	33 - 38	AR02
FBS284280	L	35.5 - 40.5	AR02
FBS284290	XL	40.5 - 46	AR02
FBS284295	XXL	46 +	AR02

Note: Measure around the knee joint.



McDavid A421 Patella Knee Support

- Delivers firm compression & soft tissue support
- Made of medical grade lightweight neoprene that retains body warmth
- Contoured design for a better fit with two medial and two lateral medical steel stays
- Top and bottom anchor straps to minimise slippage

- Available in black or beige

Code		Size	Knee (cm)	DVA:
Black	Beige			
FBS284820	FBS284852	S	30 - 33	AR02
FBS284830	FBS284853	M	33 - 38	AR02
FBS284840	FBS284854	L	38 - 43	AR02
FBS284850	FBS284855	XL	43 - 51	AR02
FBS284851	FBS284856	XXL	51 - 56	AR02

Note: Measure around the knee joint.



McDavid A425 Ligament Knee Support

- Delivering firm compression and soft tissue support
- Made of medical grade lightweight neoprene that retains body warmth
- 4 elastic straps that deliver additional support to ligaments
- Horse-shoe buttress to stabilise the knee cap
- Contoured design for a better fit with two medial and two lateral medical steel stays

- Top and bottom anchor straps minimise slippage

Code			Size	Knee (cm)
Black	DVA:	Beige		
FBS284860	AR02	FBS284890	S	30 - 33
FBS284870	AR02	FBS284895	M	33 - 38
FBS284880	AR02	FBS284900	L	38 - 43
FBS284881	AR02	FBS284901	XL	43 - 51
FBS284882	AR02	FBS284902	XXL	51 - 56

Note: Measure around the knee joint.



Rehband 8353N Genu Direxa Knee Support

- Sensorimotor performance improvements through controlled compression
- Pain relief and functional improvement of the knee joint
- Polycentric joints with extension and flexion stops of the full range of motion
- Innovative material that combines compression and breathability
- For patients with restricted grip force and mobility

Code	Size	15 cm Below Mid-Patella	15 cm Above Mid-Patella	DVA:
FBS301527	XXS	29 - 32	36 - 40	AR02
FBS301527XS	XS	32 - 35	40 - 44	AR02
FBS301527S	S	35 - 38	44 - 48	AR02
FBS301527M	M	38 - 41	48 - 52	AR02
FBS301527L	L	41 - 44	52 - 56	AR02
FBS301527XL	XL	44 - 48	56 - 61	AR02
FBS301527XXL	XXL	48 - 51	61 - 67	AR02

Note: Measure around leg, 15cm above, and 15cm below middle of knee.



Oppo 1032 Knee Support Post Operative

- Reduces swelling and pain following arthroscopic surgery and enhances protection of periarticular structures
- Provides comfortable and controlled compression with adjustable hook and loop closures
- Front opening design for easy fit and comfort

Code	Size	Knee (cm)	DVA:
FBS301634S	S	33.7 - 37.5	AR02
FBS301634M	M	37.5 - 40.6	AR02
FBS301634L	L	40.6 - 45.1	AR02
FBS301634XL	XL	45.1 - 50.2	AR02

Note: Measure around the mid-section of the knee, while in a standing position.



BodyAssist 43H Hinged Knee Brace

- Rigid hinges minimise excessive medio - lateral movement of the knee joint, while offering excellent support & stability
- Made from lightweight medical grade orthopaedic elastic & features front closures for easier application
- Top & bottom straps lock the brace in position & minimise movement

Code	Size	Knee (cm)	DVA:
FBS280770	S	Up to 32	AR02
FBS284300	M	32 - 37	AR02
FBS284310	L	37 - 42	AR02
FBS284315	XL	42 - 47	AR02
FBS284220	XXL	47+	AR02

Note: Measure around the knee joint.



Oppo 1031 Knee Stabiliser Hinged

- Neoprene circle pad around the patella provides padded protection and stabilisation to prevent patella displacement
- Two medial and lateral stays provide additional support and stability for the knee
- Two hook and loop straps insure comfortable fit and stabilisation during activities

- Helps relieve pain from tendinitis and arthritis
- Retains body heat to warm knee joint and increases blood circulation

Code	Size	Knee (cm)	DVA:
FBS301528	S	33.7 - 37.5	AR02
FBS301528M	M	37.5 - 40.6	AR02
FBS301528L	L	40.6 - 45.1	AR02
FBS301528XL	XL	45.1 - 50.2	AR02

Note: Measure around the mid-section of the knee, while in a standing position.



McDavid A429 PSII Hinged Knee

- Hinged knee with moulded condyle pads that offer superior fit, comfort & protection
- Professional quality bilateral design
- Geared rigid upright polycentric hinges provide additional medial-lateral support
- Features superior horseshoe buttress & hyper-extension stops

- Top & bottom velcro strap locks eliminate brace movement

Code	Size	Knee (cm)	DVA:
FBS301529	S	30 - 33	AR02
FBS3015292	M	33 - 38	AR02
FBS3015293	L	38 - 43	AR02
FBS3015294	XL	43 - 51	AR02
FBS3015295	XXL	51 - 56	AR02

Note: Measure around the knee joint.



Oppo 2031 Knee Open Patella Support

- May be used for stabilisation and support of the knee following strains, injury or surgery
- Elastic brace with motion hinges provides medial and lateral support for the unstable knee
- Provides firm support for the maintenance of joint structures with gentle compression over the kneecap

- For arthritic knee conditions, minor sprains and strains
- Two hook and loop straps insure comfortable fit

Code	Size	Knee (cm)	DVA:
FBS301530	S	30.5 - 34.3	AR02
FBS301530M	M	34.3 - 38.1	AR02
FBS301530L	L	38.1 - 41.9	AR02
FBS301530XL	XL	41.9 - 45.7	AR02

Note: Measure around the mid-section of the knee, while in a standing position.



Rehband 8367N Direxa Stable Knee Support

- Sensorimotor performance improvements through controlled compression
- Pain relief and functional improvement of the knee joint
- Polycentric joints with gradual extension and flexion stops (extension 0°; 10°; 20°; 30°; 45°) (flexion 0°; 10°; 20°; 30°; 45°; 60°; 75°; 90°)
- Two non-elastic, circular strap systems

Code	Size	15 cm Below Mid-Patella (cm)	15 cm Above Mid-Patella (cm)	DVA:
FBS301532	XXS	29 - 32	36 - 40	AR02
FBS301532XS	XS	32 - 35	40 - 44	AR02
FBS301532S	S	35 - 38	44 - 48	AR02
FBS301532M	M	38 - 41	48 - 52	AR02
FBS301532L	L	41 - 44	52 - 56	AR02
FBS301532XL	XL	48 - 48	56 - 61	AR02
FBS301532XXL	XXL	48 - 51	61 - 67	AR02

Note: Measure around leg, 15cm above, and 15cm below middle of knee.



Rehband 8368N Genu Direxa Stable ROM Long

- Sensorimotor performance improvements through controlled compression
- Pain relief and functional improvement of the knee joint
- Four-strap stabilisation
- Polycentric joints with gradual extension and flexion stops (extension 0°; 10°; 20°; 30°; 45°) (flexion 0°; 10°; 20°; 30°; 45°; 60°; 75°; 90°)

Code		Size	15 cm Below Mid-Patella (cm)	15 cm Above Mid-Patella (cm)	DVA:
Left	Right				
FBS301533LXXS	FBS301533	XXS	29 - 32	36 - 40	AR02
FBS301533LXS	FBS301533RXS	XS	32 - 35	40 - 44	AR02
FBS301533LS	FBS301533RS	S	35 - 38	44 - 48	AR02
FBS301533LM	FBS301533RM	M	38 - 41	48 - 52	AR02
FBS301533LL	FBS301533RL	L	41 - 44	52 - 56	AR02
FBS301533LXL	FBS301533RXL	XL	48 - 48	56 - 61	AR02
FBS301533LXXL	FBS301533RXXL	XXL	48 - 51	61 - 67	AR02

Note: Measure around leg, 15cm above, and 15cm below middle of knee.



Oppo 1036 Knee Support Multi Orthosis

- Functional 3/16" breathable neoprene knee brace
- Polycentric hinge with flexion and extension pin-stops for controlled range of motion
- Increases stability and provides support and protection
- Provides greatest comfort even in the hollow of the knee
- For patients with restricted grip force and mobility
- Lightweight, open patella design

Code	Size	Knee (cm)	DVA:
FBS301534	S	33.7 - 37.5	AR02
FBS301534M	M	37.5 - 40.6	AR02
FBS301534L	L	40.6 - 45.1	AR02
FBS301534XL	XL	45.1 - 50.2	AR02

Note: Measure around the mid-section of the knee, while in a standing position.



Oppo 4039 Post-Op Motion Control Knee Splint

- Provides stability to an injured or unstable knee by allowing normal knee joint motion with proper settings or stops with in a selected range
- For use following knee surgery, early cast removal or post-operative rehabilitation
- Settings from 0 to 120° flexion and between 0 to 75° extension in 7.5° increments
- Hook and loop closures for uniform compression and support
- Available in three lengths: 46cm, 56cm, 66cm
- Fits knee circumference 43 - 71 cm

Code			Knee (cm)	DVA:
46cm	56cm	66cm		
FBS301536	FBS30153622	FBS30153626	43 - 71cm	AR02

Note: Measure around the knee joint.



Oppo Post-Op Motion Control Cool Knee Splint

- Precision flexion and extension settings are achievable in a single motion—simply press, rotate and release
- Features a lightweight design, one size fits all
- Designed for greater breathability and patient comfort

Code	Description	Size	DVA:
FBS301537	Post-Op Motion Control Cool Knee Splint	Universal Size	AR02



Oppo Tri-Panel Knee Immobiliser

- Universal knee immobiliser
- 2 side panels velcro on to the middle panel to allow a range of circumferential sizing
- 3 panel design, length: 51cm

Code	Description	Size
FBS301539	Tri-Panel Knee Immobiliser	Universal Size



Oppo 4030 Knee Immobiliser

- Provides varying degrees of support and stabilisation for the leg and knee
- Aluminium posterior stays retain proper position for comfort
- Hook and loop closure provides quick, secure fit
- Foam padding reduces slippage and will not fray
- Available in three heights: 46cm, 51cm and 58cm

Code			Size	Thigh (cm)
46 cm	51 cm	58 cm		
FBS301538	FBS301538S51	FBS301538S58	S	33 - 38.1
FBS301538M46	FBS301538M51	FBS301538M58	M	38.1 - 45.7
FBS301538L46	FBS301538L51	FBS301538L58	L	45.7 - 53.3
FBS301538XL46	FBS301538XL51	FBS301538XL58	XL	53.3 - 61


Note: Measure the circumference of the thigh approximately 10cm above the knee.



Post-Op Knee Immobiliser

- Removable medial, lateral and posterior aluminium stays can be contoured to the shape of the leg for a comfortable and secure fit
- Available in 3 lengths: 40cm, 48cm, 58cm
- The 3 panel design (including the popliteal variation) allows adjustment of the brace
- Will fit calf circumference up to 46.5cm and thigh circumference up to 56.5cm

Code			Size
40 cm	48 cm	58 cm	
FBS301540	FBS301541	FBS301542	Universal Size



BodyAssist 480W Thigh Tendon Strap

- Rigid tendon strap that delivers firm support and assists in minimising further injury to the thigh
- Designed by a physiotherapist to assist in chronic thigh injury as a rehabilitative device
- Suitable for use in various sports

Code	Description	Size
FBS301543	480W Thigh Tendon Strap	Universal Size



Rotulax Elastic Knee - Closed Patella

- Provides support for patients with hyperlax patella
- Post-operative support after patellar operations
- Fits left or right patella

Code	Size	Thigh (cm)	DVA:
FBS303527	XS	35 - 38	AR02
FBS303528	S	39 - 42	AR02
FBS303529	M	43 - 46	AR02
FBS303530	L	47 - 49	AR02
FBS303531	XL	50 - 53	AR02
FBS303532	XXL	54 - 57	AR02

Measurement: At the thigh, 15cm above the kneecap



Rotulax Elastic Knee - Open Patella

- General stabilisation of acute patella injury
- Preventative use in case of hyperlax patella
- Fits left or right patella

Code	Size	Thigh (cm)
FBS303533	XS	35 - 38
FBS303534	S	39 - 42
FBS303535	M	43 - 46
FBS303536	L	47 - 49
FBS303537	XL	50 - 53
FBS303538	XXL	54 - 57

Measurement: At the thigh, 15cm above the kneecap



Donjoy OA Reaction Web Brace

- A unique and effective way to treat osteoarthritis knee pain utilising patented silicone web technology
- Breathable, open framework
- Low profile design - wear it discreetly under clothing
- Includes undersleeve

Code		Size	Thigh (cm)	DVA:
Right Medial/Left Lateral	Left Medial/Right Lateral			
FBS303444	FBS303451	XS	33 - 39	AR02
FBS303445	FBS303452	S	39 - 47	AR02
FBS303446	FBS303453	M	47 - 53	AR02
FBS303447	FBS303454	L	53 - 60	AR02
FBS303448	FBS303455	XL	60 - 67	AR02
FBS303449	FBS303456	XXL	67 - 75	AR02
FBS303450	FBS303457	XXXL	75 - 81	AR02

Measurement: At the thigh, 15cm above the kneecap

Bracing & Supports - Thigh & Calf

I ♥ ndis For NDIS Item Code, refer to page: 500



BodyAssist 490 Elastic Slip-On Thigh Support

- General multi purpose slip-on elastic thigh support
- Lightweight ventilated cotton elastic for improved comfort
- Ideal for holding a hot or cold pack
- Retains some body warmth & can increase blood circulation

Code	Description	Size	Thigh (cm)	DVA:
FBS284380	490 Elastic Slip-on Thigh Support	S	Up to 44	AR02
FBS284390	490 Elastic Slip-on Thigh Support	M	44 - 50	AR02
FBS284400	490 Elastic Slip-on Thigh Support	L	50 - 56	AR02
FBS284410	490 Elastic Slip-on Thigh Support	XL	56+	AR02

Note: Measure around mid-section of the thigh.



Oppo 1040 Slip-On Thigh Support

- Provides comfortable compression for weak or overstressed thigh and hamstring muscles
- Helps reduce pain from tendonitis
- Retains body heat to warm the thigh and increase blood circulation
- Minimises the potential for future injuries

Code	Description	Size	Thigh (cm)	DVA:
FBS301544	1040 Slip-on Thigh Support	S	45.1 - 50.8	AR02
FBS301544M	1040 Slip-on Thigh Support	M	50.8 - 57.2	AR02
FBS301544L	1040 Slip-on Thigh Support	L	57.2 - 62.2	AR02
FBS301544XL	1040 Slip-on Thigh Support	XL	62.2 - 68.9	AR02

Note: Measure around mid-section of the thigh.



McDavid A478 Thermal Thigh Wrap

- Provides warmth, firm compression for strains & hamstring pulls
- For minor protection & shock absorption & helps reduce muscle strain
- Simple application with Velcro closure
- One size fits most
- Retains body warmth & can increase blood circulation
- Ideal for corked thigh condition

Code	Description	Size	DVA:
FBS281375	McDavid A478 Thermal Thigh Wrap	Universal Size	AR02

DVA Contact Details

P 1300 888 052 | F 1300 787 052 | E dva@aidacare.com.au

I ♥ ndis For NDIS Item Code, refer to page: 500

aidacare
HEALTHCARE EQUIPMENT



Oppo 2420 Silicone Wrap

- Provides localised support and compression, and avoids unnecessary compression to areas where it is not required
- Tension is fully adjustable to provide from light support to total compression
- Can be used for thigh, hamstring, knee or calf
- Cohesive fabric, causing no pain or irritation to the skin whilst in use, comfortable to wear and painless to remove
- 7.5cm width

Code	Description	Size	DVA:
FBS301546	2420 Silicone Wrap	Universal Size	AR02



BodyAssist N478 Thermal Thigh Wrap

- Provides warmth, firm compression for strains & hamstring pulls
- For minor protection & shock absorption & helps reduce muscle strain
- Simple application with Velcro closure
- One size fits most
- Ideal for corked thigh condition
- Retains body warmth & can increase blood circulation

Code	Description	Size	DVA:
FBS281380	N478 Thermal Thigh Wrap	Universal Size	AR02



Oppo 1010 Slip-On Shin Support

- Designed to fit the upper portion of the tibia and calf
- Provides compression and support over the calf and shin
- Retains body heat to warm shin and increase blood circulation
- Minimises the potential for future injuries

Code	Description	Size	Shin (cm)	DVA:
FBS301547	1010 Slip-on Shin Support	S	33.0 - 36.2	AR02
FBS301547M	1010 Slip-on Shin Support	M	36.2 - 39.4	AR02
FBS301547L	1010 Slip-on Shin Support	L	39.4 - 42.5	AR02
FBS301547XL	1010 Slip-on Thigh Support	XL	42.5 - 46.4	AR02

Note: Measure around mid-section of the shin.



Bracing & Supports - Ankle

I ♥ ndis For NDIS Item Code, refer to page: 500



McDavid A431 Slip-On Thermal Ankle Support

- A general purpose slip-on neoprene ankle support
- Provides all the benefits of thermal compression

Code	Size	Shoe Size
FBS301556	S	37 - 40
FBS3015562	M	40.5 - 42.5
FBS3015563	L	43 - 45
FBS3015564	XL	45.5+

Note: Sizing by European shoe size.



McDavid A432 Thermal Ankle Support W/Straps

- Thermal slip-on ankle support with 2 orthopaedic elasticised adjustment straps for graduated support
- Provides thermal compression with additional outer support
- Lined both sides
- Outer lining is made of heavy duty nylon for greater durability
- Fits left or right
- Protects and supports without restriction

Code	Size	Shoe Size
FBS301557	S	37 - 40
FBS3015572	M	40.5 - 42.5
FBS3015573	L	43 - 45
FBS3015574	XL	45.5+

Note: Sizing by European shoe size.



McDavid A430B Open Thermal Ankle Support

- Fully adjustable thermal ankle support with 2 orthopaedic elastic adjustable straps for graduated support, lined both sides
- Outer lining is made of heavy duty nylon for greater durability
- One size fits most
- Provides all the benefits of thermal compression together with added outer support layer
- Protects and supports without restriction and fits left or right

Code	Description	Size	DVA:
FBS301550	A430B Open Thermal Ankle Support	Universal Size	AR01



McDavid A438 Thermal Ankle Wrap

- One size thermal ankle wrap providing the benefits of heat & compression
- Adjustable closures ensure an optimum fit & easy application & removal
- Allows increased blood circulation
- Fits easily into standard footwear
- Compresses and supports soft tissue while retaining therapeutic heat
- Lined both sides, one size fits most

Code	Description	Size	DVA:
FBS284965	A438 Thermal Ankle Wrap - Beige	Universal Size	AR01



BodyAssist 01A10 Elastic Ankle Wrap With Heel Loop

- This one size elastic ankle wrap with loop anchor is simple to apply & very effective
- Simply hook the loop anchor over the heel & start wrapping
- The adjustable touch tape closure assures a perfect fit and locks down anywhere over itself
- Will not stick to skin
- Pull hard & the strap locks rigid
- Orthopaedic elastic that locks like tape
- Fits easily into standard footwear

Code		Size
Black	Beige	
FBS301552B	FBS301552S	Universal Size



Oppo 1001 Slip-On Thermal Ankle Support

- Open toe and heel design allows full range of movement
- Relieves pain due to strains and stiffness
- Provides compression for weak or overstressed ankles
- Helps ankle heal faster following injury

Code	Description	Size	Ankle (cm)
FBS301553	1001 Slip-on Thermal Ankle Support	S	15.6 - 20.6
FBS301553M	1001 Slip-on Thermal Ankle Support	M	20.6 - 25.7
FBS301553L	1001 Slip-on Thermal Ankle Support	L	25.7 - 30.8
FBS301553XL	1001 Slip-on Thermal Ankle Support	XL	30.8 - 35.9

Note: Measure circumference just above ankle bone.



Oppo 2004 Slip-On Adj Elastic Ankle Support

- Ideal for weak or unstable ankles, sprains, strains, or twisting injuries
- Provides firm support to the ankle joint with gentle compression
- Increases circulation to relieve pain, speed healing and maintain flexibility
- 35% cotton blend

Code	Size	Ankle (cm)	DVA:
FBS301558	S	15.6 - 21.0	AR01
FBS301558M	M	21.5 - 26.0	AR01
FBS301558L	L	26.5 - 31.0	AR01
FBS301558XL	XL	31.1 - 36.2	AR01

Note: Measure circumference just above ankle bone.



Oppo 1003 Slip-On Thermal Ankle Support w/strap

- Removable surgical elastic strap allows adjustable compression and support for ankle joint
- Open toe and heel design allows full range of movement
- Relieves pain due to strains or stiffness
- Provides compression for weak or overstressed ankle

Code	Size	Ankle (cm)	DVA:
FBS301554	S	15.6 - 20.6	AR01
FBS301554M	M	20.6 - 25.7	AR01
FBS301554L	L	25.7 - 30.8	AR01
FBS301554XL	XL	30.8 - 35.9	AR01

Note: Measure circumference just above ankle bone.



Oppo 1008 Slip-On Thermal Ankle Support

- Open heel design allows for comfortable fit and applies less pressure on the sensitive Achilles tendon area
- Detachable strap may be removed when additional compression is not needed
- Flexible spiral springs act like tendons to reinforce the ankle
- Relieves pain due to strains or stiffness
- Provides comfortable compression for weak or overstressed ankles

Code	Size	Ankle (cm)	DVA:
FBS301555	S	15.6 - 20.6	AR01
FBS301555M	M	20.6 - 25.7	AR01
FBS301555L	L	25.7 - 30.8	AR01
FBS301555XL	XL	30.8 - 35.9	AR01

Note: Measure circumference just above ankle bone.



Oppo 2101 Elastic Ankle Wrap

- Provides comfortable and uniform compression
- Made of unique extra-strength elastic and transverse nylon fibers for long-lasting, heavy-duty support
- Wrap may be pre-positioned for controlled compression

Code	Description	Size	DVA:
FBS301551	2101 Elastic Ankle Wrap	Universal	AR01



Oppo 2204 Slip-On Adj Elastic Ankle Support

- Provides support and relieves pain from strains and sprains
- Helps minimise further injury

Code	Size	Ankle (cm)	DVA:
FBS301560	S	15.6 - 21.0	AR01
FBS301560M	M	21.5 - 26.0	AR01
FBS301560L	L	26.5 - 31.0	AR01

Note: Measure circumference just above ankle bone.



McDavid A511 Slip-On Elastic Ankle Support

- Tight weave basic tubular support
- Provides low level support & delivers protection to painful, tired ankles
- Retains some body warmth & can increase blood circulation

Code	Size	Shoe Size	DVA:
FBS301568	S	37 - 40	AR01
FBS3015682	M	40.5 - 42.5	AR01
FBS3015683	L	43 - 45	AR01
FBS3015684	XL	45.5+	AR01

Note: Sizing by European shoe size.



Oppo Air Lite 4009 Ankle Brace

- Two rigid plastic shells have dual air layers over cushioning foam
- Air is allowed to flow within the air bladder to help reduce swelling
- Heel strap is adjustable to fit a range of sizes, and two adjustable swivel straps ensure proper placement each time
- Limits inversion/eversion while permitting near-normal flexion
- The brace designed for active people resuming activity from an acute ankle sprain

Code	Description	Size	Shoe Size
FBS301569	Air Lite 4009 Ankle Brace Medium	M	40 - 41
FBS301569L	Air Lite 4009 Ankle Brace Large	L	42 - 44

Note: Sizing by European shoe size.



Fortilax Elastic Ankle

- Support and compression following acute and chronic ankle injury and mild ankle and foot osteoarthritis
- Proprioceptive and neuromuscular stimulation through compression
- Fits either left or right ankle

Code	Size	Ankle (cm)	DVA:
FBS303545	XS	19 - 20	AR01
FBS303546	S	21 - 22	AR01
FBS303547	M	23 - 24	AR01
FBS303548	L	25 - 26	AR01
FBS303549	XL	27 - 28	AR01
FBS303550	XXL	28 - 29	AR01

Measurement: Ankle circumference above the malleoli



Strapping Elastic Ankle

- Support and compression of mild ankle and foot osteoarthritis
- Also suitable following acute and chronic ankle injury and instability
- Protective and prophylactic use for sports activities
- Fits either left or right ankle

Code	Size	Ankle (cm)
FBS303539	XS	19 - 20
FBS303540	S	21 - 22
FBS303541	M	23 - 24
FBS303542	L	25 - 26
FBS303543	XL	27 - 28
FBS303544	XXL	28 - 29

Measurement: Ankle circumference above the malleoli



Stromgren S355 Stirrup Lock Ankle Support

- Designed to conform to the ankle anatomy like tape, without the pain
- Provides solid support incorporating two heel lock & two forefront lock straps
- The straps extend until they are rigid, therefore achieving a taped ankle effect without tape
- The support is locked using Velcro closure enabling easy adjustment
- Restricts inversion & eversion injuries
- No break-in time needed
- Not bulky in a shoe & washable

Code		Size	Shoe Size	DVA:
Black	White			
FBS301561	FBS301561W1	S	37 - 40	AR01
FBS301561B2	FBS301561W2	M	40.5 - 42.5	AR01
FBS301561B3	FBS301561W3	L	43 - 45	AR01
FBS301561B4	FBS301561W4	XL	45.5+	AR01

Note: Sizing by European shoe size.



BodyAssist 300 Elastic Ankle Brace

- A general purpose adjustable spiral lift elastic ankle brace
- Provides low level orthopaedic support & protection to painful, tired ankles
- Includes an adjustable closure to assist in simple application & allows variation to compression levels
- Retains some body warmth & can increase blood circulation
- Fits easily into standard footwear

Code	Size	(cm)	DVA:
FBS284080	S	Up to 45	AR01
FBS284090	M	45 - 50	AR01
FBS284100	L	50 - 55	AR01
FBS284110	XL	55+	AR01

Note: Measure Figure 8 circumference per illustration.



BodyAssist 310 Slip-On Ankle Support

- A general purpose spiral lift elastic ankle brace
- Provides low level orthopaedic support & protection to painful, tired ankles
- Retains some body warmth & can increase blood circulation
- Fits easily into standard footwear

Code	Size	(cm)
FBS301563	S	Up to 45
FBS301564	M	45 - 50
FBS301565	L	50 - 55
FBS301566	XL	50+

Note: Measure Figure 8 circumference per illustration.



BodyAssist 35 Slip-On Basic Ankle Support

- A general purpose, tubular design support
- Provides low level support & some protection to painful, tired ankles
- Retains body warmth & can increase blood circulation around the ankle
- Fits easily into standard footwear

Code		Size	(cm)
Black	Beige		
FBS301567	FBS301567B1	S	37 - 40
FBS30156702	FBS301567B2	M	40.5 - 42.5
FBS30156703	FBS301567B3	L	43 - 45
FBS30156704	FBS301567B4	XL	45.5+

Note: Measure Figure 8 circumference per illustration.



Oppo 4005 Total Stability Ankle Brace

- Unique heel-lock strapping system
- Non-stretch figure 8 straps provide firm stability to an injured ankle, preventing inversion and eversion sprains
- Removable bilateral plastic splints for additional support
- Padded heel and tongue fabric is durable and provides added comfort

Code	Description	Size	Shoe Size	DVA:
FBS301572	4005 Total Stability Ankle Brace	S	40 - 42	AR01
FBS301572M	4005 Total Stability Ankle Brace	M	43 - 45	AR01
FBS301572L	4005 Total Stability Ankle Brace	L	46 - 48	AR01
FBS301572XL	4005 Total Stability Ankle Brace	XL	49 - 51	AR01

Note: Sizing by European shoe size.



McDavid A195 Ultralight Ankle Brace

- Ultra light nylon laceup with fully adjustable figure 8 straps
- Conforms to & supports sprained or weak ankles
- Stirrup straps restrict inversion & eversion injury & simulate costly tape
- Straps fully adjustable without removing shoes or re-lacing
- Elastic heel delivers Achilles tendon support

Code		Size	Shoe Size	DVA:
Black	White			
FBS280901	FBS280906	XS	37.5 - 38.5	AR01
FBS280902	FBS280907	S	39 - 40	AR01
FBS280903	FBS280908	M	40.5 - 42.5	AR01
FBS280904	FBS280909	L	43 - 45	AR01
FBS280905	FBS280910	XL	45.5+	AR01

Note: Sizing by European shoe size.



McDavid A101 Ankle Guard

- Suitable for those returning from injury or requiring on-going support
- Medical inserts can be removed if desired, although in place they provide additional support & protection
- Replaces costly tape & supports sprained or weak ankles
- Sewn in spring steel medical stays & an elastic heel inset for achilles tendon support
- Fits right or left foot & fully washable

- Notched front allows more free plantar & dorsi flexion

Code		Size	Shoe Size	DVA:
Black	White			
FBS284677	FBS284687	XXS	Up to 37	AR01
FBS284678	FBS284688	XS	37.5 - 38.5	AR01
FBS284679	FBS284689	S	39 - 40	AR01
FBS284680	FBS284690	M	40.5 - 42.5	AR01
FBS284685	FBS284691	L	43 - 45	AR01
FBS284686	FBS284692	XL	45.5+	AR01

Note: Sizing by European shoe size.



McDavid A199 Lightweight Laced Ankle Brace

- This lightweight version of a laced ankle brace is equally designed to conform to the ankle anatomy & to restrict inversion & eversion (rolling of the ankle)
- Includes sewn in spring steel medical stays & consists of two layers of vinyl/nylon material
- Elastic heel achieves excellent results in supporting the achilles tendon
- Unique lining wicks perspiration away

- Notched front allows more free plantar & dorsi flexion movement

Code		Size	Shoe Size	DVA:
Black	White			
FBS280915	FBS280920	XS	37.5 - 38.5	AR01
FBS280916	FBS280921	S	39 - 40	AR01
FBS280917	FBS280922	M	40.5 - 42.5	AR01
FBS280918	FBS280923	L	43 - 45	AR01
FBS280919	FBS280924	XL	45.5+	AR01

Note: Sizing by European shoe size.



McDavid A189 Ankle - X Hinged Ankle Brace

- Ankle brace that offers a high level of protection & comfort for athletes of all ages
- Consists of a neoprene sleeve worn inside the rigid component
- Makes this brace very comfortable
- Pivot Hinges
- Lightweight

- Manoeuverability not restricted
- Very effective for high ankle sprains

Code	Description	Size	Shoe Size	DVA:
FBS301573L/XL	McDavid A189 Ankle	L/XL	45+	AR01

Note: Sizing by European shoe size.



McDavid A188g Ultra Graphite Hinged Ankle Brace

- Ankle support taken to a new level
- Made to the highest professional quality, this hinged ankle brace is designed to treat all ankle injuries (including high ankle sprains)
- Engineered to reduce excessive inversion, eversion and rotation ankle injuries
- Fully Padded & very quick to apply (no laces)
- Ideal for rehab or as injury prevention

- Works well with foot beds & orthotics
- Universal fitting for Left or Right

Code	Size	Shoe Size
FBS301574M	M	41 - 44
FBS301574L	L	45+

Note: Sizing by European shoe size.

Bracing & Supports - Foot & Toe



For NDIS Item Code, refer to page: 500



Darco All Purpose Boot

- Rocker Sole reduces plantar pressure under forefoot and heel by over 25%
- Square Toe Design acts as a bumper to protect the toe
- Strap seats foot firmly in shoe to reduce heel slippage and eliminate friction
- Universal Left/Right fit
- Insole can be totally removed to accommodate DARCO's innovative Peg Assist™ customisable off-loading insole
- Removable Forefoot Closure provides security while eliminating buckle pressure

Size	Shoe Size		Code	DVA:
	Female	Male		
XS	4 - 7	-	FBS280240	AJ06
S	7.5 - 10	6 - 8	FBS280241	AJ06
M	10.5 - 13	8.5 - 10	FBS280242	AJ06
L	13+	10.5 - 12	FBS280245	AJ06
XL	-	12.5 - 14	FBS280246	AJ06

Note: Sizing by US shoe size.



Darco Heelwedge Shoe

- Effective off-loading Footwear
- Reduces Heel pressure by over 26%
- Square Toe Design acts as a bumper to protect the toes and provides better fit
- Ankle Strap seats foot firmly in shoe to reduce heel slippage and eliminate friction
- Zoned Outsole has an aggressive tread under the mid-foot where traction is needed most
- Insole can be totally removed to accommodate DARCO's innovative Peg Assist™ customisable off-loading insole for an even more targeted off-load

Size	Shoe Size		Code	DVA:
	Female	Male		
XS	4 - 7	-	FBS301577	AJ06
S	7.5 - 10	6 - 8	FBS3015771B	AJ06
M	10.5 - 13	8.5 - 10	FBS3015772B	AJ06
L	13+	10.5 - 12	FBS3015773B	AJ06
XL	-	12.5 - 14	FBS3015774B	AJ06

Note: Sizing by US shoe size.



Darco Med Surg Shoe

- Square Toe Design acts as a bumper to provide extra protection and facilitates better left/right fit
- MetaShank™ protection provides more rigid control under the metatarsal heads
- Rocker Sole reduces plantar pressure on the heel and forefoot by over 25%
- Fully compatible with DARCO's Peg-Assist™ off-loading insole

• Female/Male specific for optimal fit

Size	Shoe Size		Code		
	Female	Male	Female	Male	DVA:
S	4 - 6	6 - 8	FBS301579	FBS301578	AJ06
M	6.5 - 8	8.5 - 10	FBS3015792B	FBS3015782B	AJ06
L	8.5 - 10	10.5 - 12	FBS3015793B	FBS3015783B	AJ06
XL	-	12.5 - 14	-	FBS3015784B	AJ06

Note: Sizing by US shoe size.



Darco Orthowedge Shoe

- 15° Wedge Sole effectively shifts body weight to the midfoot and heel
- Reduces forefoot pressure by as much as 57%
- Square Toe Design acts as a bumper to protect the toes when k-wires are present and provides better left/right fit
- Insole can be totally removed to accommodate DARCO's innovative Peg-Assist™ customisable off-loading insole

• Zoned Outsole has an aggressive tread under the mid-foot where traction is needed most

Size	Shoe Size		Code	DVA:
	Female	Male		
XS	4-7	-	FBS301580	AJ06
S	7.5-10	6 - 8	FBS3015801B	AJ06
M	10.5-13	8.5 - 10	FBS3015802B	AJ06
L	13+	10.5 - 12	FBS3015803B	AJ06
XL	-	12.5 - 14	FBS3015804B	AJ06

Note: Sizing by US shoe size.



Darco Soffie Shoe

- Economical, post-op shoe indicated for soft tissue trauma or surgery
- Semi-rigid metatarsal shank provides control of forefoot motion
- Tissue friendly upper material conforms to bony abnormalities and is breathable
- Easy to fasten hook and loop straps
- Female/Male specific for optimal fit
- Fits left & right foot

Size	Shoe Size		Code		
	Female	Male	Female	Male	DVA:
S	4 - 6	6 - 8	FBS280230	FBS280233	AJ06
M	6.5 - 8	8.5 - 10	FBS280231	FBS280234	AJ06
L	8.5 - 10	10.5 - 12	FBS280232	FBS280235	AJ06
XL	-	12.5 - 14	-	FBS280236	AJ06

Note: Sizing by US shoe size.



Cosmac Cast Shoe

- Curved 'rocker' sole for easy ambulation
- Canvas upper with velcro closures at front and back
- Symmetrical shape fits left & right foot

Code	Description	Size	Shoe Length (cm)	DVA:
FBS301583S	Cosmac Cast Shoe	S	23	AJ06
FBS301583M	Cosmac Cast Shoe	M	27	AJ06
FBS301583L	Cosmac Cast Shoe	L	32	AJ06



Darco Body Armor Cast Shoe

- Ultimate cast protection in a new modern design that will appeal to all patients
- Soft but durable EVA provides shock absorption while protecting the cast from the elements
- Unique bungee closure system is easy to use, easy to adjust and prevents the heel from slipping out of the cast shoe

• Soft, durable TPR
• Fits left & right foot

Size	Shoe Size		Code	DVA:
	Female	Male		
XS	Up to 5	-	FBS301582	AJ06
S	5.5 - 8	-	FBS3015821	AJ06
M	8.5 - 10	6 - 8	FBS3015822	AJ06
L	10.5 - 13	8.5 - 11	FBS3015823	AJ06
XL	13.5+	11.5 - 14	FBS3015824	AJ06

Note: Sizing by US shoe size.



Darco Slimline Cast Boot

- Be adapted to fiberglass casts, bulky compression bandages and a variety of wound care applications, making it one of the most versatile cast shoes on the market
- Forefoot Closure can be adjusted to accommodate casts or bandages
- Symmetrical shape fits left & right foot
- Higher Ankle Strap adds stability
- Rocker Sole allows for a smooth roll throughout the gait pattern

Size	Shoe Size		Code	DVA:
	Female	Male		
XS	Up to 5	-	FBS301583	AJ06
S	5.5 - 8	-	FBS3015831B	AJ06
M	8.5 - 10	6 - 8	FBS3015832B	AJ06
L	10.5 - 13	8.5 - 11	FBS3015833B	AJ06
XL	13.5+	11.5 - 14	FBS3015834B	AJ06

Note: Sizing by US shoe size.



Xceltrax Walker Boot - Standard

- Provides support and comfort to help promote healing for patients with mild to severe foot and ankle injuries
- Lightweight and ergonomic design
- Available in short or tall

Size	Shoe Size		Code		DVA:
	Female	Male	Short	Tall	
XS	3.5 - 5.5	2 - 4	FBS303437	FBS303432	AR19
S	6 - 8	4.5 - 7	FBS303438	FBS303433	AR19
M	8.5 - 11.5	7.5 - 10.5	FBS303439	FBS303434	AR19
L	11.5 - 13.5	10.5 - 12.5	FBS303440	FBS303435	AR19
XL	13.5+	12.5+	FBS303441	FBS303436	AR19

Note: Sizing by US shoe size.



Xceltrax Walker Boot - Air

- Provides support and comfort to help promote healing for patients with mild to severe foot and ankle injuries
- The soft liner has an integrated pump design making the pneumatic liner simple to inflate and adjust
- Available in short or tall

Size	Shoe Size		Code		DVA:
	Female	Male	Short	Tall	
XS	3.5 - 5.5	2 - 4	FPP303714	FPP303709	AR19
S	6 - 8	4.5 - 7	FPP303715	FPP303710	AR19
M	8.5 - 11.5	7.5 - 10.5	FPP303716	FPP303711	AR19
L	11.5 - 13.5	10.5 - 12.5	FPP303717	FPP303712	AR19
XL	13.5+	12.5+	FPP303718	FPP303713	AR19

Note: Sizing by US shoe size.



AirSelect Elite

- The Elite is the most advanced walking boot, engineered to provide the ultimate in protection, comfort and oedema control
- SoftStrike technology absorbs and dissipates shock, while the light weight and lab-tested rocker sole combine to encourage a natural gait and a continuation of everyday activities
- Suitable for use in: stable fractures of the foot/ankle/lower leg; severe ankle sprain; post-operative immobilisation

Size	Shoe Size		Code	DVA:
	Female	Male		
XS	Up To 5	Up To 4	FPP303703	AR01
S	5 - 8	4 - 7	FPP303701	AR01
M	8 - 11	7 - 10	FPP303700	AR01
L	11 - 15	10 - 13	FPP303699	AR01
XL	15+	13+	FPP303702	AR01

Note: Sizing by US shoe size.



AirSelect Short

- The Elite is the most advanced walking boot, engineered to provide the ultimate in protection, comfort and oedema control
- Suitable for use in: stable fractures of the foot/ankle/lower leg; severe ankle sprain; post-operative immobilisation
- SoftStrike technology absorbs and dissipates shock, while the light weight and lab-tested rocker sole combine to encourage a natural gait and a continuation of everyday activities

Size	Shoe Size		Code
	Female	Male	
XS	Up To 5	Up To 4	FPP303708
S	5 - 8	4 - 7	FPP303706
M	8 - 11	7 - 10	FPP303705
L	11 - 15	10 - 13	FPP303704
XL	15+	13+	FPP303707

Note: Sizing by US shoe size.



Drop Foot Ankle Foot Orthosis Open

- Limits plantar flexion, common peroneal nerve injuries
- Full length sole
- Can be adjusted for a custom fit by heating, cutting and grinding
- Available in left or right fitting

Size	Shoe Size	Code		DVA:
		Left	Right	
XS	4 - 5	FBS301588LXS	FBS301588	AR19
S	5 - 7	FBS301588LS	FBS301588RS	AR19
M	7 - 9	FBS301588LM	FBS301588RM	AR19
L	9 - 12	FBS301588LL	FBS301588RL	AR19

Note: Sizing by US shoe size.



Drop Foot Ankle Foot Orthosis

- Designed to hold the foot correctly
- Supports the foot to stop dragging
- Fits easily into shoe, heat mouldable, easily trimmed

Size	Shoe Size	Code		DVA:
		Left	Right	
S	Up to 38	FBS301589LS	FBS301589	AR19
M	39 - 43	-	FBS301589RM	AR19
L	44 - 46	FBS301589LL	FBS301589RL	AR19
XL	47+	FBS301589LXL	FBS301589RXL	AR19

Note: Sizing by European shoe size.



BodyAssist Heel Spur Heel Cups

- Silicone Medical Gel provide exceptional heel support, protection & comfort
- Ideal for sports use, gardening, factory staff, construction workers, office staff & for people who stand on hard surfaces
- Heel Spur Pads are made from a unique Medical Silicone Gel compound that assists to dissipate shock & shear, caused by impact to the heel & foot
- They emulate the fatty tissue, usually found on the underfoot

- Damaging stress forces normally absorbed by the body's musculoskeletal system, are now absorbed by the pad & away from the body

Size	Shoe Size	Code	DVA:
S	37- 40	FBS281210	AR04
M	40.5 - 42	FBS2812102	AR04
L	42.5+	FBS2812103	AR04

Note: Sizing by European shoe size.



Valgus 4050 Arch/ Metatarsal Support

- Anatomically designed moulded rubber with leather cover
- Metatarsal dome
- Available in narrow and wide fittings
- Supplied in pairs



Foot Up Drop Foot Orthosis

- A lightweight ankle-foot orthosis
- Provides visible improvements in gait by supporting the foot
- The cushioned ankle wrap is made from a breathable, 3 layer material
- The plastic inlay fits discreetly in the shoe and is practically invisible when worn
- Shoeless wrap accessory also available

Shoe Size	Code		
	Narrow	Wide	DVA:
33	FBS301590	-	AR04
34	FBS301590/100/34	-	AR04
35	FBS301590/100/35	-	AR04
36	FBS301590/100/36	-	AR04
37	FBS301590/100/37	-	AR04
38	FBS301590/100/38	FBS301590/200/38	AR04
39	FBS301590/100/39	FBS301590/200/39	AR04
40	FBS301590/100/40	FBS301590/200/40	AR04
41	FBS301590/100/41	FBS301590/200/41	AR04
42	FBS301590/100/42	FBS301590/200/42	AR04
43	-	FBS301590/200/43	AR04
44	-	FBS301590/200/44	AR04
45	-	FBS301590/200/45	AR04
46	-	FBS301590/200/46	AR04

Size	Ankle Circum (cm)	Foot Circum (cm)	Description	Code	DVA:
M	18 - 21	-	Foot Up Drop Foot Orthosis	FBS303504	AR19
L	22 - 26	-	Foot Up Drop Foot Orthosis	FBS303505	AR19
XL	27 - 33	-	Foot Up Drop Foot Orthosis	FBS303506	AR19
S	-	< 20	Wrap Accessory	FBS303507	AR19
M	-	20 - 24	Wrap Accessory	FBS303508	AR19
L	-	24 - 28	Wrap Accessory	FBS303509	AR19
XL	-	> 28	Wrap Accessory	FBS303510	AR19

Note: Sizing according to European (metric) shoe size.



BodyAssist Heel Spur Pads

- These Heel Spur Pads gently cushion the heel, while dissipating shock and shear generated by heel strike
- A removable 'tear drop' central plug allows the heel centre to be suspended around a supportive ring of silicone gel
- The pads fit snugly into the heel of a shoe

- Shock attenuation properties of this material have been proven effective in modifying heel strike trauma

Size	Shoe Size	Code	DVA:
S	Up to 38	FBS301598	AR04
M	38.5 - 42	FBS3015982	AR04
L	42.5+	FBS3015983	AR04

Note: Sizing by European shoe size.



BodyAssist Gel Shoe Insole

- Provides overall support, protection, comfort & help reduce painful shock waves to the body & lighten the load on joints
- Ideal for sports use, gardening, factory staff, construction workers, office staff & for people who stand on hard surfaces
- Made from a unique Medical Silicone Gel compound that assists to dissipate shock & shear, caused by impact to the heel & foot

- Suitable for managing minor cases of splayfoot
- Damaging stress forces normally absorbed by the body's musculoskeletal system, are now absorbed by the footbeds & away from the body

Size	Shoe Size	Code	DVA:
XS	37 - 38	FBS281200	AR04
S	39 - 41	FBS281201	AR04
M	41 - 43	FBS281202	AR04
L	44 - 46	FBS281203	AR04
XL	47+	FBS281204	AR04

Note: Sizing by European shoe size.



Trio Heel Elevators

- Trio elevators come in different heights and should be placed under the heel area of the orthotic
- Choose between sizes small, medium or large
- They are manufactured from extra firm density to prevent wearing out or crushing
- Available in 4, 6 and 8mm rises

Size	Shoe Size	Code			DVA:
		4mm	6mm	8mm	
S	Up to 37	FBS3016024S	FBS3016026S	FBS3016028S	AR04
M	38 - 43	FBS3016024M	FBS3016026M	FBS3016028M	AR04
L	44+	FBS3016024L	FBS3016026L	FBS3016028L	AR04

Note: Sizing by European shoe size.



Darco Peg-Assist

- Designed for use in most DARCO square-toe shoes, other than the DARCO Med-Surg Shoe (for which DA-PTQ Insoles should be used)
- Removable pegs allow for localised off-loading of wounds and ulcerations of the foot
- Pressure reduced by as much as 60% while allowing patient to remain ambulatory
- Poron® cover eliminates the incidence of ring oedema and edge abrasion

Size	Shoe Size		Code	DVA:
	Female	Male		
XS	Up to 5	-	FBS301699XS	AR04
S	5.5 - 8	-	FBS301699S	AR04
M	8.5 - 10	6 - 8	FBS301699M	AR04
L	10.5 - 13	8.5 - 11	FBS301699L	AR04
XL	13.5+	11.5 - 14	FBS301699XL	AR04

Note: Sizing by US shoe size.



Hallufix Bunion Aid

- A 3-in-1 aid that provides the benefits of a day and a night splint for: bunion pain relief, protection and correction
- Bunion Aid is a hinged bunion splint with adjustable straps and a comfortable gel cushion that protects irritated tissue, distributes pressure and provides corrective support
- Wearing this aid helps prevent bunion progression, overlapping toes, and hammer toes; limiting the need for surgery

Description	Size	Code	DVA:
Hallufix Bunion Aid	Universal	FBS303443	AR04
Hallufix Bunion Aid Replacement Kit	Universal	FBS303442	AR04



Footbionics Dual Density

- Provides good support for athletic activities or active patients
- A comfortable cushioning top layer of 2mm OrthoLite
- Can help with shin splints, Achilles tendon pain, stress fractures

Size	Shoe Size	Code	DVA:
XXS	32-34	FBS303503	AR04
XS	35-37	FBS303501	AR04
S	38-40	FBS303499	AR04
M	41-43	FBS303498	AR04
L	44-46	FBS303497	AR04
XL	46-48	FBS303500	AR04
XXL	49-51	FBS303502	AR04

Note: Sizing by US shoe size.



Achilles 1400 Heel Protection Sleeve

- The soft, comfortable gel pad protects the posterior of the heel from friction while gently compressing and reducing pressure
- Washable and reusable
- Slim profile design fits most footwear

Code	Description	Size	DVA:
FBS301606	Achilles 1400 Heel Protection Sleeve	One Size (Small/Medium)	AR01



Oppo 1408 Achilles Brace W/Silicone Pad

- The concentric ripples on the pad provide intermittent massage - like compression, to effectively reduce edema, relieve spasm and speed up the recovery process
- CoolPrene® material provides shock absorption, dynamic compression and better breathability on the injured ankle over prolonged use
- The butterfly shaped pad is anatomically designed to fit the Achilles tendon and provide support during motion

Size	Ankle (cm)	Code		DVA:
		Left	Right	
S	15.6 - 20.6	FBS301605LS	FBS301605	AR02
M	20.6 - 25.7	FBS301605LM	FBS301605RM	AR02
L	25.7 - 30.8	FBS301605LL	FBS301605RL	AR02
XL	30.8 - 35.9	FBS301605LXL	FBS301605RXL	AR02

Note: Measure circumference just above ankle bone.



Adjustable Night Splint

- Patented Formfit 3D molded technology allows for variable padding thickness in strategic locations for superior comfort and fit
- Soft lycra lined foam interior is plush and breathable, eliminating heat build-up
- Three lycra lined foam padded straps with user friendly centre release buckles ensure immobilisation while maintaining comfort across the flexor tendons
- Adjustable dorsi-flexion straps provide a gentle stretching of the plantar fascia and achilles tendon.

Size	Female Shoe Size	Male Shoe Size	Code	DVA:
S	Up to 7.5	Up to 7	FBS301607	AR04
M	8 - 10.5	7.5 - 10	FBS301607M	AR04
L	11+	10.5+	FBS301607L	AR04

Note: Sizing by US shoe size.



Hallux Valgus 510 Night Strap

- A light and flexible strap to unload the joints of the big toe
- Especially designed for nightly or post-operative usage
- Padded with skin-caring and breathable material
- Provides pain relief and guides the big toe
- Universal size - specify left or right foot

Code		Size	DVA:
Left	Right		
FBS301608L	FBS301608	Universal Size	AR04



Rehband Hallux Valgus 509 Day Strap

- A light and flexible strap to unload the joints of the big toe
- Padded for optimal comfort
- Can be used with or without shoes
- Gives pain relief and may be a condition for a normalised function of the muscles, guides the big toe
- Universal size - specify left or right foot

Code		Size	DVA:
Left	Right		
FBS301609L	FBS301609	Universal Size	AR04



Darco Toe Alignment Splint

- Low profile metatarsal band fits easily into most footwear and reduces slippage
- Soft toe straps will not irritate skin
- Special T-Strap gently guides toe(s) into proper alignment for post-procedural healing

- One sizes fits all
- Latex Free
- Colour - White

Code	Description	DVA:
FBS301610	Darco Toe Alignment Splint	AR02



Slip - Resistant Socks

- Provides warmth, comfort and peace of mind, with the tread on the bottom
- Available with single or double imprint
- Supplied in pairs

Size	Shoe Size	Colour	Code	DVA:
				Single Imprint
S	5 - 7	Teal	PTS642548	AR29
M	7.5 - 10	Tan	PTS642549	AR29
L	>10.5	Grey	PTS642550	AR29
XL	Bariatric	Royal Blue	PTS642551	AR29
Double Imprint				DVA:
S	5 - 7	Teal	PTS642544	AR29
M	7.5 - 10	Tan	PTS642545	AR29
L	>10.5	Grey	PTS642546	AR29
XL	Bariatric	Royal Blue	PTS642547	AR29

Note: Sizing by US shoe size.



Arthritis Sock

- Helps relieve aches, pains and stiffness associated with arthritis of the feet
- Provides mild compression for warmth and to help increase circulation
- Soft, breathable cotton material
- Supplied in pairs

Size	Female Shoe Size	Male Shoe Size	Code
S	6 - 8	5 - 7	PTS642552
M	8.5 - 11.5	7.5 - 9.5	PTS642553
L	12 - 15.5	10 - 14	PTS642554

Note: Sizing by US shoe size.



For Mobility
Walking Aids
Refer to Page 255

Therapy & Exercise

DVA Contact Details

P 1300 888 052 | F 1300 787 052 | E dva@aidacare.com.au



For NDIS Item Code,
refer to page: 500