

02/12/21

WHILL MODEL C2 - AVAILABLE SOON!

The next generation of Whill Wheelchairs is now here, introducing the Whill C2 which arrives soon! Following on from the success of the Whill C1 Wheelchair, the C2 has undergone significant design improvements to further improve its position as the coolest power wheelchair in the market. Unlike the C1, the C2 is exclusive to Aidacare, which improves our competitive ability.

IMPROVED FEATURES FROM C1 TO C2

- Size options: now available in 450mm & 500mm seat widths, suiting a larger range of users.
- Higher SWL: increased from 114kg (C1) to 136kg (C2).
- Rear wheel suspension: Absorbs the impact when driving over bump roads or obstacles.
- **Mobile phone app compatibility:** access the Whill phone app and operate your chair from anywhere, on all Apple & Android devices (previously only Apple).
- **Improved user interface:** ensures all functional controls are in one place, allowing for simple chair control by the user.
- Rear tail lights: improves safety for users in the chair when crossing roads or in high density areas.
- **Standard inclusions:** lap belt, smart key, black side covers, battery and battery charger.



New Product Launch



PRODUCT CODES

Product Code	Description
Wheelchairs	
PWS655400	Power Wheelchair - Whill - Model C2 - (Black Armrests) - 450mm - 136kg
PWS655410	Power Wheelchair - Whill - Model C2 - (Black Armrests) - 500mm - 136kg
Accessories	
WCA870400	Power Wheelchair Accessory - Whill Model C2 - Lithium-Ion Battery
WCA870410	Power Wheelchair Accessory - Whill Model C2 - Charger
WCA870420	Power Wheelchair Accessory - Whill Model C2 - Lap Belt
WCA870430	Power Wheelchair Accessory - Whill Model C2 - Side Storage Bag
WCA870440	Power Wheelchair Accessory - Whill Model C2 - Storage Basket
WCA870450	Power Wheelchair Accessory - Whill Model C2 - Walking Stick Holder
WCA870260	Power Wheelchair Accessory - Whill Model C1 & C2 - Dust Cover
WCA870460	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Black (Left)
WCA870470	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Black (Right)
WCA870480	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Gold (Left)
WCA870490	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Gold (Right)
WCA870500	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Grey (Left)
WCA870510	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Grey (Right)
WCA870520	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Light Blue (Left)
WCA870530	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Light Blue (Right)
WCA870540	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Navy (Left)
WCA870550	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Navy (Right)
WCA870560	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Pink (Left)
WCA870570	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Pink (Right)
WCA870580	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - White (Left)
WCA870590	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - White (Right)
WCA870600	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Red - Left
WCA870610	Power Wheelchair Accessory - Whill Model C2 - Arm Covers - Red - Right

