

01/12/21

ASPIRE RAPHAEL QUATTRO - NOW WITH EXTENDED LEGREST!

After extensive feedback from the market, we have increased the legrest length on our Aspire Raphael Quattro Lift Chairs. The legrest now extends further out, providing a longer support under the lower legs, which is of particular benefit to taller users.

This change will affect all Raphael Quattro lift chairs stocked colours Bronze, Nordic Stone & Marina Grey. Custom chairs (upholstered locally) will not include this change for a few months.

KEY DIFFERENCES, IMPROVEMENTS & STOCK MANAGEMENT

- o **Codes:** remain the same for the S/M/L chairs however the MAXI chairs codes will change due to a reduction in the SWL for this particular size *(previously 220kg, now 180kg)*.
 - CHP227320 is superseded by CHP227319
 - CHP227321 is superseded by CHP227322





