



# Falls Prevention Solutions

Wired & Wireless Range

Exclusively distributed by



# DELUXE FLOOR ALERTAMAT



Alarm Monitor - Wired



Wireless Option Available

Plush carpet floor mat that provides a safe falls prevention solution that suits modern care environments

Distinguished look for the care environment and is an ideal solution for situations where a patient with dementia would otherwise avoid the mat

Simple plug and play ensures carer can quickly set up and operate with ease

A robust, long life floor alarm

Long lasting welded seams for durability

Large size to prevent overstepping by the user

Available in both wired and wireless configurations to suit the necessary needs of Home, Aged and Health care environments

Code	Description	Length	Width	Warranty
BEA007615	Deluxe Floor Alertamat System (includes mat, alarm monitor and cables)	900 mm	600 mm	1 Year
BEA007570	Wireless Deluxe Floor Alertamat System (includes mat and alarm monitor)	900 mm	600 mm	1 Year



# FLOOR PLUS ALERTAMAT



Alarm Monitor - Wired



Wireless Option Available

Heavy duty floor mat provides a safe falls prevention solution with maximum longevity

Simple plug and play ensures carer can quickly set up and operate with ease

A robust, long life floor alarm

Long lasting welded seams for durability

Large size to prevent overstepping by the user

Hygienic, anti-bacterial easy-clean material

Available in both wired and wireless configurations to suit the necessary needs of Home, Aged and Health care environments

Code	Description	Length	Width	Warranty
BEA007620	Floor Plus Alertamat System (includes mat, alarm monitor and cables)	800 mm	550 mm	1 Year
BEA007575	Wireless Floor Plus Alertamat System (includes mat and alarm monitor)	800 mm	550 mm	1 Year





# BED ALERTAMAT



Alarm Monitor - Wired



Wireless option available

Provides advanced warning of falls risk if patient leaves or attempts to leave bed

Simple plug and play ensures carer can quickly set up and operate with ease

High and low volume alert

Adjustable delay to reduce false alarms

Alarm is reset automatically when patient returns to the bed

Hygienic, antibacterial easy-clean material

Available in both wired and wireless configurations to suit the necessary needs of Home, Aged and Health care environments

Code	Description	Length	Width	Warranty
BEA007605	Bed Alertamat System (includes mat, alarm monitor and cables)	760 mm	250 mm	1 Year
BEA007565	Wireless Bed Alertamat System (includes mat and alarm monitor)	760 mm	250 mm	1 Year



# CHAIR ALERTAMAT



Alarm Monitor - Wired



Wireless option available

Alerts carers if patient leaves or attempts to leave chair, helping to manage falls risk

Simple plug and play ensures carer can quickly set up and operate with ease

High and low volume alert


Adjustable delay to reduce false alarms

Alarm is reset automatically when patient returns to the chair

Hygienic, antibacterial easy-clean material

Available in both wired and wireless configurations to suit the necessary needs of Home, Aged and Health care environments

Code	Description	Length	Width	Warranty
BEA007600	Chair Alertamat System (includes mat, alarm monitor and cables)	380 mm	250 mm	1 Year
BEA007657	Wireless Chair Alertamat System (includes mat and alarm monitor)	250 mm	250 mm	1 Year





# ALERTA DETECT INFRARED SYSTEM



Alarm Monitor - Wireless



Infrared Sensor

Infrared beam senses when a patient leaves their bed, alerting carers to minimise potential falls risk

Can also be used to allow the patient to move around their room freely, only alerting carers if they leave their room

Can be used as a standalone system (sensor & alarm monitor via infrared "IR") to create an alarm within the room or directly placed into nurse call system (sensor only)

Wireless plug and play ensures carer can quickly set up and operate with ease

System includes: Sensor & Alarm Monitor, Stereo & Mono Nurse Call Cables

Wireless solution that removes excess cables and provides a streamline falls prevention solution

Large sensor range of up to 95m ensures all types of room configurations are covered

Can be battery operated (sold separately) or via a mains power cable (optional accessory)

Alarm Monitor can be carried around in the homecare space by a carer, fitting neatly into the pocket

Code	Description	Length	Width	Warranty
BEA007610	Alerta Detect Infrared System (includes sensor, alarm monitor & cables)	760 mm	250 mm	1 Year



# ALERTAMAT CRASH MAT WITH WIRELESS NURSE CALL ALERT



Alarm Monitor - Wireless



Alerta Wireless Transmitter

Soft padded mat helps to minimize serious injury in the event of a fall from bed

Alerts carers if patient impacts with mat, whether due to an accidental fall or controlled exit of the bed

Wireless connect-and-play system

No cables, extra safe and versatile

95m wireless range from transmitter to receiver

Genuine non-slip base

System includes mat, transmitter, alarm monitor & cables

Code	Description	Length	Width	Warranty
BEA007580	Alerta Crash Mat (includes mat, transmitter, alarm monitor & cables)	1800 mm	600 mm	1 Year



# WIRELESS NURSE CALL BUTTON



Provides the user with reassurance that help is at hand. With a large easy to press button it activates the alarm when a user presses the button. This alerts the caregiver to ascertain the user's need.

Connects to a nurse call system through the Alerta Alarm Monitor for centralised activation and comes with a lanyard, and a holder that can be wall mounted, so it can always be used remotely within the homecare environment.

Simple plug and play ensure carer can quickly set up and operate with ease.

High and low volume alert.



Code	Description	Warranty
BEA007631	Alerta Wireless Nurse Call Button (includes button and alarm monitor)	1 Year

## SPARES, ACCESSORIES AND PLUG IDENTIFICATION



### Wired Range

- BEA007655 Chair Alertamat (mat only)
- BEA007660 Bed Alertamat (mat only)
- BEA007616 Deluxe Alertamat (mat only)
- BEA007621 Floor + Alertamat (mat only)
- BEA007625 Nurse Call Button - Mono
- BEA007630 Nurse Call Button - Stereo
- BEA007665 Alarm Monitor - Wired

### Wireless Range

- BEA007657 Chair Alertamat (mat only)
- BEA007662 Bed Alertamat (mat only)
- BEA007618 Deluxe Floor Alertamat (mat only)
- BEA007623 Floor + Alertamat (mat only)
- BEA007632 Nurse Call Button (button only)
- BEA007670 Alarm Monitor - Wireless / AlertaDetect
- BEA007675 AlertaDetect Sensor

### Cables / Accessories

- BEA007635 Double Adaptor - Mono
- BEA007640 Double Adaptor - Stereo
- BEA007645 Extension Cable - Mono
- BEA007650 Extension Cable - Stereo
- BEA007685 Standard Cable - Stereo
- BEA007690 Standard Cable - Mono
- BEA007691 Standard Cable - Ring/Sleeve
- BEA007692 Standard Cable - Tip/Ring
- BEA007695 Dual Cable - Stereo
- BEA007700 Dual Cable - Mono
- BEA007705 Floor Mat to Alarm Cable

### Transmitters

- BEA007672 Transmitter for Wireless Systems



**MONO PLUG**



**STEREO PLUG**

Mono and Stereo cables are the most common connection types in nurse call systems in Australia. Alerta alarm systems are supplied with both mono and stereo cables when purchased as a system, to ensure there is no delay upon installation.

If you require further assistance please contact: [sales@aidacare.com.au](mailto:sales@aidacare.com.au)