

IDEA CLINIC TREATMENT CHAIR




COLOUR OPTIONS




Blue Vinyl

OPTIONAL ACCESSORIES



IV Pole CHT100720



Padded Armrest (Pair) CHT100725



BIOCOTE
Antibacterial frame coating for improved infection control



FEATURES AND BENEFITS

- Compact therapy chair that effortlessly meets the needs of oncology departments and the demands of infusion and transfusion medicine
- Gas-assisted, height and laterally adjustable swing away armrests to suit individual user needs. Armrest can be removed to facilitate side transfer
- Ideal for procedures such as chemotherapy, dialysis, blood donation, and blood transfusion
- Flip down footrest for use during portering or patient movement
- Electric adjustment allows for infinite positioning of the backrest, legrest, seat height, and Trendelenburg (CPR shock position)
- Seamless patient surface for maximised infection control compliance
- Genuine chair position with leg section at 90° angle
- Intuitive handset allows for quick and easy chair repositioning
- Ergonomic push handle with inbuilt paper roll holder
- 100mm castors with directional tracking ability
- Inbuilt IV pole receiver, additional pole available as an accessory



Code	CHT100700
Overall Length	910mm
Overall Width	730mm
Seat Height	480 - 780mm
Patient Surface Length	1720mm
Patient Surface Width	560mm
Backrest Angle	80°
Trendelenburg	6°
Castor Size	100mm