



Product Code: LSS403830 (Includes Swift & Rhino)

MANGAR Swift Transfer Slide

Transfer and reposition patients in the supine position with minimal moving and handling force needed, while providing a stable, safe experience for the patient.

Maximum User Weight: 500 kg (1102 lbs/79 st)

Inflated Height (approx): 14 cm (5.5")

Overall Length (approx): 204 cm (80")

Overall Width (approx): 100 cm (39")

Weight: 1.5 kg (3.3 lbs)

MANGAR Rhino Lifting Cushion

Lifts a person in the supine position from the floor to the height of a bed or stretcher.

Maximum User Weight: 500 kg (1102 lbs/79 st)

Inflated Height (approx): 65 cm (26")

Overall Length (approx): 204 cm (80")

Overall Width (approx): 112 cm (44")

Weight: 7.5 kg (16.5 lbs)

MANGAR SUPINE TRANSFER SYSTEM

The Mangar Patient Transfer System enables healthcare workers to lift and transfer patients with minimal moving and handling, while providing a stable, safe experience for the patient.

The system consists of the Mangar Rhino Lifting Cushion, the Mangar Swift Transfer Slide, and the Mangar Airflo Duo Compressor.

- Mains and battery operated
- Lightweight anti-rip material
- Complete Supine Transfer System weighs under 20 kg



Features

- 1. Roller Bag** - Fits the Rhino, Swift, Duo Compressor and Mains Power Supply with a total weight under 20 kg.
- 2. Lightweight, durable anti-rip material** - Easy to fold and return to bag after use.
- 3. Mains and battery powered** - Airflo Duo Compressor can be transported easily and used anywhere.
- 4. Patient Markings** - To show optimum position of patient on Rhino and Swift.




For more information visit



1300 133 120 | aidacare.com.au
 sales@aidacare.com.au



RHINO



SWIFT



AIRFLO DUO