



Comfort and performance in a proven design



Live Without Limits

The Quickie Pulse 5 is a compact, durable power wheelchair that offers a wide range of effective seating and electronics options.

QUICKIE[®] Pulse 5

Available in 2 Models:

Quickie Pulse 5 (CC): 8 Kph, 2-Pole motor, No Power Option base with Captain Seat Quickie Pulse 5 (BC): 8 Kph, 2-Pole motor, Tilt and ASAP I, Adjustable Seat Frame



POWER BASE DESIGNED FOR PERFORMANCE

- Mid wheel drive
- Reliable and compact for indoor accessibility
- Flexible suspension system for optimal outdoor performance
- Sleek base style available with shroud colour options: Candy Apple Red or Candy Blue



ELECTRONICS OPTIONS

VR-2[™] Non-Expandable Controls

- Industry standard non-expandable electronics
- Dependable, highly flexible, and simple to use
- 2 button Thru Drive for Tilt models



MULTIPLE SEAT FRAME OPTIONS

Captain Seat (CC)

- Fully adjustable seat size, legrest placement and armrest width
- Convertible seat pan for flexibility in positioning solutions
- Sleek, comfortable style
- Semi-recline and fold-down backrest
- Mounting accessories
 - O₂ Holder
 - Crutch Holder
 - Walker Holder
- Available in 16", 18" and 20" width and 15"- 22" depth
- Seat to floor height adjustment from 16.25" 18.25"









Additional Backrest Options For Asap I Seat Frames (BC)





The JAY J3[™] back's stylish design, quick-release hardware, and ability to meet a multitude of clinical needs make it an excellent complement to all Pulse models.



The SPOT Curved back's solid support system and option of Whitmyer headrests and positioning components make it a fantastic option on the Pulse.

ASAP SEAT (BC)

- ASAP seat frame with static tilt adjustment of 0° and 4°
- SPOT Tilt on ASAP seat frame provides up to 50° of tilt
- Backrests: Jay Curved back or Jay 3 back
- In-field width, depth and back height adjustability 16 -20" Wide, 16- 20" Deep



COMPLEMENT YOUR QUICKIE PULSE WITH A JAY® CUSHION

Order your JAY[®] Cushion shipped complete on-chair with any Quickie Pulse chair.

JAY Union™

JAY[®] Cushions

Comfort Cushions





JAY Soft Combi[™] P

High Needs Cushions





JAY Easy™

100375



JAY Fusion™

Quickie Pulse 5 Specifications



TGA #



100375