

Sharky Bathroom Chair

Sharky has been designed for multifunctional applications and ease of handling.

Multifunctional commode for bath and shower use



Removable seat with a range of adapters offers versatile support in different areas of the bathroom

Removable rear seat section to assist with hygiene care

Easy to clean all surfaces

Extremely slim, compact foldable frame

Detachable footrest system

BACK SUPPORT RECLINE

Simple to use system for easy and secure handling.



45° TILT

The special construction of the tilt mechanism allows to 45° tilt of seat to ensure comfort and high safety especially when showering.



HEIGHT ADJUSTMENT

Smooth movement with comfortable handles make adjustments easy and controlled.



COMPACT FOLDING

Fold down to compact height of just 28cm, even with seat attached, makes transport and storage easy.



REMOVABLE PADDING

Easy care removable padding for cleaning the Sharky.



Size	Base Code	Seat Width	Seat Depth	Back Support Height	Toilet Clearance	Seat Angle	Back Support Angle	Knee Angle	Foot Support Angle	Armrest Height	Wheel Diameter	
1	6081-0015-000	24, 27, 30cm	27-36cm	42cm	52cm	-10° to 45°	72-120°	78-138°	80-100°	14,5 - 20,5cm	75, 100cm	60kg
2	Launching Soon	29, 32, 35cm	27-40cm	42cm	52cm	-10° to 45°	72-120°	78-138°	80-100°	14,5 - 20,5cm		

Accessories

Removable Footrest System
6085-1810-000



Push Handles
6085-8600-000



Lateral Trunk Supports
1: 6085-3310-000
2: 6085-3320-000



Headrest
6085-4750-000



**Neoprene Shoulder
Anterior Chest Support**
1: 6085-6610-000
2: 6085-6620-000



**Hip Belt with Neoprene
Pad Sizes**
1: 6085-1910-000
2: 2880-1920-000



**Removable & Height
Adjustable armrest**
6085-5500-000



Abduction Block
1: 6085-1500-00



**Depth Adjustable,
removable tray table**
6085-7100-000



Neoprene Calf Strap
6085-6160-000



Neoprene Ankle Huggers
1: 6085-6101-000
2: 6085-6102-000



**Removable Commode
Pan with retainer**
6085-8210-000



Removable Splash Guard
6085-7675-000



Seat Infill pad
6085-0401-000



**Bath subframe with seat
angle adjustment**
6085-0500-000



Toilet Adapter
6085-0600-000



Seat width reducing hip supports
1: 6085-1110-000
2: 6085-1120-000
3: 6085-1130-000

