


Otto Bock®

QUALITY FOR LIFE



Sleepform

Right from the start, children want to discover the world and expand their horizons. The jillions of new impressions they receive every day are processed in dreams at night. These sleeping phases alternate with the phases of physical recovery. Children between the ages of five and twelve sleep ten to eleven hours on average, while younger children may need up to 14 hours of sleep. That is why it is very important for them to have ideal support during this time, too. The innovative positioning system called "Sleepform" will let your child sleep not only comfortably at night, but also in the therapeutically correct position.



Otto Bock Kids

LECKEY®

Designed by Leckey®

Otto Bock HealthCare GmbH

Max-Näder-Straße 15 · 37115 Duderstadt/Germany · Phone: +49 5527 848-1590 · Fax +49 5527 848-1676 · healthcare@ottobock.de · www.ottobock.com



Sleepform



- It takes only a few steps to arrange the vacuum mattress in the form desired for your child.
- It does not matter whether your child sleeps on his or her back, stomach or side, you can adapt the "Sleepform" individually to the needs of your little one.



- The top Airflow mattress provides air circulation and keeps the temperature of the child the same.
- All the components together ensure that the child will readily accept it and sleep well at night.



- The Sleepform also has hip-/chest belts, leg belts and knee positioning rollers for positioning the knees to provide the child with additional support.
- The elastic belts ensure that the child has enough freedom to move, while constantly bringing him or her back to the ideal position.



Technical Data for Sleepform

| | |
|---|--|
| Airflow mattress, size 1 & 2 | 56 × 118 cm |
| Airflow mattress, size 3 & 4 | 100 × 200 cm |
| Sleepform vacuum mattress, size 1 – 4 | (includes mattress, pump and bag) 95 × 95 cm |
| Padded chest and hip belt, size 1 | 0 – 1 years old for crib |
| Padded chest and hip belt, size 2 | 1 – 5 years old for crib |
| Padded chest and hip belt, size 3 | 1 – 5 years old for single bed (100 × 200 cm) |
| Padded chest and hip belt, size 4 | 5 – 18 years old for single bed (100 × 200 cm) |
| Padded leg and knee belt, size 1 | 0 – 1 years old for crib |
| Padded leg and knee belt, size 2 | 1 – 5 years old for crib |
| Padded leg and knee belt, size 3 | 1 – 5 years old for single bed (100 × 200 cm) |
| Padded leg and knee belt, size 4 | 5 – 18 years old for single bed (100 × 200 cm) |
| Airflow knee positioning pillow, size 1 & 2 | 23 × 33 cm |
| Airflow knee positioning pillow, size 3 & 4 | 45 × 65 cm |
| Knee roller, size 1 & 2 | 14 × 50 cm |
| Knee roller, size 3 & 4 | 19,5 × 90 cm |
| Attachment belt, size 1 & 2 | Length: 63 cm |
| Attachment belt, size 3 & 4 | Length: 102 cm |

Otto Bock®

QUALITY FOR LIFE