



82 Albion Street
P.O Box 2085, Lygon St Nth
Brunswick East, VIC 3057
Tel: (03) 9384 1846
Fax: (03) 9384 2088

ABN 30 117 449 911

'Leading national provider of healthcare equipment & services'

4 Noonan Place
Ingleburn 2565
Tel: (02) 9618 5000
Fax: (02) 9618 5111

313 Warrigal Rd
Cheltenham VIC 3192
Tel: (03) 9584 1644
Fax: (03) 9584 0562

2/603-611 Whitehorse Rd
Mitcham VIC 3132
Tel: (03) 9872 4044
Fax: (03) 9872 6433

74-76 Fyans St
Geelong VIC 3220
Tel: (03) 5223 2263
Fax: (03) 5223 2269

Shop 2 826 David St
Albury NSW 2145
Tel: (02) 6041 1037
Fax: (02) 6023 2838

Unit 3 47 Learoyd Road
Alger QLD 4115
(07) 3272 6900
(07) 3272 4930

162 Great Western Highway
Westmead NSW 2145
Tel (02) 9687 1699
Fax (02) 9687 2477


OSD WHEELCHAIR MASS ORDER FORM (JULY 2009)

Date:	
Customer:	Phone:
Therapist:	Phone:
E-mail:	Fax:

OSD Transit \$350

- 40cm (16")
 45cm (18")
 50cm (20")
- 

OSD Standard \$350

- 35cm (14")
 40cm (16")
 45cm (18")
 50cm (20")
- 

OSD Lightweight \$625

- 35cm (14")
 40cm (16")
 45cm (18")
 50cm (20")
- 

FREIGHT	<input type="checkbox"/> Q1 SEQ \$30 <input type="checkbox"/> Q2 SOUTHERN \$50 <input type="checkbox"/> Q3 CENTRAL \$60 <input type="checkbox"/> Q4 NORTHERN \$80
<input type="checkbox"/>	Elevating Legrest LEFT \$90
<input type="checkbox"/>	Elevating Legrest RIGHT \$90
<input type="checkbox"/>	Stump Support LEFT (OSD Lightweight Only) – chair set up at amputee setting \$135
<input type="checkbox"/>	Stump Support RIGHT (OSD Lightweight Only) – chair set up at amputee setting \$135
<input type="checkbox"/>	Anti-Tip Bars (OSD Lightweight Only) \$50
<input type="checkbox"/>	Set-Up at Amputee Setting (when not ordering stump supports) \$80
<input type="checkbox"/>	Solid Tyres (OSD Standard & Lightweight Only) \$110
<input type="checkbox"/>	Oxygen Bottle Carrier (C-Size) \$95
<input type="checkbox"/>	Pelvic Strap (Bodypoint Non-Padded) <input type="checkbox"/> Medium HB003A1 <input type="checkbox"/> Large HB008A1 \$80
<input type="checkbox"/>	Calf Strap (Bodypoint Padded) <input type="checkbox"/> Medium SP102M <input type="checkbox"/> Large SP102L \$95
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	

TECHNICAL SPECIFICATIONS

	OSD Transit	OSD Standard	OSD Lightweight
Seat Depth	400mm	400mm	435mm
Back Height	400mm	400mm	450mm
Back Upholstery	Standard Sling	Standard Sling	Tension Adjustable
Castors	8"	8"	6"
Wheels	16" Pneumatic Mag	24" Pneumatic Spoke	24" Pneumatic Spoke
Axles	Quick Release	Quick Release	Quick Release
Armrests	Flip Back Removable Full Length	Flip Back Removable Full Length	"T" Height Adjustable Desk Length
Weight With/ Without Wheels	15.5/13.5kg	17.5/13.5kg	15.5/11.5kg
Load Capacity	100kg	100kg	100kg

Please return to Lifehealthcare staff or fax to (07) 3272 4930