

Start^{M4} XXL
 Lightweight Wheelchair

Article number: 480F53=40000



Standard Equipment

- Aluminium tube frame with double crossbrace
- Tip-assist integrated into the frame
- Individual footrests, angle-adjustable, plastic footplate
- Swing-away footrest holders
- Desk side panel, different height settings possible; with armrest
- Seat and back upholstery of black nylon
- Modular design
- Rear wheels with reinforced spokes
- Wheel lock

Folding; max. user weight: 160 kg

| | |
|---|--|
|  <p>Lifehealthcare Hospital & Rehab</p> | Order no.: <hr style="width: 100%;"/> |
|---|--|

| | | | |
|-------------------------|------|---------------|------|
| FREIGHT CHARGES: | | | |
| Q1 (SEQ) | \$60 | Q3 (CENTRAL) | \$90 |
| Q2 (SOUTHERN) | \$80 | Q4 (NORTHERN) | \$95 |

| | |
|-----------------------------------|---|
| Invoice address / Customer number | |
| Delivery address | |
| Sales person | Name <hr style="width: 80%; margin-left: 0;"/> Date Signature |

| Variant | Price |
|-----------------|--|
| 480F53=40000*** | <input type="checkbox"/> Start ^{M4} XXL – Lightweight Wheelchair |
| | Standard Seat Width (SW) in cm |
| | <input type="checkbox"/> 50.5 <input type="checkbox"/> 55.5 <input type="checkbox"/> 58 |
| | Custom Seat Width in cm |
| 480F53=40000_K | <input type="checkbox"/> 45.5 <input type="checkbox"/> 48 <input type="checkbox"/> 53 |
| | Seat Depth incm |
| | <input type="checkbox"/> Standard 43cm |
| | <input type="checkbox"/> Custom - Please Select Depth <input type="checkbox"/> 46cm <input type="checkbox"/> 49cm |
| | Back Height |
| | <input type="checkbox"/> Standard 45cm |
| | <input type="checkbox"/> Custom - Please Select Height <input type="checkbox"/> 40cm <input type="checkbox"/> 42.5cm <input type="checkbox"/> 47cm |
| | Frame Accessories |
| MA62 | <input type="checkbox"/> Anti-Tipper Right - Pendulum Style (Weight Limit 160 kg) |
| MA63 | <input type="checkbox"/> Anti-Tipper Left - Pendulum Style (Weight Limit 160 kg) |
| MA55 | <input type="checkbox"/> Plug-on Anti-tipper - Right (Maximum Weight Limit 125 kg) |
| MA56 | <input type="checkbox"/> Plug-on Anti-tipper - Left (Maximum Weight Limit 125 kg) |
| MA57 | <input type="checkbox"/> Crutch Holder - Right |
| MA58 | <input type="checkbox"/> Crutch Holder - Left |
| | Casters |
| MF10 | <input type="checkbox"/> 8" Solid Rubber (FSH settings 47.5/49/50/51.5cm) |
| MF06 | <input type="checkbox"/> 6" Solid Rubber (FSH settings 45/46.5/47.5/49cm) |
| | Brakes |
| MH01 | <input type="checkbox"/> Wheel Lock - Standard |
| MH06 | <input type="checkbox"/> Wheel Lock for User and Attendant |

