

Lifehealthcare Queensland  
 3/47 Learoyd Road, Algester QLD 4115  
 ph: 07 3272 6900 fax: 07 3272 4930  
 www.lifehealthcare.com.au

**Otto Bock®**

QUALITY FOR LIFE

**Start<sup>M5</sup> Comfort**  
 Adaptive Wheelchair

Article number: 480F53=50000

	Order no.:
	_____

FREIGHT CHARGES:			
Q1 (SEQ)	\$60	Q3 (CENTRAL)	\$90
Q2 (SOUTHERN)	\$80	Q4 (NORTHERN)	\$95

Invoice address / Customer number
_____

Delivery address
_____

Sales person	Name
_____	_____
	.....
Date	Signature
_____	_____



**Standard Equipment**

- \* Aluminium Frame
- \* Integrated Tip assist
- \* Tension Adjust Back Upholstery
- \* Rear Wheels with Quick Release Axles
- \* Adjustable Rear Seat to Floor Height
- \* Height Adjustable Back Rest 40 - 47.5cm
- \* Adjustable Seat Depth 40 - 46cm
- \* Flip Up Armrests with 1 - 2cm Height Adjustment
- \* Swingaway Removable Individual Footrests

**Maximum User weight 125Kg**

Variant		Price
480F53=50000***	Start <sup>M5</sup> Comfort – Lightweight Wheelchair	\$1,200
	<b>Seat Width (SW) in cm</b>	
	<input type="checkbox"/> 38 <input type="checkbox"/> 40.5 <input type="checkbox"/> 43 <input type="checkbox"/> 45.5 <input type="checkbox"/> 48	
	<b>Seat Depth in cm</b>	
	<input type="checkbox"/> Standard 43cm	Incl.
	<input type="checkbox"/> Custom - Please Select Depth <input type="checkbox"/> 40cm <input type="checkbox"/> 46cm	\$120
	<b>Back Height</b>	
	<input type="checkbox"/> Standard 42.5cm	Incl.
	<input type="checkbox"/> Custom - Please Select Height <input type="checkbox"/> 40cm <input type="checkbox"/> 45cm <input type="checkbox"/> 47cm	\$120
	<b>Frame Accessories</b>	
MA55	<input type="checkbox"/> Plug-on Anti-Tipper - Right	\$29
MA56	<input type="checkbox"/> Plug-on Anti-Tipper - Left	\$29
MA57	<input type="checkbox"/> Crutch Holder - Right	\$36
MA58	<input type="checkbox"/> Crutch Holder - Left	\$36
	<b>Casters</b>	
MF06	<input type="checkbox"/> 6" Solid Rubber (FSH settings 45/46.5/47.5/49cm)	Incl.
MF04	<input type="checkbox"/> 5.5" Soft Rubber Casters (reduces FSH by 0.5cm)	\$151
MF10	<input type="checkbox"/> 8" Solid Rubber (increases FSH by 2.5cm)	\$211
	<b>Caster Adapters (Hemi) (FSH with 6" castor 39.5/40.5/42cm)</b>	
MF81	<input type="checkbox"/> Hemi Adaptor + Short Fork - Right	\$71
MF81	<input type="checkbox"/> Hemi Adaptor with Short Fork - Left	\$71
	<b>Brakes</b>	
MH01	<input type="checkbox"/> Wheel Lock - Standard	Incl.
MH06	<input type="checkbox"/> Wheel Lock for User and Attendant	\$307

