




## New Enhanced Features

- 40 degree tilt range option
- Advanced stroller handle
- Fun spoke guards
- Swing in / out hangers




## The New Zippie Iris Product Improvements

- Easier to reach tie down locations
- Ridge in rocker prevents lateral play
- Improved hub-lock design



**Please note: The new version of the Zippie Iris is available from 16th November and the old version will be available for order until 30th January 2010.**

### True Fit Growth



The next larger width strut tube kit is available at time of order at no charge.

### Zippie Iris True Fit Growth

However, when your client requires a significant change in overall size, the following items are included in the True Fit program.

The **NEW Zippie Iris True Fit program is a one-time frame size change, and includes all the parts necessary to grow or shrink the wheelchair in the first 5 years from date of purchase of the original owner.**

The program includes strut tubes, seat rails, back posts, seat pan, hangers, extension tubes, footplates and tilt cables, if required.

**Note: Items are only exchanged/grown if ordered on original chair.**

### True Fit Growth Items

- Strut tubes (5 strut tubes included)
  - Seat rails
 (if outside current seat depth growth range)
  - Back posts/Stroller handles
 (when necessary to increase back height beyond current available range)
  - Push handle width insert
  - Seat pan
 (if outside current available range)
  - One arm drive axles
 (when already supplied)
  - Hangers
  - Extension tubes
  - Leg strap
  - Footplates
  - Tilt cables