

## Lifehealthcare MASSIVE Paediatric Ex-Demo Sale

### Kimba Spring Stroller Size 1



Like New



~~New Price \$6500.00~~ Now \$4000.00

The Kimba Spring is an indoor / outdoor stroller for children requiring postural support. Built in suspension ensures a smooth ride and folds easily for transporting. With height adjustable push handle, angle adj. backrest, tilt in space, as well as a number of seating support accessories.

### Convaid Profiler Transit Tilt 14"



~~Was \$4934.00~~

Now \$1500.00

Convaid's tilt Profiler adjustable contour with patented "cat's cradle" concept is a dynamic contour that allows some yield as occupant's position is changed slightly. The contour can be adjusted for several years to accommodate growth or changing posture. No more discomfort and pressure points from continuous growth with easy access to the occupant's back to feel for perfect fit.

### Easy Seat with Castor Legs Size 1 & 2

The Easy Seat is a sturdy height adjustable chair for the school or home environment. It is specifically designed to provide postural support for children with positioning needs. Seat width and depth adjustment allows for growth in the chair and there are a large range of accessories.



~~Size 1 Was \$1495~~

~~Size 2 Was \$1578~~

~~Size 3 Was \$1780~~



Now \$900

Now \$1000.00

Now \$1100.00

### Leckey Advance Bath Chair Size 2

~~Was \$1695.00~~ Now \$999.00



The Advance Bath Chair supports and holds the child in the bath. The easy to use adjustable legs can be setup prior to use to a number of different seat heights and angles. When in the bath, a simple lever can be used to adjust the backrest recline position. The mesh fabric is breathable, soft to touch and can be easily washed. Folds flat for easy storage and transportation.

## ***Leckey Advance Seat / Therapy Chair with Hi Low Base Size 1***



**Was \$5941.00**      **Now \$1800.00**



The Leckey therapy chair is height adjustable, has tilt & recline, Adjustable seating and has many postural supports supplied. (Colour of Ex-Demo is Yellow)

## ***Leckey Totstander***



**Was \$1880**      **Now \$1100.00**

Height and width adjustable flexible lateral supports. Padded Chest plate adjusts in height and angle for multi positioning. Multi adjustable pelvic band. Height and angle adjustable knee plate with laterally adjustable padded knee cups.

## ***Powerchairs***

### ***Little Ripper Powerchair***



**Was \$4500.00**  
**Now \$1900.00**

### ***Quantum Dynamo***



**Was \$6500.00**  
**Now \$2400.00**