

## GelBodies Heel & Ankle Protector



- Combined with pressure relieving gel, which is found in high-end wheelchair cushions, to offer additional padding to points on the body where pressure sores tend to form
- Improvement of circulation provides even distribution of the body weight and gives patient added overall comfort helping them to relax
- Breathable, stretchy cool-band material
- Sizing: Measure around the ankle above the malleoli. Choose the larger size if measurement is between sizes

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Ankle Circumference
PTA533327	GelBodies - Heel & Ankle Protector - Sml	20-25cm
PTA533327M	GelBodies - Heel & Ankle Protector - Med	25-30cm
PTA533327L	GelBodies - Heel & Ankle Protector - Lge	30-36cm