

## Frontier V4 Rear Wheel Drive Power Wheelchair - Hybrid



- Excellent Stability
- Exceptionally Smooth Ride
- Great battery life/distance
- Versatile indoor and outdoor performance
- Upgradable to the Off-Road with a simple tyre change
- Splash resistant (not waterproof)
- Please remember, the wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea. Magic Mobility wheelchairs meet Australian Standards 3696 and ISO7176.

✓ Contract approved

### Support @ Home

**Function:** Powered wheelchairs

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Safe Working Load	Seat To Floor Height	Seat To Footrest	Maximum Speed
PWC655025	Frontier V4 Rear Wheel Drive - Hybrid (As per script)	930 mm	670mm (710mm with off-road wheels)	182kg (155kg with lift/tilt)	435mm	435mm	10kph