


Frontier V4 Rear Wheel Drive Power Wheelchair - Hybrid

 Contract approved



magic
Mobility™

- Excellent Stability
- Exceptionally Smooth Ride
- Great battery life/distance
- Versatile indoor and outdoor performance
- Upgradable to the Off-Road with a simple tyre change
- Splash resistant (not waterproof)
- Please remember, the wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea. Magic Mobility wheelchairs meet Australian Standards 3696 and ISO7176.

Code	Display Name	Total Length	Total Width	Safe Working Load	Seat To Floor Height	Seat To Footrest	Maximum Speed
PWC655025	Frontier V4 Rear Wheel Drive - Hybrid (As per script)	930 mm	670mm (710mm with off-road wheels)	182kg (155kg with lift/tilt)	435mm	435mm	10kph