

## MOTomed Loop p.LA



**RECK**  
MOTomed®

- Insertion Aid: The insertion aid function facilitates the insertion of the feet into the foot shells before the start of the training. Using a button on the operating panel, the user can move the foot shells continuously into the appropriate position
- MovementProtector: If muscle cramp, spasticity or movement blockage (freezing) occurs during training, this function stops the pedal movement. It does not matter how the basic muscle tension changes; the MovementProtector does not remain constant, but automatically adapts to any change
- The new 7" big colour touch display can be operated easily and intuitively and offers intelligent functions of innovative form. Naturally, the display is suitable for wipe disinfection
- Includes Therapy programs, Games and videos and Slideshows



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
FPM101149	MOTomed Loop p.LA - Leg Guides - Stand - Gas Spring - Foot & Hand Fix	700 mm	600 mm	1070 - 1220 mm	135 kg	36 kg