

K Care Shower Commode Tilt



- Assists transportation and positioning of the user for showering and can be positioned over the toilet for toileting
- Gas tilt in space mechanism allows carer to adjust the vertical position from 0 to 30 (degrees) ensuring user comfortably retains correct seating position
- Large calf pads to support the legs and provide pressure care
- Large backrest with headrest extension to support the upper body and head
- Mounted on 4 locking castors for stability
- Stainless steel frame for corrosion resistance
- Available in a standard model (pictured); a deluxe version which comes with padded armrests and side panels; and an extra wide (50 cm) version
- Note: Pictured with optional handmade padded open seat (BTC075700)
- **SEATS, PANS & PAN-CARRIERS ARE SOLD AS SEPARATE ITEMS TO THE COMMODE**



Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Warranty Information

- 1 Year

Code	Display Name	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
BTC071200	Tilt in Space Commode - Swing Away Footrest - Padded Arms - 445 mm Wide	160 kg	23 kg	430 mm	445 mm	590 mm
BTC071350	Tilt in Space Commode - Swing Away Footrest - 500 mm Wide	160 kg	23 kg	430 mm	500 mm	590 mm