

Frontier V4 Front Wheel Drive Power Wheelchair - Hybrid



magic
Mobility™

- Great battery life/distance
- Versatile indoor and outdoor performance
- Upgradable to the Off-Road with a simple tyre change
- Splash resistant (not waterproof)
- Please remember, the wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea.
- Magic Mobility wheelchairs meet Australian Standards 3696 and ISO7176.

✓ Contract approved

Code	Display Name	Total Length	Total Width	Safe Working Load	Seat To Floor Height	Seat To Footrest	Maximum Speed
PWC655020	Frontier V4 Front Wheel Drive - Hybrid (As per script)	995 mm	670mm (710mm with off-road wheels)	182kg (155kg with lift/tilt)	435mm	435mm	10kph