

Frontier V4 Front Wheel Drive Power Wheelchair - Hybrid



- Great battery life/distance
- Versatile indoor and outdoor performance
- Upgradable to the Off-Road with a simple tyre change
- Splash resistant (not waterproof)
- Please remember, the wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea.
- Magic Mobility wheelchairs meet Australian Standards 3696 and ISO7176.

✓ Contract approved

Support @ Home

Function: Powered wheelchairs

Category: Prescribed

Click [here](#) for more information

| Code | Display Name | Total Length | Total Width | Safe Working Load | Seat To Floor Height | Seat To Footrest | Maximum Speed |
|-----------|---|--------------|------------------------------------|------------------------------|----------------------|------------------|---------------|
| PWC655020 | Frontier V4 Front Wheel Drive - Hybrid (As per script) | 995 mm | 670mm (710mm with off-road wheels) | 182kg (155kg with lift/tilt) | 435mm | 435mm | 10kph |