

Skippi Plus



ottobock.

Power wheelchair skippi Plus

Easy to transport

Travelling with a power wheelchair? No problem with the Skippi. You can take the wheelchair apart quickly and easily for transportation in the boot of a regular compact car

Skippi with Squiggles seating system

Various seating systems, such as the Squiggles system, are available as options for the Skippi in order to provide the child with optimum support while sitting.

SkippiPlus

The SkippiPlus is the wheelchair model with electric seat tilt, electric back angle adjustment and electric seat height adjustment. The enAble 50 control unit has three electric options, allowing the chair to be perfectly adapted to the child

SkippiPlus with special controls

Not all children are able to operate a power wheelchair with a standard controller. This is why the SkippiPlus can be equipped with numerous alternative special and environmental controls, tailored precisely to the needs of your child.

Warranty Information

- 2 Years

 Contract approved

Support @ Home

Function: Powered wheelchairs

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122306127_0105_1_2

Wheelchair Powered –

Composite:05_419_0105_1_2

Click [here](#) for more information

Code	Display Name	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Maximum Speed	Range
PWC652268	Skippi Plus	570 mm	350 - 460 mm	50 kg	68 kg	280 - 380 mm	300 - 380 mm	430 mm	6 km/h	30 km

