

Skippi Plus



ottobock.

Power wheelchair skippi Plus

Easy to transport

Travelling with a power wheelchair? No problem with the Skippi. You can take the wheelchair apart quickly and easily for transportation in the boot of a regular compact car

Skippi with Squiggles seating system

Various seating systems, such as the Squiggles system, are available as options for the Skippi in order to provide the child with optimum support while sitting.

SkippiPlus

The SkippiPlus is the wheelchair model with electric seat tilt, electric back angle adjustment and electric seat height adjustment. The enAble 50 control unit has three electric options, allowing the chair to be perfectly adapted to the child

SkippiPlus with special controls

Not all children are able to operate a power wheelchair with a standard controller. This is why the SkippiPlus can be equipped with numerous alternative special and environmental controls, tailored precisely to the needs of your child.

Warranty Information

- 2 Years



ndis

Item Code:

05_122306127_0105_1_2

Wheelchair Powered –

Composite:05_419_0105_1_2

Click [here](#) for more information

| Code | Display Name | Seat Width | Total Width | Safe Working Load | Maximum Speed | Range |
|-----------|--------------|--------------|-------------|-------------------|---------------|-------|
| PWC652268 | Skippi Plus | 300 - 380 mm | 570 mm | 50 kg | 6 km/h | 30 km |

