

Wet & Dry Trolley



- Mobile moist linen trolley
- Transport and dispose soiled linen in the hospital or aged care environment
- Galvanized plus zinc rich epoxy coated frame
- 125 mm heavy duty nylon swivel castors
- Spring loaded base is ideal to assist in occupational health & safety
- Also available with heavy duty spring coil base
- Optional push handle and spring cover available for medium & large options
- Pictured with Hanging Coil Base

Warranty Information

• 3 Years

Ocontract appr	
Contract appr	oved

Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
LAC353380	Coil Spring Base - Small	800 mm	600 mm	800 mm
LAC354190	Coil Spring Base - Medium	1050 mm	700 mm	850 mm
LAC354210	Coil Spring Base - Large	1250 mm	700 mm	850 mm
LAC354220	Push Handle - Bolt On for Wet and Dry Linen Trolley - Medium	-	-	-
LAC354225	Push Handle - Bolt On for Wet and Dry Linen Trolley - Large	-	-	-
LAC354215	Push Handle - Small	-	-	-
LAC354227	Spring - Wet and Dry Laundry Trolley - each	-	_	-



Code	Display Name	Total Le	ength Total Wid	lth Total	Height
LAC354066 Spring cove	r - Wet and Dry Laundry Trolley	-	-	-	