

Safelife Bed Exit PIR Alarm System



- PIR Beam Detection System
- Bed Exit PIR will detect as a person swings their legs over the side of the bed breaking the invisible PIR beam
- The beam is designed to taper down, ending at the base of the bed
- Easy Battery Replacement
- Portable and compact in size

✔ Contract approved

Support @ Home

Function: Signalling, alarming and localisation

Category: Prescribed

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name
BEA007760	Safelife Bed Exit PIR Alarm System - Single
BEA007770	Safelife Bed Exit PIR Alarm System - Dual