

TiLite Aero T



permobil

TiLite

- **Pure geometry:** Our most advanced aluminium chair, the Aero T's dual-tube frame minimises weight while maximising performance. Add to that full adjustability and TiFit, and you've got something special.
- Rigid Dual-Tube frame chairs are the lightest and highest performing chairs because they are constructed from the fewest possible components and have the fewest moving parts.
- Ease of Transport and Energy Saving: no wasted energy with a swept-in dual-tube frame. Dual-tube geometry adds rigidity to the chair, transferring more of your push stroke to the wheel, meaning less push strokes throughout the day.
- The centre of gravity can be moved forward, backward, up and down, and all backrest options are adjustable to meet your specific needs. Precision Lock components make these adjustments easier and faster than ever.
- The Precision Lock Caster is a simplified and more accurate way to square your wheelchair in under a few minutes. It is the only caster system that can be adjusted by a single bolt with infinite adjustment positions.
- Precision Lock Backrest combines a minimal design with matched adjustability. Your backrest can now be adjusted in a matter of minutes. The Precision Lock Backrest offers infinite adjustments to meet your needs and can be made while the user sits in the chair.
- The Precision Lock Caster is a simplified and more accurate way to square your wheelchair in under a few minutes. It is the only caster system that can be adjusted by a single bolt with infinite adjustment positions.
- The Ergonomic Advantage means an additional unique bend in the frame creating an ergonomic seat slope sitting you lower in the chair providing greater pelvic stability and improved posture through optimal contouring.
- Zero Play Camber Plugs give a much tighter feel to the chair like sport suspension in a car. CNC machined from billet aluminium to the highest

✓ Contract approved

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122203111_0105_1_2

Wheelchair Manual –

Composite:05_417_0105_1_2

Click [here](#) for more information

tolerances, they remove any side to side movement from the rear wheels.

Code	Display Name	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
MWC427471	Aero T	120 kg	5.4 kg	30.5 - 51 cm	30.5 - 51 cm	38 - 53 cm