

## **MOTOmed Wheelchair Stabiliser**



- Recommended for severe spasticity
- Supports the wheelchair: prevents tipping and slipping backwards
- For all common wheelchairs with push handle
- With padded support for the wheelchair push handles
- Width of the handle support: 580 mm Children's stabiliser especially designed for children's wheelchairs with push handle: Width of the handle support: 330 mm



## RECK MOTOmed<sub>®</sub>

Code	Display Name

FPM102170 MOTOmed Wheelchair Stabiliser

FPM102175 MOTOmed Wheelchair Stabiliser - Children