

ROHO Adaptor Pad



permobil

ROHO

- The self-contained air cells in the adaptor Pad can be used as an interface on virtually any flat or contoured surface to protect skin and soft tissues.
- Especially useful when cut to fit various surface shapes such as sides of wheelchairs, trays, head and foot rests, inside of helmets, or for any surfaces that may need protection
- Absorbs shock
- Protects against shear and friction
- Easily cleaned and disinfected
- Durable
- Closed air cells, so product maintains efficacy even if a single cell is cut or punctured
- Made of flame-resistant, black neoprene (not made with natural rubber latex)

 Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name
------	--------------

PTA507400	ROHO Adaptor Pad
-----------	------------------