


Snooooooze Curved Body Wedge



- Velcro receptive wedges adhere firmly to the Snooooooze mattress cover and are placed alongside the body/torso when the client is lying supine
- These wedges can be used either on their own as essential support for the torso when lying supine, or as an additional part of a sleep system adding support along the body contour as required



Code	Display Name	Total Length	Total Width
PTA534205	Snooooooze Curved Body Wedge - Large	300 mm	160 mm
PTA534206	Snooooooze Curved Body Wedge - Paediatric	100 mm	80 mm
PTA534207	Snooooooze Curved Body Wedge - Small	200 mm	150 mm