

Cura Alarms - Patient Activated Monitors



- Patient activated alert monitor assists with falls prevention
- Attached to patients clothing, patient activates by pulling string which detaches magnet and sets off alarm
- Stand alone system is simple to use
- Durable and reliable
- Portable and lightweight
- Supplied with 9 volt battery

Warranty Information

• 1 Year

\oslash

ndis

Item Code: 05_222704219_0123_1_2 Personal Alarm – Composite: 05_427_0103_1_2 Click <u>here</u> for more information

Compliance

No applicable compliance standard

cura¹

Code Display Name

BEA007460 Pull String Monitor