

One Way Velvet Transfer Slide



- One Way Slides are designed to easily assist in the correct positioning of a person in the back of a chair and to prevent them slipping forward
- The sliding surfaces use a combination of low friction and strips of easy sliding one way materials
- Single layer one way slide for use to help position a person on pressure cushions if they do not have severe decubitus problems

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

Code

Display Name

TSE672521 One Way Velvet Transfer Slide