

Action 2NG



- Lightweight foldable frame.
- Cross-compatibility allows users to update and refurbish their wheelchairs with other models in the Action NG family.
- 80° footrest swings in and out for easy transfers in or out of the chair.
- Axle position can be adjusted in height by 5cm.

 Contract approved

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
MWS449196	Wheelchair - Invacare - Action 2NG Transit - 45x42.5	125 kg	8.9 - 12.1 kg	400 - 430 mm	405 - 480 mm	485 mm
MWS449199	Wheelchair - Invacare - Action 2NG Self Propelled - 45x42.5	125 kg	8.9 - 12.1 kg	400 - 430 mm	405 - 480 mm	485 mm