

## Meiko TopLine 30 Sanitiser



- Packed with powerful technology, the TopLine bedpan washer-disinfectors aim to eliminate the risk of infections stemming from hospital dirty utility rooms
- Wall Embedded
- Increased Capacity & System-redundancy if failure event
- Patented complete steam impermeability system – No risk of germ discharge, highest level of work and hygiene safety
- Controllable main telescopic rotary nozzle – Optimised internal cleaning of the care utensils
- Thorough cleaning with power nozzle system
- Air drying and cooling with filtered air (HEPA filter) produces hygienically clean results
- Replaceable and simple to change HEPA filter
- System disinfection of all water-carrying lines and protects from recontamination
- Patented control of the freely selectable A0 value (60 to 3000) and monitored by temperature sensors which crosscheck each other's readings guarantees disinfection safety
- Hygienic safety due to A0-value control system to DIN EN ISO 15883
- A0 value in the display (actual and target value)
- Supply disconnection to EN 1717 type AA
- Integrated steam generator/water tank, which means shortest possible programme cycles
- Disinfection process control with temperature sensors. The control software monitors and assesses the control process
- Maximum flexibility due to multifunctional utensils holder concept



<b>Code</b>	<b>Display Name</b>	<b>Total Length</b>	<b>Total Width</b>	<b>Total Height</b>
UPM200160	Meiko Sanitiser - TopLine 30	450 mm or 588 mm	650 mm	2210 mm
UPM200170	Meiko Sanitiser - TopLine 30 AT - with Auto Door	450 mm or 588 mm	650 mm	2210 mm
UPM200180	Meiko Sanitiser - TopLine 30 WC	450 mm or 588 mm	650 mm	2210 mm
UPM200190	Meiko Sanitiser - TopLine 30 WCAT - with Auto Door	450 mm or 588 mm	650 mm	2210 mm