

Alerta - Wireless - Chair Alertamat



- Alerts carers if patient leaves or attempts to leave chair, helping to manage falls risk
- Simple plug and play ensures carer can quickly set up and operate with ease
- High and low volume alert
- Multiple choice alarm tone
- Adjustable delay to reduce false alarms
- Alarm is reset automatically when patient returns to the chair
- Hygienic, antibacterial easy-clean material
- Works in conjunction with the Alerta wireless products and requires the AlertaDetect Receiver to work
- Requires the Alerta Wall Point Receiver to work - sold separately

 Contract approved

Support @ Home

Function: Signalling, alarming and localisation

Category: Prescribed

Click [here](#) for more information



Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width
BEA007657	Alerta - Wireless - Chair Alertamat (mat only)	380 mm	250 mm