

Positioning Slide - Foot Operated



- Designed with Care staff in mind, low friction surfaces allow for easy patient manoeuvring with minimal lift
- For patients with a lower risk of slipping while sitting in a chair
- Manoeuvre the patient easily to the back of the chair or sideways for transfer assistance
- Foot strap threads through back of chair to assist the carer with positioning
- Machine washable

Warranty Information

- 1 Year



Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123103151_0105_1_2

Minor Transfer –

Composite:05_430_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width
TSE673054	Positioning Slide	1375 mm	460 mm