


Positioning Slide - Foot Operated



- Designed with Care staff in mind, low friction surfaces allow for easy patient manoeuvring with minimal lift
- For patients with a lower risk of slipping while sitting in a chair
- Manoeuvre the patient easily to the back of the chair or sideways for transfer assistance
- Foot strap threads through back of chair to assist the carer with positioning
- Machine washable

Warranty Information

- 1 Year

 Contract approved

ndis

Item Code:

05_123103151_0105_1_2

Minor Transfer –

Composite:05_430_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

| Code | Display Name | Total Length | Total Width |
|-----------|-------------------|--------------|-------------|
| TSE673054 | Positioning Slide | 1375 mm | 460 mm |