


Alerta - Wireless - Detect Sensor



- The Wireless Alerta Detect uses infrared technology to activate the alarm when body movement is detected within the range of the infrared beam
- Most effective when mounted at the door or bedside; activating the alarm when a user leaves their safe zone and breaks the infrared beam
- Alerts the caregiver to ascertain the user's safety
- Ideal for providing users with the freedom and safety of their own space while still proactively giving the caregiver advanced warning when a user leaves their safe zone
- Can be connected to a nurse call system through the Alerta Wall Point Receiver for centralised activation, or used on its own, independent of a nurse call system to give localised alerts through the Alerta Wall Point Receiver only
- Wireless connect-and-play system
- No cables, extra safe and versatile
- Must be used with Alerta - Wireless/AlertaDetect - Alarm Monitor - sold separately

 Contract approved

Compliance

No applicable compliance standard

| Code | Display Name | RRP | Box Quantity |
|-----------|-----------------------------------|-----|--------------|
| BEA007675 | Alerta - Wireless - Detect Sensor | 359 | 1 |