

Cadex VibraPlus 8 Vibrating Watch



e-pill | CADEX VibraPlus 8 Daily Alarm Vibrating Watch The CADEX VibraPlus is a discreet and reliable medication reminder watch featuring up to 8 daily alarms with vibration, sound, or both. It looks like a regular wristwatch and is ideal for anyone who needs timely medication reminders or discreet alerts during the day.



- **Multi-Alarm Function:** Set up to 8 alarms per day with automatic daily repeat.
- **Discreet Medication Reminder:** Perfect for individuals who prefer privacy or have hearing impairments.
- **Durable Design:** Black/silver polycarbonate case (40mm × 50mm) with stainless steel back.
- **Clear Display:** Easy-to-read time digits (6.35mm high × 3mm wide) with backlight for nighttime visibility.
- **Comfortable Strap:** Replaceable black leather standard 18mm strap; fits wrists approx. 16–22 cm.
- **Strong Vibration Alerts:** Choose vibration, sound, or both for flexible alert options.
- **Countdown Interval Timer:** Adjustable from 15 minutes to 24 hours, with auto-start between set times each day.
- **Key Lock Feature:** Prevents accidental setting changes.
- **Stopwatch & Chronograph:** Includes stopwatch functionality for timing activities.
- **Battery Information:** Replaceable CR2032 3V battery with 4–6 month life.
- **Spare Battery Included:** Comes with 1 spare battery and screwdriver for easy replacement.
- **Water Resistant:** Suitable for everyday use (not for swimming or bathing).
- **User Manual:** Includes an easy-to-read instruction guide.
- **Warranty:** One-year manufacturer warranty (battery excluded).

- **The CADEX VibraPlus:** It helps prevent missed doses and medication errors, ensuring better compliance and peace of mind.
- **Battery Safety Warning:** This product contains a button or coin cell battery. Hazardous if swallowed. Keep away from children. Can cause severe or fatal injuries in 2 hours or less if swallowed or placed inside any part of the body. Seek medical attention immediately if ingestion or insertion is suspected. Contact the Australian Poisons Information Centre on 13 11 26 for 24/7 expert advice. Dispose of used button/coin batteries immediately by wrapping them in sticky tape and placing them in an outside bin, out of reach of children, or recycle safely.

Code	Display Name
DLG279000	Cadex VibraPlus 8 Alarm Vibrating Medication Reminder Watch