


Frontier V4 Rear Wheel Drive Power Wheelchair - All Terrain

 Contract approved



magic
Mobility™

- Compact with great off-road capability
- Powerful kerb climbing
- Excellent stability
- Exceptionally smooth ride
- Great battery life/distance
- Splash resistant (not waterproof)
- Please remember, the wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea. Magic Mobility wheelchairs meet Australian Standards 3696 and ISO7176.

Code	Display Name	Total Length	Total Width	Safe Working Load	Seat To Floor Height	Seat To Footrest	Maximum Speed
PWC655015	Frontier V4 Rear Wheel Drive - All Terrain (As per script)	930 mm	710mm (670mm with Hybrid wheels)	182kg (155kg with lift/tilt)	435mm	435mm	10kph