

MOTomed Letto2



RECK
MOTomed®

The MOTomed letto2 is ideal for early mobilisation from intensive care. Patients can use the MOTomed letto2 to train the lower extremities in a passive, motor-assisted, and active mode, from a hospital bed or a therapy chair. It provides a well-rounded therapy to patients in stationary care.



- Leg trainer with pedal radius 7.0 cm
- Big handle bar for an easy slide and positioning
- Chassis with big castors and manual width adjustment
- Bed mount for a stable hold and a safe training
- Pneumatic height adjustment
- Horizontally adjustable toward the patient
- rotating operating panel
- MOTomed SmoothDriveSystem with a gentle startup and rundown
- MovementProtector & SpasmControl
- Speed regulation from 0 to 60 rpm (passive)
- Regulation of motor power from 1 to 10
- Preset a therapy duration from 0 to 120 minutes
- Big, user-friendly operating panel; colour code system and big, easy-feel buttons
- Language selection between 27 different languages

Code	Display Name
FPM101118	MOTomed Letto2 (Supine / recumbent model)