

ROHO - Enhancer



permobi

ROHO

The ENHANCER Cushion uses positioning in a uniquely designed two-manifold system for enhanced midline channeling of the femurs, lateral stability and skin/soft tissue protection.

- The cellular design and the soft pliable material allow the individual to be immersed into the cushion for maximum pressure redistribution, while decreasing friction and shear
- Multiple Cell Heights: Unique contoured surface provides for channeling of the legs and enhances lateral stability and body positioning
- Adjustable and Adaptable: ROHO's DRY FLOATATION® Technology provides unparalleled performance in minimizing tissue deformation and facilitating blood flow by fitting, matching and adapting to the shape of an individual
- Positioning: The ENHANCER Cushion provides positioning with a multiple cell
 height, two-compartment system that allows the individual to adjust the
 interior bucket seat section for reduced pressure on the ischia while adjusting
 the two outer perimeter rows for enhanced stability and positioning
- Skin/Soft Tissue Protection: Proven with over 40 years of research, experience and refinement, ROHO DRY FLOATATION Technology provides an individual with the ultimate in skin and soft tissue protection. The ENHANCER Cushion provides an adjustable, low friction/shear environment that helps protect against pressure ulcers
- A unique combination of high (4"), mid (3") and low (2") profile cells in the one cushion centre the user and position the pelvis back in the wheelchair providing greater lateral stability through chanelling of the legs
- A two-valve system allows the user to (1) adjust the interior "bucket" section for reduced pressure on the ischia and coccyx; and (2) adjust the outer perimeter for enhanced stability and positioning
- Appropriate for clients with high risk of ischemic ulcers and need for greater

Contracted DVA product

⊘Contract approved

ndis

Item Code:

05_181006190_0105_1_2

Pressure Care Cushion –

Composite:05_429_0103_1_2

Click here for more information

Compliance

EN-ISO 10993-10 (1999-11) / EN 1021-1 and EN 1021-2



- positioning support
- Soft neoprene construction is flame retardant and durable
- Maximum User Weight: Unlimited
- All models include: 2 way stretch cover, Hand pump & Repair kit, Operating Instructions.

Warranty Information

2 Years

Code	Display Name	Total Height	Safe Working Load	Pressure Care Rating	Material
PTC553300	ROHO Enhancer - 15 x 15in - (8 x 8cell)	2.25 - 4.25 in	Unlimited	High Risk	Air Cells
PTC553400	ROHO Enhancer - 15 x 16.75in - (8 x 9cell)	2.25 - 4.25 in	Unlimited	High Risk	Air Cells
PTC553550	ROHO Enhancer - 16.75 x 15in - (9 x 8cell)	2.25 - 4.25 in	Unlimited	High Risk	Air Cells
PTC553600	ROHO Enhancer - 16.75 x 16.75in - (9 x 9cell)	2.25 - 4.25 in	Unlimited	High Risk	Air Cells
PTC553500	ROHO Enhancer - 16.75 x 18.5in - (9 x 10cell)	2.25 - 4.25 in	Unlimited	High Risk	Air Cells
PTC553200	ROHO Enhancer - 18 x 16.75in - (10 x 9cell)	2.25 - 4.25 in	Unlimited	High Risk	Air Cells
PTC553100	ROHO Enhancer - 18.5 x 18.5in - (10 x 10cell)	2.25 - 4.25 in	Unlimited	High Risk	Air Cells
PTC553090	ROHO Enhancer - 18.5 x 20.25in - (10 x 11cell)	2.25 - 4.25 in	Unlimited	High Risk	Air Cells
PTC553060	ROHO Enhancer - 20.25 x 18.5in - (11 x 10cell)	2.25 - 4.25 in	Unlimited	High Risk	Air Cells