

## Handi Lift Walk Belt



- The Handi Lift Walk Belt is a popular transfer belt for patient handling. You can use it for transferring people, and when walking with people who need help
- Two Carers normally work together when using the Handi Lift Walk Belt
- The 10 cm / 4 inch wide double Hook & Loop Tape belt is very supportive. You can fit the Belt around the waist on ladies, or under the arms of men and some women
- All four of the webbing handles are padded
- Wipe clean with a hospital grade disinfectant or wash up to 80°C (176°F). Air or tumble dry on a cool setting, ensuring the drum is cool. Ensure the product is dry before using or storing
- Different colours for different sizes
- Total length in specs are the measurement for the waist

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Changing body position

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Total Length
TSE673110	Handi-Lift Walk Belt - Medium	850 mm - 110 mm
TSE673130	Handi-Lift Walk Belt - Large	900 mm - 1170 mm
TSE673150	Handi-Lift Walk Belt - Small	690 mm - 930 mm
TSE673140	Handi-Lift Walk Belt - Extra Large	1080 - 1330 mm
TSE673155	Handi-Lift Walk Belt - X Small	580 - 700 mm
TSE673142	Handi-Lift Walk Belt - XX Large	-

