

Cura Alarms - Safetime Kit



cura¹

- Bed/chair pad alarm system activates when a resident or patient falls or gets out of a bed or chair
- Easily activated switch for the caregiver to switch the unit on/off or to reset
- System activates when a resident or patient removes their weight from the connected sensor pad
- Re-arms itself if the resident / patient gets back onto the sensor pad
- Ultra sensitive and durable sensors reduce false alarms to nil over the lifetime of the pad
- Anti allergic soft PVC construction for optimum comfort
- Venting facilities reduce the risk of moisture ingress
- Supplied with 9 volt battery
- Optional DV plug pack for permanent (240 volt) power source
- Available as bed sensor pad, chair pad or floor mat

Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Signalling, alarming and localisation

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_222704219_0123_1_2

Personal Alarm – Composite:

05_427_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Safe Working Load
BEA007422	Safetime - Bed Pad	510 mm	760 mm	200 kg

Code	Display Name	Total Length	Total Width	Safe Working Load
BEA007420	Bed Monitor	-	-	-