

## Action 3 Junior EVO



- The seat width of the chassis can be adjusted with an innovative telescopic crossbar over 4 seat positions from 205 to 380 mm
- As the chassis is made wider, so the upholstery on the back and on the seat of the chair can be made wider with Velcro flaps
- In addition to the standard slide-on table tray, we offer a swing-away tray which can easily be folded away
- For improved driving performance and different height options, the rear wheel can be mounted in 20 different positions
- A front stabilizer kit is available to increase forward stability for active children who are likely to climb all over the wheelchair and be especially energetic when propelling the wheelchair
- Four different types of spoke guards are offered to match the different wheelchair colours and tastes

### Warranty Information

- Frame & Crossbaraces: 2 Years
- Components: Limited 1 Year

 Contract approved

### Support @ Home

**Function:** Manual wheelchairs

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_122203114\_0105\_1\_2

**Wheelchair Manual –**

**Composite:**05\_417\_0105\_1\_2

Click [here](#) for more information

### Compliance

ISO 7176-7:1996

Code	Display Name	Safe Working Load	Seat Depth	Seat Width
MWS449180	Action 3 Junior EVO - 11in to 13in	80 kg	11 - 13 in	12 - 15 in
MWS449181	Action 3 Junior EVO - 10in to 12in	80 kg	10 - 12 in	8 - 11 in
MWS449182	Action 3 Junior EVO - 12in to 14in	80 kg	12 - 14 in	8 - 11 in

Code	Display Name	Safe Working Load	Seat Depth	Seat Width
MWS449183	Action 3 Junior EVO - 8in to 10in	80 kg	8 - 10 in	8 - 11 in