

Invacare Hammock Sling



- The sling provides support and security for lifting from a bed or the floor
- The patient must roll on to, and off, the sling
- When fitted full head support is provided
- Recommended for single or double lower limb amputees and lifting from the floor, and because weight is not taken by the patient's thighs



ndis

Item Code:

05_123603164_0105_1_2

Sling – Composite:

05_421_0105_1_2

Click [here](#) for more information

Warranty Information

- 1 Year



Code	Display Name	Safe Working Load	Material
LSS375015	Invacare Hammock Mesh - Medium	205 kg	Mesh

Compliance

AS/NZS ISO 10535:2011