

Bed Alertamat System




- Provides advanced warning of falls risk if patient leaves or attempts to leave bed
- Simple plug and play ensures carer can quickly set up and operate with ease
- High and low volume alert
- Multiple choice alarm tone
- Adjustable delay to reduce false alarms
- Alarm is reset automatically when patient returns to the chair
- · Hygienic, antibacterial easy-clean material
- System includes: Bed Pad, Alarm Monitor, Stereo & Mono Nurse Call Cables
- Mains power cable is available as an optional accessory & batteries are not included

Warranty Information

• 1 Year

Ω		
Ocontracted	DVA	product



Compliance

No applicable compliance standard



Code	Display Name	RRP	Box Quantity	Total Length	Total Width
BEA007605	Alarm - Alerta - Bed Alertamat System (includes Pad, Monitor & Mono/Stereo Cables)	356	1	760mm	250mm