

**Product Awareness Bulletin – April 2018**

**NEW PRODUCT RANGE!! 2<sup>nd</sup> Generation Aspire Wheelchairs**

The Aspire Manual Wheelchair Range that you have come to know and trust just got better.

You will find the softcopy range brochure attached or find it on our website at [www.aidacare.com.au](http://www.aidacare.com.au). This PAB covers the wheelchairs only so please see the brochure for detailed specs and accessories.

The new range supersedes the existing (Aspire Evoke, Assist, Transit & Relax) but we have retained the same names for Evoke, Assist, and Transit with “2” added for 2nd generation. The new TIS chairs are very different to the Relax so we have a new name - The Rehab RX and Rehab RS... Details below.

These changes bring feature enhancement and range expansion.

**Aspire Evoke2 & Evoke2 HD (DVA Contracted)**

These are highly adjustable, lightweight, aluminium wheelchairs designed to enable a personalised fit for users.

A fantastic accessory to the new Evoke2 is the **Contoured Padded Overlay**, please refer to the attached flyer for details.

Important product improvements that have been made to this product are:

- **NOW FULLY CRASH TESTED** – can be used in approved vehicles for patient transport.
- SWL upgrade in non-HD from 130kg to 140kg.
- SWL for HD chairs is 180kg – includes crash tested, greater width range and more features and accessories.
- Increased levels of adjustment in the seat height, depth, backrest height & angle.
- Expanded range of performance enhancing accessories including the new contoured padded overlay.

Codes and Sizes below – including comparable old chairs being superseded:

<b>New Code</b>	<b>Description</b>	<b>Old Code</b>	<b>Description</b>
MWS449500	Aspire Evoke2 - 350mm - 140kg - Silver	MWS440100	Aspire Evoke - 350mm - 130kg - Silver
MWS449510	Aspire Evoke2 - 400mm - 140kg - Silver	MWS440200	Aspire Evoke - 400mm - 130kg - Silver
MWS449520	Aspire Evoke2 - 450mm - 140kg - Silver	MWS440300	Aspire Evoke - 450mm - 130kg - Silver
MWS449530	Aspire Evoke2 - 500mm - 140kg - Silver	MWS440500	Aspire Evoke - 500mm - 130kg - Silver
MWS449540	Aspire Evoke2 HD - 570mm - 180kg - Silver	MWS440525	Aspire Evoke HD - 560mm - 190kg - Silver
MWS449550	Aspire Evoke2 HD - 610mm - 180kg - Silver	MWS440550	Aspire Evoke HD - 580mm - 190kg - Silver



## Aspire Evoke2 JNR (NEW!)

The Aspire Evoke2 JNR (Junior) is a paediatric version of the adult Evoke2 chairs. The frame is reduced in length and width while all the adjustable features of the Evoke2 are retained to suit the unique needs of children, with key features such as:

- Smaller 20in rear self propelling wheels
- Lightweight design promotes independent manoeuvring and handling by juniors
- Lightweight also makes handling a breeze for parents and caregivers
- **FULLY CRASH TESTED** – can be used in approved vehicles for patient transport

The Aspire Evoke2 JNR is available in one size:

New Code	Description
MWS440100	Aspire Evoke2 JNR - 300mm - 140kg - White



## Aspire Rehab - Tilt in Space Wheelchairs

Three new Tilt in Space wheelchairs in multiple sizes replace the old Aspire Relax TIS chairs.

1. The **Aspire Rehab RX** and **Aspire Rehab RX JNR** are very highly configurable chairs geared towards Complex Rehab clients.
2. The **Aspire Rehab RS** is a medium configurable chair designed for Aged Care and Community use

### 1. Aspire Rehab RX and Rehab RX JNR (NEW!) (DVA Contracted)

The **Aspire Rehab RX** (in black frame) is geared towards the high needs rehabilitation environment for users with more complex needs. **The Rehab RX JNR** (in white frame) has all the features of the Rehab RX reduced in size to suit the unique needs of children.

On top of extensive adjustability, **Rehab RX** chairs also come with specially contoured pressure care surfaces and the most advanced articulating elevated legrests for optimum pressure, positioning and postural support.

Key features:

- **FULLY CRASH TESTED** – can be used in approved vehicles for patient transport.
- Increased SWL from 136kg (old Aspire Relax) to 180kg!
- Multilayer contoured support surfaces with castellated top layer for high level comfort and pressure reduction.
- Breathable multi stretch covers for enhance pressure redistribution and skin microclimate.
- Advanced articulating Elevating Legrests (ELR's) that mirror the anatomical range of movement at the knee joint, ensuring thighs do not elevate as legs are raised.
- Extensive accessories and options as well as the ability to configure after market backrests and cushions.

Codes and Sizes below – including comparable old chairs being superseded:

New Code	Description	Old Code	Description
MWS449710	Aspire Rehab RX - 410-460mm - Black	MWS442025	Aspire RELAX - 440-460mm - Silver
MWS449720	Aspire Rehab RX - 460-510mm - Black	MWS442050	Aspire RELAX - 490-510mm - Silver
MWS449740	Aspire Rehab RX - 510-560mm - Black		
MWS449630	Aspire Rehab RX JUNIOR - 320-370mm - White		



## **2. Aspire Rehab RS (NEW!) (DVA Contracted)**

The Aspire Rehab RS is an ideal solution for users requiring fully adjustable postural support and regular bodyweight redistribution.

The Gunmetal Grey sub-frame is the same as the Rehab RX, but the Rehab RS comes with less advanced support surfaces and ELR's allowing it to achieve a pricepoint for Aged Care, Community use.

Key features:

- **FULLY CRASH TESTED** – can be used in approved vehicles for patient transport.
- Increased SWL from 136kg (old Aspire Relax) to 180kg!
- Standard foam comfort cushions on back and seat
- Elevating Legrests with calf pad
- Extensive accessories and options as well as the ability to fit after market backrests and cushions
- 3 sizes – covering seat widths from 410mm to 560mm

Codes and Sizes below – including comparable old chairs being superseded:

New Code	Description	Old Code	Description
MWS449600	Aspire Rehab RS - 410-460mm - Gunmetal Grey	MWS442025	Aspire RELAX - 440-460mm - Silver
MWS449610	Aspire Rehab RS - 460-510mm - Gunmetal Grey	MWS442050	Aspire RELAX - 490-510mm - Silver
MWS449620	Aspire Rehab RS - 510-560mm - Gunmetal Grey		

## **Aspire Assist2 (DVA Contracted)**

The Aspire Assist2 is a versatile and robust self propelling wheelchair, designed to meet the rigorous demands of Community, Hire, Hospital & Aged Care environments.

Key features:

- Increased SWL to 140kg
- Caters to a wide variety of environments with large puncture-proof tyres that handle varied terrain with ease
- Strong steel frame, removable swing-away legrests and quick release self-propelling wheels
- Low maintenance folding design.

Codes and Sizes below – including comparable old chairs being superseded:

<b>New Code</b>	<b>Description</b>	<b>Old Code</b>	<b>Description</b>
MWS449300	Aspire Assist2 - 350mm - 140kg - Silver	MWS441400	Aspire ASSIST - 350mm - 125kg - Silver
MWS449310	Aspire Assist2 - 400mm - 140kg - Silver	MWS441500	Aspire ASSIST - 400mm - 125kg - Silver
MWS449320	Aspire Assist2 - 450mm - 140kg - Silver	MWS441600	Aspire ASSIST - 450mm - 125kg - Silver
MWS449330	Aspire Assist2 - 500mm - 140kg - Silver	MWS441700	Aspire ASSIST - 500mm - 125kg - Silver

## **Aspire Assist2 Ward Vinyl (NEW!) (DVA Contracted)**

The Aspire Assist2 Ward Vinyl wheelchair has all the features of the Assist2, but with durable ward grade vinyl, suitable for the Hospital & Clinical environment. Fully welded seams and concealed mounting of upholstery ensure that rigorous infection control standards can be met.

Key features:

- Larger range of seat widths available (350-610mm) to cater for Hospital Ward environments
- Increased SWL (140kg for 350-500mm & 180kg for 570-610mm)
- Infection control focus:
  - Durable hospital vinyl that is easy to wipe clean
  - Fully welded seams to inhibit bacterial growth
  - Concealed mounting of upholstery to prevent fluid ingress

New codes are as follows:

<b>New Code</b>	<b>Description</b>
MWS449340	Aspire Assist2 - Ward Vinyl - 350mm – 140kg - Silver
MWS449350	Aspire Assist2 - Ward Vinyl - 400mm – 140kg - Silver
MWS449360	Aspire Assist2 - Ward Vinyl - 450mm - 140kg - Silver
MWS449370	Aspire Assist2 - Ward Vinyl - 500mm - 140kg - Silver
MWS449380	Aspire Assist2 - Ward Vinyl - 570mm - 180kg - Silver
MWS449390	Aspire Assist2 - Ward Vinyl - 610mm - 180kg - Silver

**Aspire Assist2 and Assist2 Ward Vinyl**



## **Aspire Transit2 (DVA Contracted)**

The Aspire Transit2 wheelchair is a versatile and robust attendant propelled transport wheelchair. It is a reliable and comfortable solution, with various accessories available to suit the needs of individuals and carers in the Community, Hospital or Aged Care environments. Key features are:

- Increased SWL to 140kg
- Strong steel frame, removable swing-away legrests and quick release transit wheels
- Low maintenance folding design
- Comfortable padded nylon upholstery



Codes and Sizes below – including comparable old chairs being superseded:

<b>New Code</b>	<b>Description</b>	<b>Old Code</b>	<b>Description</b>
MWS449410	Aspire Transit2 - 400mm - 140kg - Silver	MWS441800	Aspire TRANSIT - 400mm - 125kg - Silver
MWS449420	Aspire Transit2 - 450mm - 140kg - Silver	MWS441900	Aspire TRANSIT - 450mm - 125kg - Silver
MWS449430	Aspire Transit2 - 500mm - 140kg - Silver	MWS442000	Aspire TRANSIT - 500mm - 125kg - Silver

## **Aspire Lite & Lite Transit (DVA Contracted)**

There has been no change to the two Aspire LITE wheelchairs and both are still available and in stock:

<b>Code</b>	<b>Description</b>
MWS442300	Wheelchair Folding - Aspire LITE - 460mm - 120kg - Blue
MWS442301	Wheelchair Folding - Aspire LITE TRANSIT - 460mm - 120kg - Green

